

Product datasheet for **RC401975**

NF2 (NM_000268) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	NF2 (NM_000268) Human Mutant ORF Clone
Mutation Description:	Q428X
Affected Codon#:	428
Affected NT#:	1282
Nucleotide Mutation:	NF2 Mutant (Q428X), Myc-DDK-tagged ORF clone of Homo sapiens neurofibromin 2 (merlin) (NF2), transcript variant 1 as transfection-ready DNA
Effect:	Neurofibromosis 2
Symbol:	NF2
Synonyms:	ACN; BANF; merlin-1; SCH
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_000268
ORF Size:	1281 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC401975 representing NM_000268
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCGGGGCCATCGCTTCCCGCATGAGCTTTCAGCTCTCTCAAGAGGAAGCAACCCAAGACGTTACCCG
 TGAGGATCGTACCATGGACGCCGAGATGGAGTTCAATTGCGAGATGAAGTGAAAGGGAAGGACCTCTT
 TGATTTGGTGTGCCGACTCTGGGGCTCCGAGAAACCTGGTTCTTTGGACTGCAGTACACAATCAAGGAC
 ACAGTGGCCTGGCTCAAATGGACAAGAAGGACTGGATCATGATGTTTCAAAGGAAGAACCAGTCACCT
 TTCACTTCTGGCCAAATTTTATCCTGAGAATGCTGAAGAGGAGCTGGTTCAGGAGATCACACAACATTT
 ATTCTTCTACAGGTAAGAAGCAGATTTAGATGAAAAGATCTACTGCCCTCCTGAGGCTTCTGTGCTC
 CTGGCTTCTTACGCCGTCAGGCCAAGTATGGTACTACGACCCAGTGTTCACAAGCGGGGATTTTGG
 CCCAAGAGGAATTGCTTCCAAAAGGGTAATAAATCTGTATCAGATGACTCCGGAATGTTGGGAGGAGAG
 AATTACTGCTTGGTACGCAGAGCACCGAGGCCAGGGATGAAGCTGAAATGGAATATCTGAAGATA
 GCTCAGGACCTGGAGATGTACGGTGTGAACTACTTTGCAATCCGGAATAAAAAGGGCACAGAGCTGCTGC
 TTGGAGTGGATGCCCTGGGGCTTACATTTATGACCCTGAGAACAGACTGACCCCAAGATCTCCTTCCC
 GTGGAATGAAATCCGAAACATCTCGTACAGTGACAAGGAGTTTACTATTAACCCTGGATAAGAAAATT
 GATGTCTTCAAGTTTAACTCCTCAAAGCTTCGTGTTAATAAGCTGATTCTCCAGCTATGTATCGGGAACC
 ATGATCTATTTATGAGGAGAAGGAAAGCCGATTCTTTGGAAGTTCAGCAGATGAAAGCCCAGGCCAGGGA
 GGAGAAGGCTAGAAAGCAGATGGAGCGGCAGCGCTCGCTCGAGAGAAGCAGATGAGGGAGGAGGCTGAA
 CGCAGGAGGATGAGTTGGAGAGGAGGCTGCTGCAGATGAAAGAAGAAGCAACAATGGCCAACGAAGCAC
 TGATGCGGTCTGAGGAGACAGTGACCTGTTGGCTGAAAAGGCCAGATCACCGAGGAGGAGGCAAAACT
 TCTGGCCAGAAAGGCCGAGAGGCTGAGCAGGAAATGCAGCGCATCAAGGCCACAGCGATTTCGCACGGAG
 GAGGAGAAGCGCCTGATGGAG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:

>RC401975 representing NM_000268
 Red=Cloning site Green=Tags(s)

MAGAIASRMSFSSLKRKQPKFTVIRIVTMDAEMEFNCEMKWKGDLDLVCRTLGLRETWFFGLQYTIKD
 TVAWLKMDKKVLDHDVSKEEPVTFHFLAKFYPENAEELVQEITQHLFFLQVKKQILDEKIYCPPEASVL
 LASAVQAKYGDYDPSVHKRGFLAQEELLPKRVINLYQMPPEMWEERITAWYAEHRGRARDEAEEMEYLKI
 AQDLEMYGVNYFAIRNKKGTLELLGVDALGLHIYDPENRLTPKISFPWNEIRNISYSKFEFTIKPLDKKI
 DVFKFNSKLRVNLILQLCIGNHDLFMRRRKADSLEVQMKQAAREEKARKQMERQRLAREKQMRREEAE
 RTRDELERLLQMKEEATMANEALMRSEETADLLAEKAQITEEEAKLLAQKAAEAEQEMQRIKATAIRTE
 EEKRLME

SGPTRTRRLE**QKLI**SEEDLAANDILDYKDDDDK**V**

Restriction Sites:

SgfI-MluI

Cloning Scheme:

OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq:

[NP_000259](#)

RefSeq Size:

1281 bp

RefSeq ORF:

1788 bp

Locus ID:

4771

Cytogenetics:

22q12.2

Domains:

B41, ERM

Protein Families:

Druggable Genome

MW:

47 kDa

Gene Summary:

This gene encodes a protein that is similar to some members of the ERM (ezrin, radixin, moesin) family of proteins that are thought to link cytoskeletal components with proteins in the cell membrane. This gene product has been shown to interact with cell-surface proteins, proteins involved in cytoskeletal dynamics and proteins involved in regulating ion transport. This gene is expressed at high levels during embryonic development; in adults, significant expression is found in Schwann cells, meningeal cells, lens and nerve. Mutations in this gene are associated with neurofibromatosis type II which is characterized by nervous system and skin tumors and ocular abnormalities. Two predominant isoforms and a number of minor isoforms are produced by alternatively spliced transcripts. [provided by RefSeq, Jul 2008]