

## Product datasheet for **RC402229**

### LKB1 (STK11) (NM\_000455) Human Mutant ORF Clone

#### Product data:

Product Type:	Mutant ORF Clones
Product Name:	LKB1 (STK11) (NM_000455) Human Mutant ORF Clone
Mutation Description:	Q220X
Affected Codon#:	220
Affected NT#:	658
Nucleotide Mutation:	STK11 Mutant (Q220X), Myc-DDK-tagged ORF clone of Homo sapiens serine/threonine kinase 11 (STK11) as transfection-ready DNA
Effect:	Peutz-Jeghers syndrome
Symbol:	STK11
Synonyms:	hLKB1; LKB1; PJS
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_000455
ORF Size:	657 bp
Restriction Sites:	Sgfi-MluI



[View online »](#)

**ORF Nucleotide Sequence:**

>RC402229 representing NM\_000455  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGAGGTGGTGGACCCGAGCAGCTGGGCATGTTACGGAGGGCGAGCTGATGTGGTGGTATGGACA  
 CGTTCATCCACCGCATCGACTCCACCGAGGTCACTACCAGCCGCGCCGAAGCGGGCCAAGCTCATCGG  
 CAAGTACCTGATGGGGACCTGCTGGGGGAAGGCTCTTACGGCAAGGTGAAGGAGGTGCTGGACTCGGAG  
 ACGCTGTGCAGGAGGGCCGCAAGATCCTCAAGAAGAAGAAGTTGCAAGGATCCCCAACGGGGAGGCCA  
 ACGTGAAGAAGGAAATCAACTACTGAGGAGGTTACGGCACAAAAATGTCATCCAGCTGGTGGATGTGT  
 ATACAACGAAGAGAAGCAGAAAATGTATATGGTATGGAGTACTGCGTGTGGCATGCAGGAAATGCTG  
 GACAGCGTGCCGGAAGCGTTCCAGTGTGCCAGGCCACGGGTACTTCTGTCAGCTGATTGACGGCC  
 TGGAGTACCTGCATAGCCAGGGCATTGTGCACAAGGACATCAAGCCGGGAACCTGCTGCTCACCACCGG  
 TGGACCCTCAAATCTCGACCTGGGCGTGGCCGAGGACTGCACCCGTTTCGGCGCGGACGACACTGC  
 CGGACCAGCCAGGGCTCCCCGGCTTTC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

**Protein Sequence:**


>RC402229 representing NM\_000455  
 Red=Cloning site Green=Tags(s)

MEVVDPPQLGMFTEGELMSVGMDFIHRIDSTEVYQPRRKRAKLIGKYLMDLLGEGSYGKVKVELDSE  
 TLCRRAYKILKKKLRIPNGEANKKEIQLLRRLRHKNVIQLVDVLYNEEKQKMYMMEYCVGMQEML  
 DSVPEKRFPVCQAHGYFCQLIDGLEYLHSQIVHKDIKPGNLLTTGGTLKISDLGVAEALHPFAADTCT  
 RTSQGSFAP

**SGP**TRRRLE**QKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**


\* The last codon before the Stop codon of the ORF

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>RefSeq:</b>	<a href="#">NP_000446</a>
<b>RefSeq Size:</b>	657 bp
<b>RefSeq ORF:</b>	1302 bp
<b>Locus ID:</b>	6794
<b>Cytogenetics:</b>	19p13.3
<b>Domains:</b>	pkinase, TyrKc, S_TKc
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	Adipocytokine signaling pathway, mTOR signaling pathway
<b>MW:</b>	24.1 kDa
<b>Gene Summary:</b>	<p>This gene, which encodes a member of the serine/threonine kinase family, regulates cell polarity and functions as a tumor suppressor. Mutations in this gene have been associated with Peutz-Jeghers syndrome, an autosomal dominant disorder characterized by the growth of polyps in the gastrointestinal tract, pigmented macules on the skin and mouth, and other neoplasms. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized. [provided by RefSeq, Jul 2008]</p>