

## Product datasheet for **RC402379**

### Tuberin (TSC2) (NM\_000548) Human Mutant ORF Clone

#### Product data:

|                           |                                                                                                                                           |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Mutant ORF Clones                                                                                                                         |
| Product Name:             | Tuberin (TSC2) (NM_000548) Human Mutant ORF Clone                                                                                         |
| Mutation Description:     | E700K                                                                                                                                     |
| Affected Codon#:          | 700                                                                                                                                       |
| Affected NT#:             | 2098                                                                                                                                      |
| Nucleotide Mutation:      | TSC2 Mutant (E700K), Myc-DDK-tagged ORF clone of Homo sapiens tuberous sclerosis 2 (TSC2), transcript variant 1 as transfection-ready DNA |
| Effect:                   | Tuberous sclerosis                                                                                                                        |
| Symbol:                   | Tuberin                                                                                                                                   |
| Synonyms:                 | LAM; PPP1R160; TSC4                                                                                                                       |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                      |
| Mammalian Cell Selection: | Neomycin                                                                                                                                  |
| Vector:                   | pCMV6-Entry (PS100001)                                                                                                                    |
| Tag:                      | Myc-DDK                                                                                                                                   |
| ACCN:                     | NM_000548                                                                                                                                 |
| ORF Size:                 | 5421 bp                                                                                                                                   |
| Restriction Sites:        | SgfI-XhoI                                                                                                                                 |
| ORF Nucleotide Sequence:  | >RC402379 representing NM_000548<br>Red=Cloning site Blue=ORF Green=Tags(s)                                                               |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGC**C

ATGGCCAAACCAACAAGCAAAGATTGAGGCTTGAAGGAGAAGTTTAAGATTCTGTTGGGACTGGGAACAC  
 CGAGGCCAAATCCCAGGTCTGCAGAGGGTAAACAGACGGAGTTTATCATCACCGCGGAAATACTGAGAGA  
 ACTGAGCATGGAATGTGGCCTCAACAATCGCATCCGGATGATAGGCGAGATTTGTGAAGTCGCAAAAC  
 AAGAAATTTGAAGAGCACGCAGTGGAAGCACTCTGGAAGGCGGTCGCGGATCTGTTGCAGCCGAGCGGC  
 CGCTGGAGGCCCGGCACGCGGTGCTGGCTCTGCTGAAGGCCATCGTGCAGGGGCAGGGCGAGCGTTTGGG  
 GGTCTCAGAGCCCTCTCTTTAAGGTCATCAAGGATTACCCTTCCAACGAAGACCTTCACGAAAGGCTG  
 GAGGTTTTCAAGGCCCTCACAGACAATGGGAGACACATCACCTACTTGGAGGAAGAGCTGGCTGACTTTG


[View online »](#)

TCCTGCAGTGGATGGATGTTGGCTTGCTCCTCGGAATTCCTTCTGGTGCTGGTGAACCTGGTCAAATTCAA  
TAGCTGTTACCTCGACGAGTACATCGCAAGGATGGTTCAGATGATCTGTCTGCTGTGCGTCCGGACCGCG  
TCCTCTGTGGACATAGAGGTCTCCCTGCAGGTGCTGGACGCCGTGGTCTGCTACAACCTGCCTGCCGGCTG  
AGAGCCTCCCGCTGTTTCATCGTTACCTCTGTGCGACCATCAACGTCAAGGAGCTCTGCGAGCCTTGCTG  
GAAGCTGATGCGGAACCTCCTTGGCACCACCTGGGCCACAGCGCCATCTACAACATGTGCCACCTCATG  
GAGGACAGAGCCTACATGGAGGACGCGCCCTGCTGAGAGGAGCCGTGTTTTTGGGCGATGGCTCTCT  
GGGAGGCCACCGGCTCTATTCTCTCAGGAACCTGCGGACATCTGTGTTGCCATCATTTTACCAGGCCAT  
GGCATGTCCGAACGAGGTGGTGTCTATGAGATCGTCTGTCCATCACCAGGCTCATCAAGAAGTATAGG  
AAGGAGCTCCAGGTGGTGGCGTGGGACATTCTGCTGAACATCATCGAACGGCTCCTTACGACGCTCCAGA  
CCTTGGACAGCCCGGAGCTCAGGACCATCGTCCATGACCTGTTGACCACGGTGGAGGAGCTGTGTGACCA  
GAACGAGTTCACGGGTCTCAGGAGAGATACTTGAACCTGGTGGAGAGATGTGCGGACCAGAGGCCTGAG  
TCCTCCCTCCTGAACCTGATCTCTATAGAGCGCAGTCCATCCACCCGGCCAAGGACGGCTGGATTGAGA  
ACCTGCAGGCGCTGATGGAGAGATTCTCAGGAGCGAGTCCGAGGCGCCGTGCGCATCAAGGTGCTGGA  
CGTGCTGTCTTTGTGCTGCTCATCAACAGGCAGTTCTATGAGGAGGAGCTGATTAACCTAGTGGTCATC  
TCGAGCTCTCCACATCCCGAGGATAAAGACCACAGGTCCGAAAGCTGGCCACCCAGTTGCTGGTGG  
ACCTGGCAGAGGGCTGCCACACACACCACTTCAACAGCCTGCTGGACATCATCGAGAAGGTGATGGCCCG  
CTCCCTCTCCCCACCCCGGAGCTGGAAGAAAGGGATGTGGCCGCATACTCGGCCCTCCTTGGAGGATGTG  
AAGACAGCCGTCTGGGGCTTCTGGTCACTCTTACAGCAAGCTGTACACCTGCCTGCAAGCCACGCCA  
CGCGTGTGATAGAGTGTGTCAGCCACATTCAGTCCACTACAAGCACAGCTACACCTGCCAATCGC  
GAGCAGCATCCGGCTGCAGGCCCTTCTGACTTCTGTTGCTGCTGCGGGCCGACTCACTGCACCCGCTGGGC  
CTGCCCCAACAAGGATGGAGTCGTGCGGTTACGCCCCCTACTGCGTCTGCGACTACATGGAGCCAGAGAG  
GCTCTGAGAAGAAGACCAGCGGCCCTTCTCTCTCCACAGGGCCTCCTGGCCCGCGCCTGCAGGCC  
CGCCGTGCGGCTGGGGTCCGTGCCCTACTCCCTGCTCTTCCGCGTCTGCTGCAGTGTGAAGCAGAAG  
TCTGACTGGAAGGTGCTGAAGCTGGTTCTGGGACGGCTGCCTGAGTCCCTGCGCTATAAAGTGCTCATCT  
TTACTTCCCCTTGCAGTGTGGACCAGCTGTGCTGCTCTGCTCTGCTCCATGCTTTACGGCCAAAGACACT  
GGAGCGGCTCCGAGGCGCCCCAGAAGGCTTCTCCAGAAGTGTGACCTGGCCGTGGTTCCAGTGCTG  
ACAGCATTAATCTCTTACCATAACTACCTGGACAAAACCAAACAGCGGAGATGGTCTACTGCCTGGAGC  
AGGGCCTCATCCACCGCTGTGCCAGCCAGTGCCTGCTGGCCTTGTCCATCTGCAGCGTGGAGATGCCTGA  
CATCATCATCAAGGCGCTGCCTGTTCTGGTGGTGAAGCTCACGCACATCTCAGCCACAGCCAGCATGGCC  
GTCCCACTGCTGGAGTTCCTGTCCACTCTGCCAGGCTGCCGACCTCTACAGGAACCTTGGCCGGGAGC  
AGTATGCCAGTGTGTTGCCATCTCCCTGCCGTACACCAACCCCTCCAAGTTTAATCAGTACATCGTGTG  
TCTGGCCCATCAGTCATAGCCATGTGGTTTCATCAGGTGCCGCTGCCCTTCCGGAAGGATTTTGTCCCT  
TTCATCACTAAGGGCCTGCGGTCCAATGTCTCTTGTCTTTGATGACACCCCGAGAAGGACAGCTTCA  
GGGCCCCGAGTACTAGTCTCAACGAGAGACCAAGAGTCTGAGGATAGCCAGACCCCCCAACAAGGCTT  
GAATAACTCTCCACCCGTGAAAGAATTCAGGAGAGCTCTGCAGCCGAGGCCCTTCCGGTGCCGCAGCATC  
AGTGTGTCTGAACATGTGGTCCGAGCAGGATACAGACGTCCCTCACCAGTGCCAGCTTGGGGTCTGCAG  
ATGAGAAGTCCGTGGCCAGGCTGACGATAGCCTGAAAAACCTCCACCTGGAGCTCACGGAAACCTGTCT  
GGACATGATGGCTCGATACGTCTTCTCCAATTACGGCTGTCCCGAAGAGGTCTCCTGTGGGCGAGTTC  
CTCTAGCGGGTGGCAGGACAAAACCTGGCTGGTTGGGAACAAGCTTGTCACTGTGACGACAAGCGTGG  
GAACCGGACCCGGTCTTACTAGGCCCTGGACTCGGGGGAGCTGCAAGTCCGGCCCGGAGTCGAGCTCCAG  
CCCCGGGTGCTGTGAGACAGACCAAGGAGGCGCCGCGCAAGCTGGAGTCCAGGCTGGGCAGCAGGTG  
TCCCGTGGGGCCCGGATCGGGTCCGTTCCATGTGCGGGGGCCATGGTCTTGCAGTTGGCGCCCTGGACG  
TGCCGGCTCCAGTTCTGGGAGTGCCACTTCTCCAGGACCAGGACTGCACAGCCGCGAAACCTGA  
GAAGGCCTCAGCTGGACCCGGGTTCTGTGACGAGAGAAGACGAACCTGGCGGCCTATGTGCCCTGCTG  
ACCCAGGGCTGGGCGGAGATCCTGGTCCGAGGCCCCACAGGGAACACCAGCTGGCTGATGAGCCTGGAGA  
ACCCGCTCAGCCCTTCTCTCGGACATCAACAACATGCCCTGCAGGAGCTGTCTAACGCCCTCATGGC  
GGCTGAGCGCTTCAAGGAGCACCGGGACACAGCCCTGTACAAGTCACTGTGGTGGCGGACGCCAGCAGC  
GCCAAACCCCTCCTCTGCCTCGCTCCAACACAGTGGCTCTTTCTCCTCCTGTACCAGTCCAGCTGCC  
AAGGACAGCTGCACAGGAGCGTTTCTGGGCGAGCTCCGCCGTGGTCACTGGAGGAGGGAAGTCCGGGCGA  
GGTTCCTGTGCTGGTGGAGCCCCAGGGTTGGAGGACGTTGAGGCAGCGCTAGGCATGGACAGGCGCAGC  
GATGCCTACAGCAGGTGCTCCTCAGTCTCCAGCCAGGAGGAGAAGTGCCTCCACGCGGAGGAGCTGGTTG  
GCAGGGGCATCCCCATCAGCGAGTGTCTCCTCGAGGGTGGCCGGCCCTGTGTGACCTCTCCTTCCA

GCCCTCGCAGCCCCTGAGCAAGTCCAGCTCCTCTCCCGAGCTGCAGACTCTGCAGGACATCCTCGGGGAC  
 CCTGGGGACAAGGCCGACGTGGGCCGGCTGAGCCCTGAGGTTAAGGCCCGGTACAGTCAGGGACCCTGG  
 ACGGGGAAAGTGTGCTGGTGGCCTCGGGCGAAGACAGTCGGGGCCAGCCGAGGGTCCCTTGCCCTC  
 CAGCTCCCCCGCTCGCCAGTGGCCTCCGGCCCCGAGGTTACACCATCTCCGACTCGGCCCCATCACGC  
 AGGGGCAAGAGAGTAGAGAGGGACGCCCTTAAGAGCAGAGCCACAGCCTCCAATGCAGAGAAAGTGCCAG  
 GCATCAACCCAGTTTCGTGTTCTGCAGCTCTACCATTCCCCCTTCTTTGGCGACGAGTCAAACAAGCC  
 AATCCTGCTGCCAATGAGTCACAGTCCTTTGAGCGGTCGGTGCAGCTCCTCGACCAGATCCCATCATAC  
 GACACCCACAAGATCGCCGTCTGTATGTTGGAGAAGGCCAGAGCAACAGCGAGCTCGCCATCCTGTCCA  
 ATGAGCATGGCTCCTACAGGTACACGGAGTTCTTGACGGGCTGGGCCGGCTCATCGAGCTGAAGGACTG  
 CCAGCCGGACAAGGTGTACCTGGGAGGCCTGGACGTGTGTGGTGAAGACGGCCAGTTCACCTACTGCTGG  
 CACGATGACATCATGAAGCCGTCTTCCACATCGCCACCCTGATGCCACCAAGGACGTGGACAAGCACC  
 GCTGCGACAAGAAGCGCCACCTGGGCAACGACTTTGTGTCCATTGTCTACAATGACTCCGGTGAAGACTT  
 CAAGCTTGGCACCATCAAGGGCCAGTTCAACTTTGTCCACGTGATCGTACCCCGCTGGACTACGAGTGC  
 AACCTGGTGTCCCTGCAGTGCAGGAAAGACATGGAGGGCCTTGTGGACACCAGCGTGGCCAAGATCGTGT  
 CTGACCGCAACCTGCCCTTCGTGGCCCGCCAGATGGCCCTGCACGCAAAATATGGCCTCACAGGTGCATCA  
 TAGCCGCTCCAACCCACCGATATCTACCCCTCCAAGTGGATTGCCCGGCTCCGCCACATCAAGCGGCTC  
 CGCCAGCGGATCTGCGAGGAAGCCGCCTACTCCAACCCAGCCTACCTCTGGTGCACCCTCCGTCCCAT  
 GCAAAGCCCCTGCACAGACTCCAGCCGAGCCACACCTGGCTATGAGGTGGGCCAGCGGAAGCGCCTCAT  
 CTCCTCGGTGGAGGACTTCACCGAGTTTGTG

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

**Protein Sequence:**

>RC402379 representing NM\_000548

Red=Cloning site Green=Tags(s)


MAKPTSKDSGLKEKFKILLGLGTRPNRPSAEGKQTEFIITAEILRELSMECGLNNRIRMIGQICEVAK  
KKFEEHAVEALWKAVADLLQPERPLEARHAVLALLKAIVQGQGERLGVLRALFFVKIKDYPSEDNLHERL  
EVFKALTDNGRHITYLEELADFVLQWMDVGLSSEFLLVLVNLVKFNSCYLDEYIARMVMICLLCVRTA  
SSVDIEVSLQVLDAVVCYNCLPAESLPLFIVTLCRTINVKELCEPCWKLNRNLGTHLGHSATYNMCHLM  
EDRAYMEDAPLLRGAVFFVGMALWGHAHRLYSLRNSPTSVLPSFYQAMACPNEVVSYEIVLSITRLIKKY  
KELQVVAWDILLNIIERLLQQLQTLDSPELRTIVHDLLTTVEELCDQNEFHGSQERYFELVERCADQRP  
SSLLNLSISYRAQSIHPAKDGIWQNLQALMERFFRSESRGAVRIKVLDVLSFVLLINRQFYEEELINSVVI  
SQLSHIPEDKDHQVRKLATQLLDLAEGCHTHFNSLLDIEKVMARSLPPPELEERDVAAYSASLEDV  
KTAVLGLLVILQTKLYTLPAASHATRVYEMLVSHIQLHYKHSYTLPIASSIRLQAFDFLLLRADSLHRLG  
LPNKDGVVRFSPYCVCDYMERPERGSEKKTSGPLSPPTGPPGAPAGPAVRLGSVPYSLFRVLLCQLKQK  
SDWKVLKLVGLRPLESRLYKVLIFTPSCPSVDQLCSALCSLSPGKTPLEIRGAPEGFSRTDLHLAAVVPVL  
TALISYHNYLDKTKQREMYCYCLEQGLIHRCASQCVAALSCMSVEMPDIIRKALPVLVVKLTHSATSMA  
VPLLEFLSTLARPLHYRNFAAEQYASVFAISLPTYNPSKFNQYIVCLAHHVIAWFIIRCRLPFRKDFVP  
FITKGLRSNVLLSFDTPKDSFRARSTSLNERPKSLRIARPPKQGLNNSPPVKEFKESSAAEAFRCRSI  
SVSEHVRSRIQTSLSASLGSADENVAQADDSLKNLHLELTETCLDMMARYVFSNFTAVPKRSPVGEF  
LLAGGRKTWLGNKLVTVTTSVGTGTRSLLGLDSGELQSGPESSSSPGVHVRQTKEAPAKLESQAGQV  
SRGARDRVSRMSGGHGLRVGALDVPASQFLGSATSPGPRTAPAAKPEKASAGTRVPVQEKTNLAAYVPL  
TQGWAEILVRRPTGNTSWLMSLENPLSPFSSDINNMPQLQELSNALMAAERFKEHRDTALYKLSVPAAST  
AKPPPLPRSNTVASFSSLYQSSCQQLHRSVSWADSAVVMEEGSGPEVPVLEPPGLEDVEAALGMDRRT  
DAYSRSSSVSSQEEKSLHAEELVGRGIPIERVVSSEGGRPVDSLQFQSPQLSKSSSSPELQTLQDILGD  
PGDKADVGRLSPEVKARSQSGTLDGESAASGSGEDSRGQPEGLPSSSPRSPSGLRPRGYTISDSAPSR  
RGKRVERDALKSRATASNAEKVPGINPSFVFLQLYHSPFFGDESINKPILLPNESQSFSERSVQLLDQIPSY  
DTHKIAVLVYGEQGSNSLAILENEHGSYRYTEFLTGLGRLIELKDCQPKVYLLGGLDVCGEDGQFTYCW  
HDDIMQAVFHIAELMPTKDVAKHKRCKKRHLGNDFVSIYVNSDGEDFKLGTIKGFNFVHVIVTPLDYEC  
NLVSLQCRKDMEGLVDTVSXKIVSDRNLPFVARQMALHANMASQVHHSRNSPTDIYPSKWIARLRIKRL  
RORICEFAAYSNSPLPLVPPSHSKAPAOPTPAEPTPGYEVGORKRLISSVDEFTFV

SGPTRRRL**E**QKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

Sgfl-XhoI

### Cloning Scheme:



\* The last codon before the Stop codon of the ORF

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>   | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| <b>OTI Annotation:</b>   | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| <b>Components:</b>       | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                               |
| <b>Note:</b>             | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq:</b>           | <a href="#">NP_000539</a>                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq Size:</b>      | 5421 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>RefSeq ORF:</b>       | 5424 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>         | 7249                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Cytogenetics:</b>     | 16p13.3                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Domains:</b>          | Rap_GAP, Tuberin                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Families:</b> | Druggable Genome                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Pathways:</b> | Insulin signaling pathway, mTOR signaling pathway, p53 signaling pathway                                                                                                                                                                                                                                                                                                                                                                           |
| <b>MW:</b>               | 198.8 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Gene Summary:</b>     | Mutations in this gene lead to tuberous sclerosis complex. Its gene product is believed to be a tumor suppressor and is able to stimulate specific GTPases. The protein associates with hamartin in a cytosolic complex, possibly acting as a chaperone for hamartin. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]                                                     |