

## Product datasheet for **RC402561**

### **BubR1 (BUB1B) (NM\_001211) Human Mutant ORF Clone**

#### **Product data:**

Product Type:	Mutant ORF Clones
Product Name:	BubR1 (BUB1B) (NM_001211) Human Mutant ORF Clone
Mutation Description:	R194X
Affected Codon#:	194
Affected NT#:	580
Nucleotide Mutation:	BUB1B Mutant (R194X), Myc-DDK-tagged ORF clone of Homo sapiens budding uninhibited by benzimidazoles 1 homolog beta (yeast) (BUB1B) as transfection-ready DNA
Effect:	Mosi vrieded neuploidy
Symbol:	BUB1B
Synonyms:	Bub1A; BUB1beta; BUBR1; hBUBR1; MAD3L; MVA1; SSK1
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_001211
ORF Size:	579 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)



<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
<b>Components:</b>	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
<b>RefSeq:</b>	<p><a href="#">NP_001202</a></p>
<b>RefSeq Size:</b>	<p>579 bp</p>
<b>RefSeq ORF:</b>	<p>3153 bp</p>
<b>Locus ID:</b>	<p>701</p>
<b>Cytogenetics:</b>	<p>15q15.1</p>
<b>Protein Families:</b>	<p>Druggable Genome, Protein Kinase, Stem cell - Pluripotency</p>
<b>Protein Pathways:</b>	<p>Cell cycle</p>
<b>MW:</b>	<p>21.2 kDa</p>
<b>Gene Summary:</b>	<p>This gene encodes a kinase involved in spindle checkpoint function. The protein has been localized to the kinetochore and plays a role in the inhibition of the anaphase-promoting complex/cyclosome (APC/C), delaying the onset of anaphase and ensuring proper chromosome segregation. Impaired spindle checkpoint function has been found in many forms of cancer. [provided by RefSeq, Jul 2008]</p>