

Product datasheet for **RC402622**

c-Myc (MYC) (NM_002467) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	c-Myc (MYC) (NM_002467) Human Mutant ORF Clone
Mutation Description:	N26S
Affected Codon#:	26
Affected NT#:	77
Nucleotide Mutation:	MYC Mutant (N26S), Myc-DDK-tagged ORF clone of Homo sapiens v-myc myelocytomatosis viral oncogene homolog (avian) (MYC) as transfection-ready DNA
Effect:	Breast cancer, increased risk, association with
Symbol:	c-Myc
Synonyms:	bHLHe39; c-Myc; MRTL; MYCC
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_002467
ORF Size:	1362 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC402622 representing NM_002467
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

CTGGATTTTTTCGGGTAGTGGAAAACCAGCAGCCTCCCGCAGCATGCCCTCAACGTTAGCTTACCA
 ACAGGAGCTATGACCTCGACTACGACTCGGTGCAGCCGATTTTCTACTGCGACGAGGAGAGAATTCTA
 CCAGCAGCAGCAGCAGAGCGAGCTGCAGCCCCGGCGCCAGCGAGGATATCTGGAAGAAATTCGAGCTG
 CTGCCACCCCGCCCTGTCCCCTAGCCGCCGCTCCGGGCTCTGCTCGCCCTCCTACGTTGCGGTACAC
 CTTTCTCCCTCGGGGAGACAACGACGGCGGTGGCGGGAGCTTCTCCACGGCCGACCAGCTGGAGATGGT
 GACCGAGCTGCTGGGAGGAGACATGGTGAACCAGAGTTTCATCTGCGACCCGGACGACGAGACCTTCATC
 AAAACATCATCATCCAGGACTGTATGTGGAGCGGCTTCTCGGCCCGCCAAGCTCGTCTCAGAGAAGC
 TGGCTCTACCAGGCTGCGCGCAAAGACAGCGGCAGCCGAACCCCGCCGCGCCACAGCGTCTGCTC
 CACCTCCAGCTTGTACCTGCAGGATCTGAGCGCCGCCCTCAGAGTGCATCGACCCCTCGGTGGTCTTC
 CCTACCTCTCAACGACAGCAGCTCGCCAAGTCTGCGCCTCGCAAGACTCCAGCGCCTTCTCTCCGT
 CCTCGGATTCTGCTCTCTCTCGACGGAGTCTCCCGCAGGGCAGCCCGAGCCCTGGTGTCTCATGA
 GGAGACACCGCCACCACCAGCAGCGACTCTGAGGAGGAACAAGAAGATGAGGAAGAAATCGATGTTGTT
 TCTGTGAAAAGAGGCAGGCTCTGGCAAAAGGTGAGAGTCTGGATCACCTTCTGCTGGAGGCCACAGCA
 AACCTCCTCACAGCCACTGGTCTCAAGAGGTGCCACGTCTCCACACATCAGCACAACACGACGCGCC
 TCCCTCCACTCGGAAGGACTATCCTGCTGCCAAGAGGGTCAAGTTGGACAGTGTGAGAGTCTGAGACAG
 ATCAGCAACAACCGAAAATGCACCAGCCCCAGGCTCCTCGGACACCGAGGAGAATGTCAAGAGGCGAACAC
 ACAACGCTTGGAGCGCCAGAGGAGGAACGAGCTAAAACGGAGCTTTTTTGGCCCTGCGTGACCACTCC
 GGAGTTGAAAACAATGAAAAGGCCCAAGGTAGTTATCCTTAAAAAGCCACAGCATACATCCTGTCC
 GTCCAAGCAGAGGAGCAAAAGCTCATTCTGAAGAGGACTTGTTCGCGAAAACGACGAGAACAGTTGAAAC
 ACAAACTGAACAGCTACGGAACCTTGTGCG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:

>RC402622 representing NM_002467
 Red=Cloning site Green=Tags(s)


LDFFRVVENQPPATMPLNVSFTNRSYDLDYDSVQPYFYCDEEENFYQQQQSELQPPAPSEDIWKKFEL
 LPTPPLSPRRSGLCSPSYVAVTPFLRGDNDGGGSFSTADQLEMVTELLGGDMVNQSFICDPDDETFI
 KNII IQDCMWSGF SAAAKLVSEKLASYQAARKDSGSPNPARGHSVCSTSSLYLQDL SAAASECIDPSVVF
 PYPLNDSSSPKSCASQDSSAFSPSSDLLSSTESSPQGSPEPLVLHEETPPTTSSDSEEEQEDEEEIDVV
 SVEKRQAPGKRSESGSPSAGGHSKPPHSPLVLKRCHVSTHQHNYAAPSTRKDYPAAKRVKLDsvrvlrq
 ISNNRKCTSPRSSDTEENVKRRTHNVLERQRRNELKRSFFALRDQIPELENNEKAPKVVILKkatayils
 VQAEQKLI SEEDLLRKRREQLKHKLEQLRNSCA

SGPTRRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NP_002458](#)

RefSeq Size:

1362 bp

RefSeq ORF:

1365 bp

Locus ID:

4609

Cytogenetics:

8q24.21

Domains:

HLH, Myc_N_term, Myc-LZ

Protein Families:	Druggable Genome, Embryonic stem cells, Induced pluripotent stem cells, Stem cell - Pluripotency, Stem cell relevant signaling - JAK/STAT signaling pathway, Stem cell relevant signaling - TGFb/BMP signaling pathway, Stem cell relevant signaling - Wnt Signaling pathway, Transcription Factors
Protein Pathways:	Acute myeloid leukemia, Bladder cancer, Cell cycle, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Jak-STAT signaling pathway, MAPK signaling pathway, Pathways in cancer, Small cell lung cancer, TGF-beta signaling pathway, Thyroid cancer, Wnt signaling pathway
MW:	49.9 kDa
Gene Summary:	This gene is a proto-oncogene and encodes a nuclear phosphoprotein that plays a role in cell cycle progression, apoptosis and cellular transformation. The encoded protein forms a heterodimer with the related transcription factor MAX. This complex binds to the E box DNA consensus sequence and regulates the transcription of specific target genes. Amplification of this gene is frequently observed in numerous human cancers. Translocations involving this gene are associated with Burkitt lymphoma and multiple myeloma in human patients. There is evidence to show that translation initiates both from an upstream, in-frame non-AUG (CUG) and a downstream AUG start site, resulting in the production of two isoforms with distinct N-termini. [provided by RefSeq, Aug 2017]