

## Product datasheet for **RC403210**

### AMHR2 (NM\_020547) Human Mutant ORF Clone

#### Product data:

Product Type:	Mutant ORF Clones
Product Name:	AMHR2 (NM_020547) Human Mutant ORF Clone
Mutation Description:	R80X
Affected Codon#:	80
Affected NT#:	238
Nucleotide Mutation:	AMHR2 Mutant (R80X), Myc-DDK-tagged ORF clone of Homo sapiens anti-Mullerian hormone receptor, type II (AMHR2), transcript variant 1 as transfection-ready DNA
Effect:	Persisen Mullerin dus syndrome
Symbol:	AMHR2
Synonyms:	AMHR; MISR2; MISRII; MRII
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_020547
ORF Size:	237 bp
Restriction Sites:	Sgfl-Mlul
ORF Nucleotide Sequence:	>RC403210 representing NM_020547 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCTAGGGTCTTTGGGCTTTGGGCATTACTTCCCACAGCTGTGGAAGCACCCCAAACAGGCGAACCT  
GTGTGTTCTTTGAGGCCCTGGAGTGCGGGGAAGCACAAAGACTGGGAGAGCTGCTAGATACAGGCAC  
AGAGCTCCCCAGAGCTATCCGCTGCCTTACAGCCGCTGCTGTTGGGATCTGGAACCTGACCCAAGAC  
CGGGCACAGGTGGAATGCAAGGATGC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGA TAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC403210 representing NM\_020547  
 Red=Cloning site Green=Tags(s)

MLGSLGLWALLPTAVEAPPNRRTCVFFEAPGVRGSKTLGELLDGTGELPRAIRCLYSRCCFGIWNLTQD  
 RAQVEMQGC

SGP TRTRRL EQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**RefSeq:** [NP\\_065434](#)

**RefSeq Size:** 237 bp

**RefSeq ORF:** 1722 bp

**Locus ID:** 269

**Cytogenetics:** 12q13.13

**Protein Families:** Druggable Genome, Protein Kinase, Transmembrane

**Protein Pathways:** Cytokine-cytokine receptor interaction, TGF-beta signaling pathway

**MW:** 8.7 kDa

**Gene Summary:**

This gene encodes the receptor for the anti-Mullerian hormone (AMH) which, in addition to testosterone, results in male sex differentiation. AMH and testosterone are produced in the testes by different cells and have different effects. Testosterone promotes the development of male genitalia while the binding of AMH to the encoded receptor prevents the development of the mullerian ducts into uterus and Fallopian tubes. Mutations in this gene are associated with persistent Mullerian duct syndrome type II. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Sep 2009]