

Product datasheet for **RC403445**

Wilms Tumor Protein (WT1) (NM_024426) Human Mutant ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Mutant ORF Clones |
| Product Name: | Wilms Tumor Protein (WT1) (NM_024426) Human Mutant ORF Clone |
| Mutation Description: | G321A |
| Affected Codon#: | 321 |
| Affected NT#: | 962 |
| Nucleotide Mutation: | WT1 Mutant (G321A), Myc-DDK-tagged ORF clone of Homo sapiens Wilms tumor 1 (WT1), transcript variant D as transfection-ready DNA |
| Effect: | Wilms tumour |
| Symbol: | Wilms Tumor Protein |
| Synonyms: | AWT1; GUD; NPHS4; WAGR; WIT-2; WT33 |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| Tag: | Myc-DDK |
| ACCN: | NM_024426 |
| ORF Size: | 1551 bp |
| Restriction Sites: | Sgfl-Mlul |



[View online »](#)

ORF Nucleotide Sequence:

>RC403445 representing NM_024426
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCAGGACCCGGCTTCCACGTGTGTCCCGGAGCCGGCTCTCAGCACACGCTCCGCTCCGGGCCTGGT
 GCCTACAGCAGCCAGAGCAGCAGGGAGTCCGGGACCCGGCGGCATCTGGGCCAAGTTAGCGCCGCCGA
 GGCCAGCGCTGAACGTCTCCAGGGCCGAGGAGCCGCGGGCGTCCGGGTCTGAGCCGACGAAATGGGC
 TCCGACGTGCGGGACCTGAACGCGCTGCTGCCGCCGTCCCTCCCTGGGTGGCGGCGCGGCTGTGCC
 TGCCTGTGAGCGGCGCGCAGTGGCGCCGGTGTGGACTTTGCGCCCCGGGCGCTTCGGCTTACGG
 GTCGTTGGGCGGCCCGCGCCGACCCGGCTCCGCCGCCACCCCGCGCGCGCCCTACTCCTTCATC
 AAACAGGAGCCGAGCTGGGGCGCGGAGCCGCACGAGGAGCAGTGCCTGAGCGCCTTACTGTCCACT
 TTTCCGGCCAGTTCACTGGCACAGCCGAGCCTGTGCTACGGGCCCTTCGGTCTCTCCGCCAGCCA
 GCGTCATCCGGCCAGGCCAGGATGTTTCTAACGCGCCCTACCTGCCAGCTGCCTCGAGAGCCAGCCC
 GCTATTCGCAATCAGGGTTACAGCACGGTCACCTTCGACGGGACGCCAGCTACGGTACACGCCCTCGC
 ACCATGCGGCGCAGTTCCCAACCACTATTCAAGCATGAGGATCCCATGGGCCAGCAGGGCTCGTGGG
 TGAGCAGCAGTACTCGGTGCCGCCCGGTCTATGGTGCCACACCCCAACCGACAGCTGCACCGGCAGC
 CAGGCTTTGCTGCTGAGGACGCCCTACAGCAGTGACAATTTATACCAATGACATCCCAGCTTGAATGCA
 TGACCTGGAATCAGATGAACTTAGGAGCCACCTTAAAGGGAGTTGCTGCTGCGAGCTCCAGCTCAGTGAA
 ATGGACAGAAGGGCAGAGCAACCACAGCACAGGATACGAGAGCGATAACCACACAACGCCATCCTCTGC
 GGAGCCCAATACAGAATACACACGCACGGTGTCTTACAGGGCATTACAGGATGTGCGACGTGTGCCGGAG
 TAGCCCCGACTTTGTACGGTCCGATCTGAGACCAGTGAAGAACGCCCTTTCATGTGCTTACCCAGG
 CTGCAATAAGAGATATTTTAAGCTGTCCCACTTACAGATGCACAGCAGGAAGCACACTGGTGAGAAACCA
 TACCAGTGTGACTTCAAGGACTGTGAACGAAGTTTTCTCGTTACAGCAGCTCAAAAGACACCAAGGA
 GACATACAGGTGTGAAACCATTCCAGTGTAAAACCTTGTGAGCGAAAGTTCTCCCGGTCCGACCACCTGAA
 GACCCACACCAGGACTCATAAGGTAACAAGTGAAGGCCCTTACAGTGTGCGTGGCCAAAGTTGTGAG
 AAAAGTTTCCCGGTGAGTGAATAGTCCGCCATCACAAATGCATCAGAGAAACATGACCAAACTCC
 AGCTGGCGCTT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:


>RC403445 representing NM_024426
 Red=Cloning site Green=Tags(s)

MQDPASTCVPEPASQHTLRSGPGCLQQPEQQVDRDPGGIWAKLGAAEASAERLQRRSRGASGSEPPQMG
 SDVDRDLNALLPAVPSLGGGGCALPVSGAAQWAPVLDFAAPPASAYGSLGGPAPPAPPPPPPPPHSFI
 KQEPSWGAEPHEEQCLSAFTVHFSGQFTGTAGACRYGPFPPPPSQASSGQARMFPNAPYLPSCLSEQP
 AIRNQYSTVTFDGTSPSYGHTPSHAAQFPNHSFKHEDPMGQQGSLGEQQYSVPPPVYGCHTPTDSTGS
 QALLLRTPYSSDNLQMTSQLCMTWNQMLGATLKGVAASSSSVKWTEGQSNHSTGYESDNHTPTILC
 GAQYRIHTHGVRGIQDVRVPGVAPTLVRSASETSEKRPFMCAYPGCNKRYFKLSHLQMSRKHTGEKP
 YQCDFKDCERRFSRSDQLKRHQRRTGVKPFQCKTCQRKFSRSDHLKTHTRHTGKTSEKPFSCRWPSCQ
 KKFARSDELVRHHNMHRNMTKLQAL

SGPTRRRLEQKLI**SEEDLAANDILDYKDDDDKV**

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NP_077744](#)

RefSeq Size:

1551 bp

RefSeq ORF:

1569 bp

Locus ID:

7490

Cytogenetics:

11p13

Domains:

WT1, zf-C2H2

Protein Families: Druggable Genome, Transcription Factors

MW: 56.9 kDa

Gene Summary: This gene encodes a transcription factor that contains four zinc-finger motifs at the C-terminus and a proline/glutamine-rich DNA-binding domain at the N-terminus. It has an essential role in the normal development of the urogenital system, and it is mutated in a small subset of patients with Wilms tumor. This gene exhibits complex tissue-specific and polymorphic imprinting pattern, with biallelic, and monoallelic expression from the maternal and paternal alleles in different tissues. Multiple transcript variants have been described. In several variants, there is evidence for the use of a non-AUG (CUG) translation initiation codon upstream of, and in-frame with the first AUG. Authors of PMID:7926762 also provide evidence that WT1 mRNA undergoes RNA editing in human and rat, and that this process is tissue-restricted and developmentally regulated. [provided by RefSeq, Mar 2015]