

Product datasheet for **RC403526**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	Y77X
Affected Codon#:	77
Affected NT#:	231
Nucleotide Mutation:	MEN1 Mutant (Y77X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	MEN1
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	228 bp
Restriction Sites:	SgfI-MluI
ORF Nucleotide Sequence:	>RC403526 representing NM_130799 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCGCCTGTTTG
CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGCTCCTTTCCCTGGTGGCTGGGCTTCGTGGAGCATTG
TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTCCAGCCAGCCCGCCCGCCGAC
CCGCTGGCGGCTCACC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGA TAAGTTTAA



[View online »](#)


Protein Sequence: >RC403526 representing NM_130799
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDDVVRLFAAELGREEPDLVLLSLVLGFVEHFLAVNRVIPTNPVPELTFQPSPAPD
 PPGGLT

SGP TRTRRL EQKL I SEEDLA AND I LDYK DDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq: [NP_570711](#)

RefSeq Size: 228 bp

RefSeq ORF: 1833 bp

Locus ID: 4221

Cytogenetics: 11q13.1

Domains: Menin

Protein Families: Druggable Genome, Transcription Factors

MW: 8.4 kDa

Gene Summary:

This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]