

## Product datasheet for **RC403577**

### Menin (MEN1) (NM\_130799) Human Mutant ORF Clone

#### Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	W183X
Affected Codon#:	183
Affected NT#:	548
Nucleotide Mutation:	MEN1 Mutant (W183X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	Menin
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	546 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

**ORF Nucleotide Sequence:**

>RC403577 representing NM\_130799  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCCTGTTTGC  
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGCTCCTTTCCCTGGTGTGGCTTCGTGGAGCATTT  
 TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCAGCCCGCCCCGAC  
 CCGCTGGCGCCTCACCTACTTTCCCGTGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGCTTCA  
 CCGCCAGATCCGAGGCGCCGTGCACCTGTCCCTCTATCCTCGAGAAGGGGGTGTCTCCAGCCGTGAGCT  
 GGTGAAGAAGGTCTCCGATGTATGGAACAGCCTCAGCCGCTCTACTTCAAGGATCGGGCCACATC  
 CAGTCCCTTTCAGTTCATCACAGGCACAAATTGGACAGCTCCGGTGTGCCCTTGTGTGGTTGGGG  
 CCTGCCAGGCCCTGGGTCTCCGGATGTCCACCTGCCTGTCTGAGGATCATGCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGA TAAGGTTAA

**Protein Sequence:**


>RC403577 representing NM\_130799  
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDDVVRFLFAAELGREEDLVLLSLVLGFVEHFLAVNRVIPTNVPCLTFQSPAPD  
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDLSLYPREGGVSSRELVKKVSVDVIWNSLSRSYFKDRAHI  
 QSLFSFITGKLDSSGVAFVVGACQALGLRDVHLALSEDHA

SGP**TRTRRLEQKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**

**OTI Disclaimer:**

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<a href="#">NP_570711</a>
<b>RefSeq Size:</b>	546 bp
<b>RefSeq ORF:</b>	1833 bp
<b>Locus ID:</b>	4221
<b>Cytogenetics:</b>	11q13.1
<b>Domains:</b>	Menin
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>MW:</b>	20 kDa
<b>Gene Summary:</b>	This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]