

## Product datasheet for **RC403580**

### Menin (MEN1) (NM\_130799) Human Mutant ORF Clone

#### Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	E191X
Affected Codon#:	191
Affected NT#:	571
Nucleotide Mutation:	MEN1 Mutant (E191X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	Menin
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	570 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

**ORF Nucleotide Sequence:**

>RC403580 representing NM\_130799  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCGCCCTGTTTG  
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGCTCCTTTCCCTGGTCTGGGCTTCGTGGAGCATTT  
 TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCAGCCCCGCCCGCCGAC  
 CCGCTGGCGGCCTCACCTACTTCCCGTGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGCTTCA  
 CCGCCAGATCCGAGGCGCCGTGCACCTGTCCCTCTATCCTCGAGAAGGGGGTCTCCAGCCGTGAGCT  
 GGTGAAGAAGTCTCCGATGTATGGAACAGCCTCAGCCGCTCCTACTTCAAGGATCGGGCCACATC  
 CAGTCCCTTTCAGTTCATCACAGGCACAAATTGGACAGCTCCGGTGTGCCCTTGTGTGGTGGGG  
 CCTGCCAGGCCCTGGTCTCCGGATGTCCACCTGCCTGTCTGAGGATCATGCCTGGGTAGTGTGGTGG  
 GCCAATGGG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

**Protein Sequence:**


>RC403580 representing NM\_130799  
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDVVRFLFAAELGREEDLVLLSLVLFVHFVFLAVNRVIPTNVPDLTFQSPAPD  
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDSLYPEGGVSSRELVKVSDVIWNSLSRSYFKDRAHI  
 QLSLFSFITGKLDSSGVAVFVGGACQALGLRDVHLALSEDHAWVFGPNG

**SGP**TRRRLE**QKLI**SEEDLAANDILDYKDDDDK**V**

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**


\* The last codon before the Stop codon of the ORF

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<a href="#">NP_570711</a>
<b>RefSeq Size:</b>	570 bp
<b>RefSeq ORF:</b>	1833 bp
<b>Locus ID:</b>	4221
<b>Cytogenetics:</b>	11q13.1
<b>Domains:</b>	Menin
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>MW:</b>	20.9 kDa
<b>Gene Summary:</b>	This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]