

Product datasheet for **RC403612**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	Q258X
Affected Codon#:	258
Affected NT#:	772
Nucleotide Mutation:	MEN1 Mutant (Q258X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	Menin
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	771 bp
Restriction Sites:	Sgfi-MluI



[View online »](#)

ORF Nucleotide Sequence:

>RC403612 representing NM_130799
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCGCCCTGTTTG
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTCTCCTTTCCCTGGTCTGGCTTCGTGGAGCATTT
 TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCAGCCCCGCCCCGAC
 CCGCTGGCGGCCTCACCTACTTTCCCGTGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGCTTCA
 CCGCCAGATCCGAGGCGCCGTGACCTGTCCCTCTATCCTCGAGAAGGGGGTGTCTCCAGCCGTGAGCT
 GGTGAAGAAGGTCTCCGATGTCATATGGAACAGCCTCAGCCGCTCCTACTTCAAGGATCGGGCCACATC
 CAGTCCCTTTCAGTTCATCACAGGCACAAATTGGACAGCTCCGGTGTGCCCTTTGCTGTGGTGGGG
 CCTGCCAGGCCCTGGGTCTCCGGATGTCCACCTGCCCTGTCTGAGGATCATGCCTGGGTAGTGTTTGG
 GCCCAATGGGGAGCAGACAGCTGAGGTACCTGGCAGCGCAAGGGCAACGAGGACCGAGGGCCAGACA
 GTCAATGCCGGTGTGGCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCCTGTGACCGCAAGA
 TGGAGGTGGCGTTCATGGTGTGTGCCATCAACCCTTCCATTGACCTGCACACCGACTCGCTGGAGCTTCT
 G

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:


>RC403612 representing NM_130799
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDDVVRLF AAE L GREPDLVLLSLVLGFVEHFLAVNRVIPTNPPELTFQPSPAPD
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDLSLYPREGGVSSREL VKKVSDVIWNSLSRSYFKDRAHI
 QSLFSFITGKLDSSGVAFVVGACQALGLRDVHLAL SEDHAWVVF GPNGEQTA EVTWHGKGNEDRRGQT
 VNAGVAERSWL YLKGSYMRCDRKMEVAFMVCAINPSIDLHTDSLELL

SGPTRRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NP_570711
RefSeq Size:	771 bp
RefSeq ORF:	1833 bp
Locus ID:	4221
Cytogenetics:	11q13.1
Domains:	Menin
Protein Families:	Druggable Genome, Transcription Factors
MW:	28.3 kDa
Gene Summary:	This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]