

Product datasheet for **RC403627**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	Y276X
Affected Codon#:	276
Affected NT#:	828
Nucleotide Mutation:	MEN1 Mutant (Y276X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	Menin
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	825 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC403627 representing NM_130799
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCGCCCTGTTTG
 CTGCCGAGCTGGGCCGAGAGGCCGGACCTGGTGCTCCTTTCCCTGGTCTGGCTTCGTGGAGCATTT
 TCTGGCTGTCAACCGCGTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCCAGCCCCGCCCCGAC
 CCGCTGGCGGCCTCACCTACTTTCCCGTGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGCTTCA
 CCGCCAGATCCGAGGCGCCGTGACCTGTCCCTCTATCCTCGAGAAGGGGGTCTCCAGCCGTGAGCT
 GGTGAAGAAGGTCTCCGATGTCATATGGAACAGCCTCAGCCGCTCCTACTTCAAGGATCGGGCCACATC
 CAGTCCCTTTCAGTTCATCACAGGCACAAATTGGACAGCTCCGGTGTGCCCTTGTCTGGTGGTGGG
 CCTGCCAGGCCCTGGGTCTCCGGATGTCCACCTGCCCTGTCTGAGGATCATGCCTGGTGTGTTGG
 GCCCAATGGGGAGCAGACAGCTGAGGTACCTGGCAGCGCAAGGGCAACGAGGACCGAGGGCCAGACA
 GTCAATGCCGGTGTGGCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCCTGTGACCCGAAGA
 TGGAGGTGGCGTTCATGGTGTGTGCCATCAACCCTCCATTGACCTGCACACCGACTCGCTGGAGCTTCT
 GCAGCTGCAGCAGAAGCTGCTCTGGCTGCTCTATGACCTGGGACATCTGGAAAGG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:


>RC403627 representing NM_130799
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDDVVRLFAAELGREPDLVLLSLVLGFVEHFLAVNRVIPTNPPELTFQPSPAPD
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDLSLYPREGGVSSREL VKKVSDVIWNSLSRSYFKDRAHI
 QSLFSFITGKLDSSGVAFVVGACQALGLRDVHLAL SEDHAWVVF GPNGEQTA EVTWHGKGNEDRRRGQT
 VNAGVAERSWL YLKGSYMRCDRKMEVAFMVCAINPSIDLHTDSLELLQLQKLLWLLYDLGHLER

SGPTRRRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


* The last codon before the Stop codon of the ORF

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NP_570711
RefSeq Size:	825 bp
RefSeq ORF:	1833 bp
Locus ID:	4221
Cytogenetics:	11q13.1
Domains:	Menin
Protein Families:	Druggable Genome, Transcription Factors
MW:	30.3 kDa
Gene Summary:	This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]