

Product datasheet for **RC403703**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	W423X
Affected Codon#:	423
Affected NT#:	1268
Nucleotide Mutation:	MEN1 Mutant (W423X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	MEN1
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	1266 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC403703 representing NM_130799
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCCTGTTTGC
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGCTCCTTTCCCTGGTCTGGCTTCGTGGAGCATTT
 TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCAGCCCCGCCCGCCGAC
 CCGCTGGCGGCCTCACCTACTTTCCCGTGGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGTTCA
 CCGCCAGATCCGAGGCGCGTGCACCTGTCCCTCTATCCTCGAGAAGGGGGTCTCCAGCCGTGAGCT
 GGTGAAGAAGGTCTCCGATGTCATATGGAACAGCCTCAGCCGCTCCTACTTCAAGGATCGGGCCACATC
 CAGTCCCTTTCAGTTCATCACAGGCACCAAATTGGACAGCTCCGGTGTGCCCTTTGCTGTGGTTGGGG
 CCTGCCAGGCCCTGGGTCTCCGGATGTCCACCTCGCCCTGTCTGAGGATCATGCCTGGGTAGTGTTTGG
 GCCCAATGGGGAGCAGACAGCTGAGGTACCTGGCAGGCAAGGGCAACGAGGACCGCAGGGGCCAGACA
 GTCAATGCCGGTGTGGCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCCTGTGACCCGAAGA
 TGGAGGTGGCGTTCATGGTGTGTGCCATCAACCCTCCATTGACCTGCACACCGACTCGCTGGAGCTTCT
 GCAGCTGCAGCAGAAGCTGCTCTGGCTGCTCTATGACCTGGGACATCTGGAAAGGTACCCCATGGCCTTA
 GGGAACTGGCAGATCTAGAGGAGCTGGAGCCACCCCTGGCCGGCCAGACCCACTCACCTCTACCACA
 AGGGCATTGCCTCAGCCAAGACCTACTATCGGGATGAACACATCTACCCCTACATGTACCTGGCTGGCTA
 CCACTGTCGAACCGCAATGTGCGGGAAGCCCTGCAGGCTGGGCGGACACGGCCACTGTATCCAGGAC
 TACAACCTACTGCCGGGAAGACGAGGAGATCTACAAGGAGTTCTTTGAAGTAGCCAATGATGTCATCCCA
 ACCTGCTGAAGGAGGCAGCCAGCTTGTGGAGCGGGCGAGGAGCGCCGGGGGAGCAAAGCCAGGGCAC
 CCAGAGCAAGGTTCCGCCCTCAGGACCTGAGTGCTTCGCCACCTGCTGCGATTCTACGACGGCATC
 TGCAAA

AG**GCGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGTTTAA

Protein Sequence:


>RC403703 representing NM_130799
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDVVRFLFAAELGREEPDLVLLSLVLFVEHFLAVNRVIPTNPVPELTFQPSAPD
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDSLYPREGVSSREL VKKVSDVIWNSLSRSYFKDRAHI
 QSLFSFITGKLDSSGVAFVAVGACQALGLRDVHLALSEDHAWVVFPGNGEQTA EVTWHGKGNEDRRGQT
 VNAGVAERSWL YLKGSYMRCDRKMEVAFMVCAINPSIDLHTDSLELLQLQKLLWLLYDLGHLERYPMAL
 GNLADLEELEPTPGRPDPLTL YHKGIASAKTYRDEHIYPYMYLAGYHCRNRNVREALQAWADTATVIQD
 YNYCREDEEIIYKEFFEYANDVIPNLLKEASLLEAGEERPGEQSQGTQSQGSALQDPECF AHLLRFYDGI
 CK

SGP**TRTRRLEQKLI SEEDLAANDILDYKDDDDK**V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq:

[NP_570711](#)

RefSeq Size:

1266 bp

RefSeq ORF:

1833 bp

Locus ID:

4221

Cytogenetics:

11q13.1

Domains:

Menin

Protein Families:

Druggable Genome, Transcription Factors

MW:

46.4 kDa

Gene Summary:

This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]