

Product datasheet for **RC403717**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Mutant ORF Clones |
| Product Name: | Menin (MEN1) (NM_130799) Human Mutant ORF Clone |
| Mutation Description: | Q450X |
| Affected Codon#: | 450 |
| Affected NT#: | 1348 |
| Nucleotide Mutation: | MEN1 Mutant (Q450X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA |
| Effect: | Multiple endocrine neoplasia 1 |
| Symbol: | Menin |
| Synonyms: | MEAI; SCG2 |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| Tag: | Myc-DDK |
| ACCN: | NM_130799 |
| ORF Size: | 1347 bp |
| Restriction Sites: | Sgfi-MluI |



[View online »](#)

ORF Nucleotide Sequence:

>RC403717 representing NM_130799
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCGCCCTGTTTG
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGCTCCTTTCCCTGGTCTGGCTTCGTGGAGCATTT
 TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCCAGCCCCGCCCGAC
 CCGCTGGCGGCCTCACCTACTTTCCCGTGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGTTCA
 CCGCCAGATCCGAGGCGCCGTGACCTGTCCCTCTATCCTCGAGAAGGGGGTCTCCAGCCGTGAGCT
 GGTGAAGAAGTCTCCGATGTCATATGGAACAGCCTCAGCCGCTCCTACTTCAAGGATCGGGCCACATC
 CAGTCCCTTTCAGTTCATCACAGGCACCAAATTGGACAGCTCCGGTGTGCCCTTGTGTGGTTGGGG
 CCTGCCAGGCCCTGGGTCTCCGGATGTCCACCTGCCTGTCTGAGGATCATGCCTGGGTAGTGTGGG
 GCCCAATGGGGAGCAGACAGCTGAGGTACCTGGCACGGCAAGGGCAACGAGGACCGCAGGGGCCAGACA
 GTCAATGCCGGTGTGGCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCCTGTGACCCGAAGA
 TGGAGGTGGCGTTCATGGTGTGTGCCATCAACCCTCCATTGACCTGCACACCGACTCGCTGGAGCTTCT
 GCAGCTGCAGCAGAAGCTGCTCTGGTCTCTATGACCTGGGACATCTGGAAAGGTACCCCATGGCCTTA
 GGGAACTGGCAGATCTAGAGGAGCTGGAGCCACCCCTGGCCGGCCAGACCCACTCACCTCTACCACA
 AGGGCATTGCCTCAGCCAAGACCTACTATCGGGATGAACACATCTACCCCTACATGTACCTGGCTGGCTA
 CCACTGTCGAACCGCAATGTGCGGGAAGCCCTGCAGGCTGGCGGACACGGCCACTGTATCCAGGAC
 TACAACCTACTGCCGGGAAGACGAGGAGATCTACAAGGAGTTCTTTGAAGTAGCCAATGATGTCATCCCA
 ACCTGCTGAAGGAGGCAGCCAGCTTGTGGAGGCGGGCAGGAGCGGCCGGGGGAGCAAAGCCAGGGCAC
 CCAGAGCCAAGGTTCCGCCCTCCAGGACCTGAGTGCTTCGCCACCTGCTGCGATTCTACGACGGCATC
 TGCAAATGGGAGGAGGGCAGTCCACGCCTGTGCTGCAGTGGGCTGGGCCACCTTTCTGTGCAGTCCC
 TAGGCCGTTTTGAGGGA

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:


>RC403717 representing NM_130799
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDVVRLLFAAELGREPDVLLSLVLFVEHFLAVNRVIPTNVPELTFQPSAPD
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDLSLYPREGGVSSRELVKKVSVDIWNLSRSYFKDRAHI
 QSLFSFITGKLDSSGVAFVAVGACQALGLRDVHLALSEDHAWVVFPGNGEQTAEVTHWGKGNEDRRGQT
 VNAGVAERSWL YLKGSYMRCDRKMEVAFMVCAINPSIDLHTDSLELLQLQKLLWLLYDLGHLERYPMAL
 GNLADLEELEPTPGRPDPLTYHKGIASAKTYRDEHIYPMYLAGYHCRNRNVREALQAWADTATVIQD
 YNYCREDEEIIYKEFFEYANDVIPNLLKEAASLLEAGEERPGEQSQGTQSQGSALQDPECFALLRFYDGI
 CKWEEGSPTPVLHVGWATFLVQSLGRFEG

SGPTRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NP_570711](#)

RefSeq Size:

1347 bp

RefSeq ORF:

1833 bp

Locus ID:

4221

Cytogenetics:

11q13.1

Domains:

Menin

Protein Families:

Druggable Genome, Transcription Factors

MW:

49.4 kDa

Gene Summary:

This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]