

Product datasheet for **RC403729**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	R527X
Affected Codon#:	527
Affected NT#:	1579
Nucleotide Mutation:	MEN1 Mutant (R527X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	Menin
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	1578 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC403729 representing NM_130799
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCGCCCTGTTT
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGTCTCTTCCCTGGTGTGGCTTCGTGGAGCATTT
 TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCCAGCCCCGCCCCGAC
 CCGCTGGCGGCCTCACCTACTTCCCGTGGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGCTTCA
 CCGCCAGATCCGAGGCGCCGTGACCTGTCCCTCTATCTCGAGAAGGGGTGTCTCCAGCCGTGAGCT
 GGTGAAGAAGGTCTCCGATGTATATGGAACAGCCTCAGCCGCTCTACTTCAAGGATCGGGCCACATC
 CAGTCCCTTTCAGTTCATCACAGGCACCAAATTGGACAGCTCCGGTGTGGCCTTTGCTGTGGTTGGG
 CCTGCCAGGCCCTGGGTCTCCGGATGTCCACCTCGCCCTGTCTGAGGATCATGCCTGGGTAGTGTGG
 GCCAATGGGGAGCAGACAGCTGAGGTACCTGGCAGGCAAGGGCAACGAGGACCGCAGGGGCCAGACA
 GTCAATGCCGGTGTGGCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCCTGTGACCGCAAGA
 TGGAGGTGGCGTTCATGGTGTGTGCCATCAACCCTCCATTGACCTGCACACCGACTCGCTGGAGCTTCT
 GCAGCTGCAGCAGAAGCTGCTCTGGTGTCTATGACCTGGGACATCTGGAAAGGTACCCCATGGCCTTA
 GGGAACTGGCAGATCTAGAGGAGCTGGAGCCACCCCTGGCCGGCCAGACCCACTCACCTCTACCACA
 AGGGCATTGCCTCAGCCAAGACCTACTATCGGGATGAACACATCTACCCCTACATGTACCTGGCTGGTA
 CCACTGTGCAACCGCAATGTGCGGGAAGCCCTGCAGGCTGGGCGGACAGGCCACTGTATCCAGGAC
 TACAACCTGCCGGGAAGACGAGGAGATCTACAAGGAGTCTTTGAAGTAGCCAATGATGTCATCCCCA
 ACCTGCTGAAGGAGGAGCCAGCTTGTGGAGCGGGGAGGAGCGGCCGGGGGAGCAAAGCCAGGCAC
 CCAGAGCCAAGGTTCCGCCCTCCAGGACCTGAGTGTCTCGCCACCTGCTGCGATTCTACGACGGCATC
 TGCAAATGGGAGGAGGGCAGTCCACGCTGTGCTGCAGTGGGCTGGGCCACCTTTCTGTGCAGTCCC
 TAGGCCGTTTTGAGGGACAGGTGCGGCAGAAGGTGCGCATAGTGAAGGAGGCGGAGGCGGCCGAGGC
 CGAGGAGCCGTGGGGCAGGAAGCCCGGAAGGCCGGCGGGGCCACGGCGGGAGTCCAAGCCAGAG
 GAGCCCCCGCCCAAGAAGCCAGCACTGGACAAGGGCTGGGCACCGCCAGGGTGCAGTGTGAGGAC
 CCCCCGGAAGCTCTGGGACTGTGCTGGCACAGCC

AGCGGACCGACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGTTTAA

Protein Sequence:


>RC403729 representing NM_130799
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDVVRFLFAAELGREPDLVLLSLVLGFVEHFLAVNRVIPTNPVELTFQSPAPD
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDLSLYPREGGVSSRELVKKVSVDIWNLSRSYFKDRAHI
 QSLFSFITGTKLDSSGVAFVVGACQALGLRDVHLALSEDHAWVFGPNGEQTAEVTHWGKGNEDRRGQT
 VNAGVAERSWLYLKGSYMRCDRMEVAFMVCAINPSIDLHTDSLELLQLQKLLWLLYDLGHLERYPMAL
 GNLADLEELEPTGRPDPLTYHKGIASAKTYRDEHIYPMYLAGYHCRNRNVREALQAWADTATVIQD
 YNYCREDEEIIYKEFFEVDVIPNLLKEAASLLEAGEERPGEQSQGTQSQGSALQDPECFHLLRFYDGI
 CKWEEGSPTPVLHVGWATFLVQSLGRFEGQVRQKVRIVSREAAEAEEPPWGEAREGRRRGPRESKPE
 EPPPPKPPALDKGLGTGQAVSGPPRPPGTVAGTA

SGPTRRRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NP_570711](#)

RefSeq Size:

1578 bp

RefSeq ORF:

1833 bp

Locus ID:

4221

Cytogenetics:

11q13.1

Domains:

Menin

Protein Families:

Druggable Genome, Transcription Factors

MW:

57.9 kDa

Gene Summary:

This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]