

Product datasheet for **RC403730**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	E530X
Affected Codon#:	530
Affected NT#:	1588
Nucleotide Mutation:	MEN1 Mutant (E530X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	Menin
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	1587 bp
Restriction Sites:	Sgfl-MluI



[View online »](#)

ORF Nucleotide Sequence:

>RC403730 representing NM_130799
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCCTGTTTG
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGCTCCTTTCCCTGGTGCTGGGCTTCGTGGAGCATTT
 TCTGGCTGTCAACCGCGTCATCCCTACCAACGTTCCCGAGCTCACCTTCAGCCAGCCCGCCCCGAC
 CCGCTGGCGGCTCACCTACTTTCCCGTGCCGACCTGTCTATCATCGCCGCCCTCTATGCCGCTTCA
 CCGCCAGATCCGAGGCGCCGTCGACCTGTCCCTCTATCCTCGAGAAGGGGTGTCTCCAGCCGTGAGCT
 GGTGAAGAAGGTCTCCGATGTCATATGGAACAGCCTCAGCCGCTCCTACTTCAAGGATCGGGCCACATC
 CAGTCCCTCTCAGTTCATCACAGGCACCAATTGGACAGCTCCGGTGTGGCCTTTGCTGTGGTTGGG
 CCTGCCAGGCCCTGGGTCTCCGGGATGTCCACCTCGCCTGTCTGAGGATCATGCCTGGGTAGTGTGG
 GCCAATGGGGAGCAGACAGCTGAGGTACCTGGCAGGCAAGGGCAACGAGGACCGAGGGGCCAGACA
 GTCAATGCCGGTGTGGCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCCTGTGACCGCAAGA
 TGGAGGTGGCGTTCATGGTGTGTGCCATCAACCCTCCATTGACCTGCACACCGACTCGCTGGAGCTTCT
 GCAGCTGCAGCAGAAGCTGCTCTGGTGCTCTATGACCTGGGACATCTGGAAAGGTACCCCATGGCCTTA
 GGGAACCTGGCAGATCTAGAGGAGCTGGAGCCACCCCTGGCCGCCAGACCCACTCACCTCTACCACA
 AGGGCATTGCCTCAGCCAAGACCTACTATCGGGATGAACACATCTACCCCTACATGTACCTGGCTGGCTA
 CCACTGTCGCAACCGCAATGTGCGGGAAGCCCTGCAGGCTGGCGGACACGGCCACTGTATCCAGGAC
 TACAATACTGCCGGAAGACGAGGAGATCTACAAGGAGTTCTTTGAAGTAGCCAATGATGTCATCCCCA
 ACCTGCTGAAGGAGGAGCCAGCTTGTGGAGCGGGCGAGGAGCGGCCGGGGGAGCAAGCCAGGGCAC
 CCAGAGCCAAGGTTCCGCCCTCCAGGACCTGAGTGCTTCGCCACCTGCTGCGATTCTACGACGGCATC
 TGCAAATGGGAGGAGGAGTCCACGCTGTGCTGCACGTGGGCTGGGCCACCTTTCTGTGCAGTCCC
 TAGGCCGTTTTGAGGGACAGGTGCGGCAGAAGGTGCGCATAGTGAGCCGAGAGGCCGAGGCGGCCAGGC
 CGAGGAGCCGTGGGGCGAGGAAGCCCGGAAGGCCGCGGCGGGGCCACGGCGGAGTCCAAGCCAGAG
 GAGCCCCCGCCCAAGAAGCCAGCACTGGACAAGGGCTGGGCACCGGCCAGGCTGCAGTGTGAGGAC
 CCCCCGAAGCTCTGGGACTGTGCTGGCACAGCCGAGGCCCT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA TAAGGTTTAA

Protein Sequence:

>RC403730 representing NM_130799
 Red=Cloning site Green=Tags(s)


MGLKAAQKTLFPLRSIDVVRLLFAELGREEDLVLLSLVLGFVEHFLAVNRVIPTNPVELTFQPSAPD
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDLSLYPREGGVSSRELVKKVSVDIWNLSRSYFKDRAHI
 QSLFSFITGKLDSSGVAFVVGACQALGLRDVHLALSEDHAWVFGPNGEQTAEVTHWGKGNEDRRGQT
 VNAGVAERSWLYLKSYMRCRMEVAFMVCAINPSIDLHTDSLELLQLQQKLLWLLYDLGHLERYPMAL
 GNLADLEELEPTPRPDPLTYHKGIASAKTYRDEHIYPMYLAGYHCRNRNVREALQAWADTATVIQD
 YNYCREDEEIIYKEFFEANDVIPNLLKEASLLEAGEERPGEQSQGTQSQGSALQDPECFAHLLRFYDGI
 CKWEEGSPTPVLHVGWATFLVQSLGRFEGQVRQKVRIVSREAEAAEAEPWGEEAREGRRRGPRRESKPE
 EPPPPKPPALDKGLGTGQGAVSGPPRPPGTVAGTARGP

SGPTRRRLE**QKLI**SEED**LA**AND**ILDYKDDDDKV**

Restriction Sites:

SgfI-MluI

Cloning Scheme:



OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NP_570711](#)

RefSeq Size:

1587 bp

RefSeq ORF:

1833 bp

Locus ID:

4221

Cytogenetics:

11q13.1

Domains:

Menin

Protein Families:

Druggable Genome, Transcription Factors

MW:

58.2 kDa

Gene Summary:

This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]