

Product datasheet for **RC403742**

Menin (MEN1) (NM_130799) Human Mutant ORF Clone

Product data:

Product Type:	Mutant ORF Clones
Product Name:	Menin (MEN1) (NM_130799) Human Mutant ORF Clone
Mutation Description:	K559X
Affected Codon#:	559
Affected NT#:	1675
Nucleotide Mutation:	MEN1 Mutant (K559X), Myc-DDK-tagged ORF clone of Homo sapiens multiple endocrine neoplasia I (MEN1), transcript variant 2 as transfection-ready DNA
Effect:	Multiple endocrine neoplasia 1
Symbol:	MEN1
Synonyms:	MEAI; SCG2
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
Tag:	Myc-DDK
ACCN:	NM_130799
ORF Size:	1674 bp
Restriction Sites:	Sgfl-Mlul



[View online »](#)

ORF Nucleotide Sequence:

>RC403742 representing NM_130799
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGGCTGAAGGCCGCCAGAAGACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTGCGCCCTGTTT
 CTGCCGAGCTGGGCCGAGAGGAGCCGGACCTGGTGCTCCTTTCCCTGGTCTGGCTTCGTGGAGCATTT
 TCTGGCTGTCAACCGCTCATCCCTACCAACGTTCCCGAGCTCACCTTCCAGCCCAGCCCCGCCCCGAC
 CCGCTGGCGGCCTCACCTACTTTCCCGTGCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGTTCA
 CCGCCAGATCCGAGGCGCCGTGACCTGTCCCTCTATCTCGAGAAGGGGTGTCTCCAGCCGTGAGCT
 GGTGAAGAAGGTCTCCGATGTATGGAACAGCCTCAGCCGCTCTACTTCAAGGATCGGGCCACATC
 CAGTCCCTTTCAGTTCATCACAGGCACAAATTGGACAGCTCCGGTGTGCCCTTGTGTGGTGGGG
 CCTGCCAGGCCCTGGGTCTCCGGATGTCCACCTCGCCCTGTCTGAGGATCATGCCTGGGTAGTGTGG
 GCCAATGGGGAGCAGACAGCTGAGGTACCTGGCACGGCAAGGCAACGAGGACCGCAGGGCCAGACA
 GTCAATGCCGGTGTGGCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCCTGTGACCGCAAGA
 TGGAGGTGGCGTTCATGGTGTGTGCCATCAACCCTCCATTGACCTGCACACCGACTCGCTGGAGCTTCT
 GCAGCTGCAGCAGAAGCTGCTCTGGCTGCTCTATGACCTGGGACATCTGGAAAGGTACCCCATGGCCTTA
 GGGAACTGGCAGATCTAGAGGAGCTGGAGCCACCCCTGGCCGGCCAGACCCACTCACCTCTACCACA
 AGGGCATTGCCTCAGCCAAGACCTACTATCGGGATGAACACATCTACCCCTACATGTACCTGGCTGGTA
 CCACTGTCGAACCGCAATGTGCGGGAAGCCCTGCAGGCTGGCGGACACGGCCACTGTATCCAGGAC
 TACAACCTACTGCCGGGAAGACGAGGAGATCTACAAGGAGTCTTTGAAGTAGCCAATGATGTCATCCCA
 ACCTGCTGAAGGAGCCAGCTTGTGGAGCGGGGAGGAGCGGCCGGGGGAGCAAGCCAGGGCAC
 CCAGAGCCAAGGTTCCGCCCTCCAGGACCTGAGTGTTCGCCACCTGCTGCGATTCTACGACGGCATC
 TGCAAATGGGAGGAGGGCAGTCCACGCTGTGCTGCAGTGGCTGGGCCACCTTTCTGTGCAGTCCC
 TAGGCCGTTTTGAGGGACAGGTGCGGCAGAAGGTGCGCATAGTGAAGGAGGCGGAGGCGGCCAGGC
 CGAGGAGCCGTGGGGCAGGAAGCCCGGAAGGCCGGCGGGGCCACGGCGGGAGTCCAAGCCAGAG
 GAGCCCCCGCCCAAGAAGCCAGCACTGGACAAGGGCTGGGCACCGCCAGGGTGCAGTGTGAGGAC
 CCCCCCGAAGCCTCTGGGACTGTGCTGGCACAGCCCGAGGCCCTGAAGGTGGCAGCACGGCTCAGGT
 GCCAGCACCCGACATACCACCGCCGAGGGTCCAGTGCTCACTTTCCAGAGTGAGAAGATG

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGA **TAAGTTTAA**

Protein Sequence:

>RC403742 representing NM_130799
 Red=Cloning site Green=Tags(s)

MGLKAAQKTLFPLRSIDDVVRLF AAELGREPDLVLLSLVLFVVEHFLAVNRVIPTNPVPELTFQPSAPD
 PPGGLTYFPVADLSIIAALYARFTAQIRGAVDLSLYPREGGVSSRELVKKVSVDVIWNSLSRSYFKDRAHI
 QSLFSFITGKLDSSGVAFAVVGACQALGLRDVHLALSEDHAWVVFPGNGEQTA EVTWHGKGNEDRRGQT
 VNAGVAERSWL YLKGSYMRCDRMEVAFMVCAINPSIDLHTDSLELLQLQKLLWLLYDLGHLERYPMAL
 GNLADLEELPTPGRPDPLTYHKGIASAKTYRDEHIYPYMYLAGYHCRNRNVREALQAWADTATVIQD
 YNYCREDEEIIYKEFFEIVANDVIPNLLKEASLLEAGEERPGEQSQGTQSQGSALQDPECFHLLRFYDGI
 CKWEEGSPTPVLHVGWATFLVQSLGRFEGVQRQKVRIVSREAEAAEAEPEWGEAREGRRRGPRESKPE
 EPPPPKPKALDKGLGTGQAVSGPPRPPGTVAGTARGPEGGSTAQVPAPAASPPPEGPVLTQSEKM

SGPTRRRLEQKLI**SEEDLAANDILDYKDDDDKV**

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

RefSeq:

[NP_570711](#)

RefSeq Size:

1674 bp

RefSeq ORF:	1833 bp
Locus ID:	4221
Cytogenetics:	11q13.1
Domains:	Menin
Protein Families:	Druggable Genome, Transcription Factors
MW:	61.4 kDa
Gene Summary:	This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]