

Product datasheet for **RC600018**

FGFR1 (NM_023110) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FGFR1 (NM_023110) Human Tagged ORF Clone
Tag:	DDK-His
Symbol:	FGFR1
Synonyms:	bFGF-R-1; BFGFR; CD331; CEK; ECCL; FGFBR; FGFR-1; FLG; FLT-2; FLT2; HBGFR; HH2; HRTFDS; KAL2; N-SAM; OGD
Mammalian Cell Selection:	None
Vector:	pCMV6-XL5-DDK-His (PS100068)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC600018 representing leader sequence plus the extracellular domain region of NM_023110

Red=Cloning site Blue=ORF Green=Tags(s)

GTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCTGGTACCGAGGAGATCCGCCGCCG
CGATCGCC

ATGTGGAGCTGGAAGTGCCTCCTTCTGGGCTGTGCTGGTACAGCCACTCTGCACCCTAGGCCGT
CCCCGACCTTGCCTGAACAAGCCAGCCCTGGGGAGCCCTGTGGAAGTGGAGTCCTTCTGGTCCACCC
CGGTGACCTGCTGCAGCTTCGTGTGGCTGCGGGACGATGTGCAGAGCATCAACTGGTGCAGGACGGG
GTGCAGCTGGCGAAAGCAACCGCACCCGCATCACAGGGGAGGAGGTGGAGGTGCAGGACTCCGTGCCCC
CAGACTCCGGCCTCTATGCTTGCCTAACAGCAGCCCTCGGGCAGTGACACCACCTACTTCTCCGTCAA
TGTTTCAGATGCTCTCCCTCCTCGGAGGATGATGATGATGATGACTCCTTTCAGAGGAGAAAGAA
ACAGATAACACCAACCAAAACCGTATGCCCGTAGCTCCATATTGGACATCCCAGAAAAGATGAAAAGA
AATTGCATGCAGTGCCGGCTGCCAAGACAGTGAAGTTCAAATGCCCTTCCAGTGGGACCCAAACCCAC
ACTGCGCTGGTTGAAAAATGGCAAAGAATTCAAACCTGACCACAGAATTGGAGGCTACAAGTCCGTAT
GCCACCTGGAGCATCATAATGGACTCTGTGGTGCCTCTGACAAGGGCAACTACACCTGCATTGTGGAGA
ATGAGTACGGCAGCATCAACCACACATACCAGTGGATGTCGTGGAGCGGTCCCCTACCGGCCATCCT
GCAAGCAGGGTTGCCCGCAACAAAACAGTGGCCCTGGGTAGCAACGTGGAGTTCATGTGAAGGTGAC
AGTGACCCGCAGCCGCACATCCAGTGGCTAAAGCACATCGAGGTGAATGGGAGCAAGATTGGCCAGACA
ACCTGCCTTATGTCCAGATCTTGAAGACTGCTGGAGTTAATACCACCGACAAAGAGATGGAGGTGCTTCA
CTTAAGAAATGTCTCCTTTGAGGACGCAGGGGAGTATACGTGCTTGGCGGGTAACTCTATCGGACTCTCC
CATCACTCTGCATGGTTGACCGTCTTGGAAAGCCCTGGAAGAGAGGCCGGCAGTGATGACCTCGCCCTGT
ACCTGGAG

ACGCGTTCAGGGCGACTACAAGGATGACGACGATAAGGGATCTCATCATCACCATCACCATTAATGAGATC
TGGTACCGATATCAAGCTTGTGACTCTAGA

Protein Sequence:

>RC600018 representing signal peptide plus the extracellular domain region of NM_023110

Red=Cloning sites Green= DDK and 6XHIS Tags


MWSWKCLLFWAVLVLTATLCTARPSPTLPEQAQPWVAPVEVESFLVHPGDLQLRCRLRDDVQSINWLRDG
VQLAESNRTRITGEEVEVQDSVPADSGLYACVTSSPSGSDTTYFSVNVSDALPSEDDDDDDSSSEEKE
TDNTKPNRMPVAPYWTSPEKMEKLLHAVPAAKTVKFKCPSGTPNPTLRWLKNGKEFKPDHRIGGYKVRY
ATWSIIMDSVVPVSDKGNVYTCIVENEYGSINHTYQLDVVERSPHRPILQAGLPANKTVALGNSVEFMCKVY
SDPQPHIQWLKHIEVNGSKIGPDNLPYVQILKTAGVNTTDKEMEVLHLRNVSFEDAGEYTCLAGNSIGLS
HHSAWLTVLEALEERPAMTSPLYLE

TRSGTRSGDYKDDDDKGSHHHHHH

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_023110

ORF Size: 1128 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the extra cellular domain of the protein with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_023110.2](#), [NP_075598.2](#)

RefSeq Size: 5917 bp

RefSeq ORF: 2469 bp

Locus ID: 2260

UniProt ID: [P11362](#)

Cytogenetics: 8p11.23

Domains:	ig, IGc2, IG
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Adherens junction, MAPK signaling pathway, Melanoma, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton
MW:	41.8 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the fibroblast growth factor receptor (FGFR) family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds both acidic and basic fibroblast growth factors and is involved in limb induction. Mutations in this gene have been associated with Pfeiffer syndrome, Jackson-Weiss syndrome, Antley-Bixler syndrome, osteoglophonic dysplasia, and autosomal dominant Kallmann syndrome 2. Chromosomal aberrations involving this gene are associated with stem cell myeloproliferative disorder and stem cell leukemia lymphoma syndrome. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008]</p>