

Product datasheet for **RG200138**

RBM23 (NM_018107) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RBM23 (NM_018107) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RBM23
Synonyms:	CAPERbeta; PP239; RNPC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG200138 representing NM_018107 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCATCTGATGACTTTGACATAGTGATTGAGGCCATGCTGGAAGCTCCCTATAAAAAAGAAGAGGATG
 AGCAACAAAGGAAAGAAGTTAAAAAGGATTATCCTAGCAATACCACCAGCAGCACCAGCAACAGTGGCAA
 TGAGACCAAGTGAAGCAGCACCATCGGGGAGACAAGCAATCGTAGTCGAGATCGGGATCGGTATAGACGG
 AGAAATAGTCGGAGCCGAAGTCCAGGTCGGCAGTGTCTCACCAGTACCGGTAGCTGGGATCGTCGACATG
 GTAGTGAGTCGCGAAGTCCGGACCATCGTCGTGAGGATCGTGTGCATTACAGGAGTCTCCACTTGCCAC
 TGGTTATAGATATGGACACAGTAAGAGTCCTCATTTTCAGAGAGAAGAGCCAGTCAGGGAGCCAGTTGAT
 AATCTGAGTCCTGAGGAGCGTGATGCCCGCACAGTTTTCTGTATGCAGTTAGCTGCCCGAATTCGGCCTC
 GAGATCTGGAGGACTTTTTCTCTGCTGTAGGCAAGGTTTCGCGATGTACGTATCATTTTCAGATCGGAATC
 ACGTCGTTCTAAGGGCATTGCCTACGTGGAATTCGTGAAATCCAGTCTGTGCCACTGGCCATTGGGCTG
 ACTGGGCAGCGGTTGCTGGGAGTGCCTATCATTGTACAGGCTTCACAGGCAGAGAAAAACCGACTGGCAG
 CCATGGCCAACAACCTGCAAAAGGGCAATGGTGGACCAATGCGCCTCTATGTGGGTTCCCTGCACCTCAA
 TACTACTGAAGACATGCTCCGGGCGCATCTTTGAGCCCTTTGGTAAATTGATAATATTGCTCTGATGAAG
 GACTCAGATACAGGCCGCTCTAAAGGTTATGGTTTCATCAGTTCTCTGATTCTGAGTGTGCCGGCGGG
 CCCTGGAACAGTTGAATGGGTTTGAGCTTGCTGGTCGACCTATGAGGGTTGGCCATGTGACTGAGCGACT
 GGATGGTGGCACAGACATCACTTTCTGATGGGACCAGGAGCTGGATCTGGGATCAGCAGGTGGACGT
 TTTTCAGCTCATGGCAAACTGGCAGAAGGCGCTGGAATCCAAGTCCCAAGCACTGCTGCTGCTGCTG
 CCGCCGCCGCCAGGCTGCTGCCTTGCAACTGAATGGAGCAGTTCCTTGGGGGCCCTGAATCCAGC
 AGCTCTGACTGCTCTGAGTCCAGCCCTGAACCTTGCTCCAGTGTTCAGCTCTCCAGCCTCTTACC
 CCCAGACCATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG200138 representing NM_018107
 Red=Cloning site Green=Tags(s)

MASDDFDIVIEAMLEAPYKKEEDEQQRKEVKKDYPSTTSSTNSGNETSGSSTIGETSNRSRDRDRYRR
 RNSRSRSPGRQCRHRSRSDRRHGSESRSDHRREDRVHYRSPPLATGYRYGHSKSPHFREKSPVREPVD
 NLSPEEDARTVFCMQLAARIRPRDLEDFFSAVGKVRDVRIISDRNSRRSKGIAYVEFCEIQSVPLAIGL
 TGQRLLGVPIIVQASQAENRLAAMANNLQKNGGPMRLYVGSLSLHFNITEDMLRGIFEPFGKIDNIVLMK
 DSDTGRSKGYGFIITFSDSECARRALEQLNGFELAGRPMRVGHVTERLDGGTDITFPDGDQELDLGSAGGR
 FQLMAKLAEGAGIQLPSTAAAAAAAAAQAALQLNGAVPLGALNPAALTALSPALNLSAQCFQLSSLFT
 PQTM

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_018107

ORF Size: 1272 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_018107.3](#), [NP_060577.2](#)

RefSeq Size: 2576 bp

RefSeq ORF: 1272 bp

Locus ID: 55147


UniProt ID: [Q86U06](#)

Cytogenetics: 14q11.2

Domains: RRM

Gene Summary: This gene encodes a member of the U2AF-like family of RNA binding proteins. This protein interacts with some steroid nuclear receptors, localizes to the promoter of a steroid-responsive gene, and increases transcription of steroid-responsive transcriptional reporters in a hormone-dependent manner. It is also implicated in the steroid receptor-dependent regulation of alternative splicing. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG200138