

Product datasheet for **RG200324**

ASF1A (NM_014034) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ASF1A (NM_014034) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ASF1A
Synonyms: CGI-98; CIA; HSPC146
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG200324 representing NM_014034
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCAAAGTTTCAGGTGAACAATGTAGTGGTCTGGATAACCCTTCTCCTTTCTACAACCCGTTCCAGT
 TCGAGATCACCTTCGAGTGCATCGAGGACCTGTCTGAAGACTTGGAAATGGAAAATATCTATGTGGGCTC
 TGCAGAAAGTGAAGAATACGATCAAGTTTTAGACTCTGTTTTAGTGGTCTGTTCCCGCAGGAAGGCAT
 ATGTTTGTATTTTCAGGCTGATGCACCTAATCCAGGACTCATTCCAGATGCAGATGCAGTAGGCGTAACTG
 TTGTGCTAATTACTTGTACCTATCGAGGACAAGAATTTATTAGAGTTGGCTATTATGTAATAATGAATA
 TACTGAGACAGAATTAAGGGAAAATCCACCAGTAAAACCAGACTTTTCTAAGCTTCAAAGGAATATTTTG
 GCATCTAATCCCAGGTCACAAGATCCACATTAATTGGGAAGATAACACAGAAAACTGGAAGATGCAG
 AGAGCAGTAATCCAAATCTACAGTCACTTCTTTCAACAGATGCATTACCTTCAGCATCAAAGGGATGGTC
 CACATCAGAAAATCACTAAATGTCATGTTAGAATCCCACATGGACTGCATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG200324 representing NM_014034
 Red=Cloning site Green=Tags(s)

MAKVQVNNVVLDNPSPFYNPFFQFEITFECIEDLSEDLWKIIYVGSASEEYDQVLDVSLVGPVPAGRH
 MFVVFQADAPNPLIPDADAVGVTVVLITCTYRQQEFIRVGYVNNYETETELRENPPVKPDFSKLQRNII
 ASNPRVTRFHINWEDNTEKLEDAESSNPNLQSLSTDALPSASKGWSTSENSLNVMLESHMDCM

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_014034

ORF Size: 612 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014034.3](#)

RefSeq Size: 2525 bp

RefSeq ORF: 615 bp

Locus ID: 25842


UniProt ID: [Q9Y294](#)

Cytogenetics: 6q22.31

Domains: Anti-silence

Gene Summary: This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG200324