

Product datasheet for **RG200367**

ERH (NM_004450) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ERH (NM_004450) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ERH
Synonyms: DROER
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG200367 representing NM_004450
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTCACACCATTTTGTGGTACAGCCTACCAAGAGGCCAGAAGGCAGAACTTATGCTGACTACGAAT
CTGTGAATGAATGCATGGAAGGTGTTTGTAAAATGTATGAAGAACATCTGAAAAGAATGAATCCCAACAG
TCCCTCTATCACATATGACATCAGTCAGTTGTTTGATTCATCGATGATCTGGCAGACCTCAGCTGCCTG
GTTTACCGAGCTGATACCCAGACATACCAGCCTTATAACAAAGACTGGATTAAGAGAAGATCTACGTGC
TCCTTCGTCGGCAGGCCAACAGGCTGGGAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG200367 representing NM_004450
Red=Cloning site Green=Tags(s)
MSHTILLVQPTKRPEGRYADYESVNECMEGVCKMYEEHLKRMNPNSPITYDISQLFDFIDDLADLSCL
VYRADTQTYQPYNKDWIKEKIYVLLRRQAQAGK

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_004450

ORF Size: 312 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_004450.3](#)

RefSeq Size: 815 bp

RefSeq ORF: 315 bp

Locus ID: 2079
UniProt ID: [P84090](#)
Cytogenetics: 14q24.1
Domains: ER
Gene Summary: May have a role in the cell cycle.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG200367