

Product datasheet for **RG200469**

GRB2 (NM_002086) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GRB2 (NM_002086) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	GRB2
Synonyms:	ASH; EGFRBP-GRB2; Grb3-3; MST084; MSTP084; NCKAP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG200469 representing NM_002086 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGCCATCGCCAAATATGACTTCAAAGCTACTGCAGACGACGAGCTGAGCTTCAAAGGGGGGACA
TCCTCAAGGTTTTGAACGAAGAATGTGATCAGAAGTGGTACAAGGCAGAGCTTAATGGAAAAGACGGCTT
CATTCCCAAGAACTACATAGAAATGAAACCACATCCGTGGTTTTTTGGCAAATCCCCAGAGCCAAGGCA
GAAGAAATGCTTAGCAAACAGCGGCACGATGGGGCCTTTCTTATCCGAGAGAGTGAGAGCGCTCCTGGGG
ACTTCTCCCTCTCTGTCAAGTTTGGAAACGATGTGCAGCACTTCAAGGTGCTCCGAGATGGAGCCGGGAA
GTACTTCTCTGGGTGGTGAAGTTCAATTTTGAATGAGCTGGTGGATTATCACAGATCTACATCTGTC
TCCAGAAACCAGCAGATATTCCTGCGGGACATAGAACAGGTGCCACAGCAGCCGACATACGTCCAGGCC
TCTTTGACTTTGATCCCCAGGAGGATGGAGAGCTGGGCTTCCGCCGGGGAGATTTTATCCATGTCATGGA
TAACTCAGACCCCAACTGGTGAAAGGAGCTTGCCACGGGCAGACCGGCATGTTTCCCCGAATTATGTC
ACCCCGTGAACCGAACGTC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG200469 representing NM_002086
 Red=Cloning site Green=Tags(s)

MEAIKYDFKATADDELSFKRGDILKVLNEECDQNWYKAELNGKDGFIKPNYIEMKPHPWFFGKIPRAKA
 EEMLSKQRHDGAFLIRESESAPGDFSLSVKFGNDVQHFKVLRDGAGKYFLWVVKFNSLNELVDYHRSTSV
 SRNQIQIFLRDIEQVPQQPTYVQALFDFDPQEDGELGFRRGDFIHVMDNSDPNWWKGACHGQTMGFPRNYV
 TPVNRNV

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_002086

ORF Size: 651 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_002086.5
RefSeq Size:	3317 bp
RefSeq ORF:	654 bp
Locus ID:	2885
UniProt ID:	P62993
Cytogenetics:	17q25.1
Domains:	SH2, SH3
Protein Families:	Druggable Genome
Protein Pathways:	Acute myeloid leukemia, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Dorso-ventral axis formation, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Focal adhesion, Gap junction, Glioma, GnRH signaling pathway, Insulin signaling pathway, Jak-STAT signaling pathway, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pathways in cancer, Prostate cancer, Renal cell carcinoma, T cell receptor signaling pathway
Gene Summary:	The protein encoded by this gene binds the epidermal growth factor receptor and contains one SH2 domain and two SH3 domains. Its two SH3 domains direct complex formation with proline-rich regions of other proteins, and its SH2 domain binds tyrosine phosphorylated sequences. This gene is similar to the Sem5 gene of C.elegans, which is involved in the signal transduction pathway. Two alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG200469