

Product datasheet for RG200482

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Casein Kinase 1 delta (CSNK1D) (NM_001893) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Casein Kinase 1 delta (CSNK1D) (NM_001893) Human Tagged ORF Clone

Tag: TurboGFP

Symbol: Casein Kinase 1 delta

Synonyms: ASPS; CKI-delta; CKId; CKIdelta; FASPS2; HCKID

Mammalian Cell

Selection:

Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide >RG200482 representing NM_001893

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGGAGCTGAGAGTCGGGAACAGGTACCGGCTGGGCCGGAAGATCGGCAGCGGCTCCTTCGGAGACATCT ATCTCGGTACGGACATTGCTGCAGGAGAAGAGGTTGCCATCAAGCTTGAATGTGTCAAAACCAAACACCC TGCGGGGCAGAGGGGGACTACAACGTCATGGTGATGGAGCTGCTGGGGCCAAGCCTGGAGGACCTCTTCA ACTTCTGCTCCAGGAAATTCAGCCTCAAAACCGTCCTGCTGCTGCTGACCAAATGATCAGTCGCATCGA ATACATTCATCAAAGAACTTCATCCACCGGGATGTGAAGCCAGACAACTTCCTCATGGGCCTGGGGAAG ACATCCCCTATCGTGAGAACAAGAACCTCACGGGGACGGCGCGGTACGCCTCCATCAACACGCACCTTGG AATTGAACAATCCCGAAGAGATGACTTGGAGTCTCTGGGCTACGTGCTAATGTACTTCAACCTGGGCTCT CTCCCCTGGCAGGGGCTGAAGGCTGCCACCAAGAGAAATACGAAAGGATTAGCGAGAAGAAAATGT CCACCCCATCGAAGTGTTGTGTAAAGGCTACCCTTCCGAATTTGCCACATACCTGAATTTCTGCCGTTC CTTGCGTTTTGACGACAAGCCTGACTACTCGTACCTGCGGCAGCTTTTCCGGAATCTGTTCCATCGCCAG GGCTTCTCCTATGACTACGTGTTCGACTGGAACATGCTCAAATTTGGTGCCAGCCGGGCCGCCGATGACG CCGAGCGGGAGCGCAGGGACCGAGAGGAGCGGCTGAGACACTCGCGGAACCCGGCTACCCGCGGCCTCCC TTCCACAGCCTCCGGCCGCCTGCGGGGGACGCAGGAAGTGGCTCCCCCCACACCCCTCACCCCACACCTCA ACCGCGGGGCCCCCGTCAACATCTCCTCGTCCGACCTCACAGGCCGACAAGATACCTCTCGCATGTCCAC CTCACAGATTCCTGGTCGGGTGGCTTCCAGTGGTCTTCAGTCTGTCGTGCACCGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA





Protein Sequence: >RG200482 representing NM_001893

Red=Cloning site Green=Tags(s)


MELRVGNRYRLGRKIGSGSFGDIYLGTDIAAGEEVAIKLECVKTKHPQLHIESKIYKMMQGGVGIPTIRW CGAEGDYNVMVMELLGPSLEDLFNFCSRKFSLKTVLLLADQMISRIEYIHSKNFIHRDVKPDNFLMGLGK KGNLVYIIDFGLAKKYRDARTHQHIPYRENKNLTGTARYASINTHLGIEQSRRDDLESLGYVLMYFNLGS LPWQGLKAATKRQKYERISEKKMSTPIEVLCKGYPSEFATYLNFCRSLRFDDKPDYSYLRQLFRNLFHRQ GFSYDYVFDWNMLKFGASRAADDAERERRDREERLRHSRNPATRGLPSTASGRLRGTQEVAPPTPLTPTS HTANTSPRPVSGMERERKVSMRLHRGAPVNISSSDLTGRQDTSRMSTSQIPGRVASSGLQSVVHR


TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:





ACCN: NM_001893

ORF Size: 1245 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts

of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at customercom or by

calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Casein Kinase 1 delta (CSNK1D) (NM_001893) Human Tagged ORF Clone - RG200482

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001893.2</u>

 RefSeq Size:
 2030 bp

 RefSeq ORF:
 1248 bp

 Locus ID:
 1453

 UniProt ID:
 P48730

 Cytogenetics:
 17q25.3

Domains: pkinase, TyrKc, S_TKc

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Circadian rhythm - mammal, Gap junction, Hedgehog signaling pathway


Gene Summary: This gene is a member of the casein kinase I (CKI) gene family whose members have been

implicated in the control of cytoplasmic and nuclear processes, including DNA replication and repair. The encoded protein may also be involved in the regulation of apoptosis, circadian rhythm, microtubule dynamics, chromosome segregation, and p53-mediated effects on growth. The encoded protein is highly similar to the mouse and rat CK1 delta homologs. Three transcript variants encoding different isoforms have been found for this gene.

[provided by RefSeq, Feb 2014]



Product images:



Circular map for RG200482