

Product datasheet for **RG200485**

CPT1A (NM_001031847) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CPT1A (NM_001031847) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CPT1A
Synonyms:	CPT1; CPT1-L; L-CPT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG200485 representing NM_001031847
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGAAGCTACCAAGCTGTGGCCTTTCAGTTCACGGTCACTCCGGACGGGATTGACCTGCGGCTGA
 GCCATGAAGCTCTTAGACAAATCTATCTCTGGACTTATTCCCTGGAAAAAGAAGTTCATCAGATTCAA
 GAACGGCATCATCACTGGCGTGTACCCGCAAGCCCTCCAGTTGGCTTATCGTGGTGGTGGCGTGATG
 ACAACGATGTACGCCAAGATCGACCCCTCGTTAGGAATAATTGCAAAAATCAATCGGACTCTGAAACGG
 CCAACTGCATGTCCAGCCAGACGAAGAAGCTGGTCAAGCGCGTGTGTTGGCACCGGCCTGTGGGTGGC
 CCTCATCGTCACCATGCGCTACTCCCTGAAAGTGTGCTCTCTACCACGGGTGGATGTTCACTGAGCAC
 GGCAAGATGAGTCGTGCCACCAAGATCTGGATGGGTATGGTCAAGATCTTTTCAGGCCGAAAACCCATGT
 TGTACAGCTTCCAGACATCGCTGCCTCGCTGCCGGTCCCGGCTGTCAAAGACACTGTGAACAGGTATCT
 ACAGTCGGTGAGGCTCTTATGAAGGAAGAAGACTTCAAACGGATGACAGCACTTGTCAAGATTTTGTCT
 GTCGGTCTTGACCAAGATTACAGTGGTATTTGAAGTTAAAAATCCCTGGTGGGCTACAAAATACGTGAGCG
 ACTGGTGGGAGGAGTACATCTACCTCCGAGGACGAGGGCCGCTCATGGTGAACAGCAACTATTATGCCAT
 GGATCTGCTGTATATCTTCCAACCTCACATTACGGCAGCAAGAGCCGGCAACGCCATCCATGCCATCCTG
 CTTTACAGGCGCAAACCTGGACCGGGAGGAAATCAAACCAATTCGTCTTTTGGGATCCACGATCCACTCT
 GCTCCGCTCAGTGGGAGCGGATGTTTAATACTTCCCGGATCCCAGGAGAGGAGACAGACACCATCCAGCA
 CATGAGAGACAGCAAGCACATCGTCGTGTACCATCGAGGACGCTACTTCAAGGTCTGGCTCTACCATGAT
 GGGCGGCTGTGAAGCCCCGGGAGATGGAGCAGCAGATGCAGAGGATCCTGGACAATACCTCGGAGCCTC
 AGCCCCGGGAGGCCAGGCTGGCAGCCCTCACCGCAGGAGACAGAGTTCCCTGGCCAGGTGTCTCAGGC
 CTATTTTGGACGTGGGAAAAATAAGCAGTCTCTTGATGCTGTGGAGAAAGCAGCGTTCTTCTGTGAGTTA
 GATGAAACTGAAGAAGGATACAGAAGTGAAGACCCGGATACGTCAATGGACAGCTACGCCAAATCTCTAC
 TACACGGCCGATGTTACGACAGGTGGTTTGACAAGTCTGTTACGTTTGTGTTCTTCAAAAACGGGAAGAT
 GGGCCTCAACGCTGAACACTCCTGGGCAGATGCGCCGATCGTGGCCACCTTTGGGAGTACGTATGTCC
 ATTGACAGCCTCCAGCTGGGCTATGCGGAGGATGGGCACTGCAAAGGCGACATCAATCCGAACATCCGT
 ACCCCACCAGGCTGCAGTGGGACATCCCGGGGGAATGTCAAGAGGTTATAGAGACCTCCCTGAACACCGC
 AAATCTTCTGGCAAACGACGTGGATTTCCATTCTTCCATTCTGAGCCTTTGGTAAAGGAATCATCAAG
 AAATGTGCGCAGGAGCCAGACGCTTTGTGACGCTGGCCCTCCAGCTGGCGCACTACAAGGACATGGGCA
 AGTTTTGCCTCACATACGAGGCCTCCATGACCCGGCTTCCGAGAGGGGAGGACGGAGACCGTGCCTC
 CTGCAACACTGAGTCATGCGACTTCGTGCGGGCCATGGTGGACCCGGCCAGACGGTGAACAGAGGCTG
 AAGTTGTTCAAGTTGGCGTCTGAGAAGCATCAGCATATGTATCGCCTCGCCATGACCCGGCTCTGGGATCG
 ATCGTCACCTCTTCTGCCTTACGTGGTGTCTAAATATCTCGCTGTGGAGTCCCCCTTCTTAAGGAAGT
 TTTATCTGAGCCTTGGAGATTATCAACAAGCCAGACCCCTCAGCAGCAAGTGGAGCTGTTTGACTTGGAG
 AATAACCCAGAGTACGTGTCCAGCGGAGGGGCTTTGGACCGGTGTGTGATGACGGCTATGGTGTGTCGT
 ACATCCTTGTGGGAGAGAACCTCATCAATTTCCACATTTCTTCCAAGTTCTTGGCCTGAGACGGGGAT
 TATAAGTCAAGGACCAAGTTCAGATACT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG200485 representing NM_001031847
 Red=Cloning site Green=Tags(s)

MAEAHQAVAFQFTVTPDGIDLRLSHEALRQIYLSGLHSWKKKFIKFKNGIITGVYPASPSSWLIVVVGVM
 TTYAKIDPSLGIIAKINRTLETANCMSSQTKNVVSGVLFGTGLWVALIVTMRYSKLVLLSYHGWMFTEH
 GKMSRATKIWMGMVKIFSGRKPLYSFQTSPLRPLVPAVKDVTNRYLQSVRPLMKKEEDFKRMTALAQDFA
 VGLGPRLQWYLKLSWWATNYVSDWWEYIYLRGRGPLMVNSNYAMDLLYILPHTIQAARAGNAIHAIL
 LYRRKLDREEIKPIRLLGSTIPLCSAQWERMFNTRIPGEETDTIQHMRDSKHIVVYHRGRYFKVWLYHD
 GRLLKPREMEQQMQRILDNTSEPPGEARLAALTAGDRVPWARCRQAYFGRGKNKQSLDAVEKAFFVTL
 DETEEGYRSEDPDTSMDSYAKSLLHGRCYDRWFDKSFYVVFKNKMGMLNAEHSWADAPIVAHLWEYVMS
 IDSLQLGYAEDGHCKGDINPNIPYPTRLQWDIPGECQEVETSLNTANLLANDVDFHSFPFVAFGKGIK
 KCRTSPDAFVQLALQLAHYKDMGKFCLEYEASMTRLFREGRTETVRSCTTESCDFVRAMVDAQTVEQRL
 KLFKLASEKHQHYRLAMTGSIDRHLFCLYVVS KYLAVESPFLEKVLSEPWRLSTSQTPTQQVELFDLE
 NNPEYVSSGGGFGPVADDGYGVSYILVGENLINFHISSEKFCPETGIISQGPSSDT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001031847

ORF Size: 2268 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001031847.2
RefSeq Size:	2546 bp
RefSeq ORF:	2271 bp
Locus ID:	1374
UniProt ID:	P50416
Cytogenetics:	11q13.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Adipocytokine signaling pathway, Fatty acid metabolism, PPAR signaling pathway
Gene Summary:	<p>The mitochondrial oxidation of long-chain fatty acids is initiated by the sequential action of carnitine palmitoyltransferase I (which is located in the outer membrane and is detergent-labile) and carnitine palmitoyltransferase II (which is located in the inner membrane and is detergent-stable), together with a carnitine-acylcarnitine translocase. CPT I is the key enzyme in the carnitine-dependent transport across the mitochondrial inner membrane and its deficiency results in a decreased rate of fatty acid beta-oxidation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p>

Product images:



Circular map for RG200485