

Product datasheet for **RG200506**

AKAP1 (NM_003488) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AKAP1 (NM_003488) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AKAP1
Synonyms:	AKAP; AKAP84; AKAP121; AKAP149; D-AKAP1; PPP1R43; PRKA1; SAKAP84; TDRD17
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG200506 representing NM_003488
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGCAATCCAGTTCGTTTCGCTCTTCCCCTTGGCATTGCCTGGGATGCTGGCGCTCCTCGGCTGGTGGT
 GGTTTTTCTCTCGTAAAAAGGCCATGTCAGCAGCCATGATGAGCAGCAGGTGGAGGCTGGTGTGTGCA
 GCTGAGGGCTGACCTGCCATCAAGGAACCTCTCCCCGTGGAAGACGTCTGTCCCAAAGTAGTGTCCACA
 CCCCCAGTGTACAGAGCCTCCAGAAAAGAACTGTCCACCGTGAGCAAGCTGCCTGCAGAGCCCCAG
 CATTGCTCCAGACACCCACCTTGCCGAAGATCAGAGTCTCGGGCATTCTTCTAACACCACAGACAT
 GAGATTGCGACCAGAACACGCAGAGATGACAGTACAAAGCTGGAGCTAGCCCTGACAGGTGGTGAAGCC
 AAATCGATTCTCTAGAGTGCCCCCTTTCATCCCCAAAGGGTGTACTATTCTCCAGCAAATCAGCTGAGG
 TGTGTAAGCAAGATTCCCCCTTTCAGCAGGGTGCCAAGGAAGGTCCAGCCAGGCTACCCCGTAGTCCCCGC
 AGAGAAGCGTAGCTCTGGGGAGAGGGCAAGAGAGACAGGTGGGGCCGAAGGGACTGGTGTGCGCGTTG
 GGGGAAAAGGTGCTTGAAGAAGCTCTGTTGTCTCGGGAGCATGTCTTGAATTGGAGAACAGCAAGGGCC
 CCAGCCTGGCCTCTTTAGAGGGGGAAGAAGATAAGGGGAAGAGCAGCTCATCCCAGGTGGTGGGGCCAGT
 GCAGGAGGAAGAGTATGTAGCAGAGAAGTTGCCAAGTAGGTTTCATCGAGTCGGCTCACACAGAGCTGGCA
 AAGGACGATGCGGCGCCAGCACCCCCAGTCGCAGACGCCAAAGCCAGGATAGAGGTGTCGAGGGAGAAC
 TGGGCAATGAGGAGAGCTTGGATAGAAATGAGGAGGGCTTGGATAGAAATGAGGAGGGCTTGGATAGAAA
 TGAGGAGAGCTTGGATAGAAATGAGGAGGGCTTGGATAGAAATGAGGAGATTAAAGCGGGCTGCCTTCCAG
 ATAATCTCCCAAGTGATCTCAGAAGCAACCGAACAGGTGCTGGCCACCACGGTTGGCAAGGTTGCAGGTG
 GTGTGTGTCAGGCCAGTCAGCTCCAAGGGCAGAAGGAAGAGAGCTGTGTCCAGTTTACCAGAAAAGTGT
 CTTGGGGCCAGACACTGCGGAGCCTGCCACAGCAGAGGCAGCTGTTGCCCGCCGGATGCTGGCCTCCCC
 TTGCCAGGCCTACCAGCAGAGGGCTCACCACCACCAAAGACCTACGTGAGCTGCCTGAAGAGCCTTCTGT
 CCAGCCCCACCAAGGACAGTAAGCCAAATATCTCTGCACACCACATCTCCCTGGCCTCCTGCCTGGCACT
 GACCACCCCCAGTGAAGAGTTGCCGACCGGGCAGGCATCCTGGTGAAGATGCCACCTGTGTACCTGC
 ATGTCAGACAGCAGCCAAAGTGTCCCTTGGTGGCTTCTCCAGGACACTGCTCAGATTCTTTCAGCACTT
 CAGGGCTTGAAGACTCTTGCACAGAGACCAGCTCGAGCCCCAGGACAAGGCCATACCCCGCCACTGCC
 AGAAAGTACTGTGCCCTTCAAGCAATGGGGTGTGAAGGGGGAGTTGTGAGACTTGGGGGCTGAGGATGGA
 TGGACCATGGATGCGGAAGCAGATCATTAGGAGGTTCTGACAGGAACAGCATGGATTCCGTGGATAGCT
 GTTGACAGTCTCAAGAAGACTGAGAGCTTCCAAAATGCCCAGGCAGGCTCCAACCTAAGAAGGTCGACCT
 CATCATCTGGGAGATCGAGGTGCCAAAGCACTTAGTCGGTCGGCTAATTGGCAAGCAGGGGCGCTATGTG
 AGTTTTCTGAAGCAAACATCTGGTGCCAAGATCTACATTTCAACCCTGCCTTACACCCAGAGCGTCCAGA
 TCTGCCACATAGAAGGCTCTCAACATCATGTAGACAAAGCGCTGAACCTGATTGGGAAGAAGTTCAAAGA
 GCTGAACCTCACCAATATCTACGCTCCCCATTGCCTTCACTGGCACTGCCTTCTCTGCCGATGACATCC
 TGGCTCATGCTGCCTGATGGCATACCGTGGAGGTCATTGTGGTCAACCAGGTCAATGCCGGGCACCTGT
 TCGTGCAGCAGCACACACCCCTACCTTCCACGCGCTGCGCAGCCTCGACCAGCAGATGTACCTCTGTTA
 CTCTCAGCCTGGAATCCCACCTTGCCACCCAGTGGAAATAACGGTCATCTGTGCCGCCCTGGTGGC
 GACGGGGCTGGTGGCGAGCCCAAGTGGTGCCTCTACGAGGAGACCAACGAAGTGGAGATTTCGATACG
 TGCACTACGGCGGATATAAGAGGGTGAAAGTAGACGTGCTCCGGCAAATCAGGTCTGACTTTGTACCCCT
 GCCGTTTCAGGGAGCAGAAGTCCTTCTGGACAGTGTGATGCCCTGTGAGACGATGACCAGTTTTACCG
 GAAGCAGATGCCGCCATGAGCGAGATGACGGGAATACAGCACTGCTTGCTCAGGTGACAAGTTACAGTC
 CAACTGGTCTTCTCTGATTAGCTGTGGAGTGTGGTTGGAGATGAAGTGGTGTGATAAACCGGTCCCT
 GGTGGAGCGAGGCCTTGCCAGTGGGTAGACAGCTACTACACAAGCCTT

ACGCGTACGCGGCGGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG200506 representing NM_003488
Red=Cloning site Green=Tags(s)


[illegible]

TRTRPLE - GFP Tag - V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:



ACCN:	NM_003488
ORF Size:	2709 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003488.3 , NP_003479.1
RefSeq Size:	3968 bp
RefSeq ORF:	2712 bp
Locus ID:	8165
UniProt ID:	Q92667
Cytogenetics:	17q22
Domains:	TUDOR, KH, TUDOR
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	The A-kinase anchor proteins (AKAPs) are a group of structurally diverse proteins, which have the common function of binding to the regulatory subunit of protein kinase A (PKA) and confining the holoenzyme to discrete locations within the cell. This gene encodes a member of the AKAP family. The encoded protein binds to type I and type II regulatory subunits of PKA and anchors them to the mitochondrion. This protein is speculated to be involved in the cAMP-dependent signal transduction pathway and in directing RNA to a specific cellular compartment. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG200506