

Product datasheet for **RG200540**

RIMS3 (NM_014747) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	RIMS3
Synonyms:	NIM3; RIM 3; RIM3
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG200540 representing NM_014747
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGTTTAACGGGGAGCCAGGTCCTGCCTCATCTGGGGCTCCAGGAATGTGGTGCGGAGCTCCAGCATT
GCGGTGAAATCTGCGGATCCCAGCAAGCCGGGGCGGGCTGGGACCACCACCGCAAGAAGCGGCGGAG
CAGCCTGGGTGCCAAGATGGTGGCCATCGTGGGCTGACTCAGTGGAGCAAGAGCACACTCCAGCTTCCG
CAGCCTGAAGGGGCCACCAAGAAGCTGCGCAGCAACATCCGCCGGAGCAGGAGACAGGCATCGCGGTGG
AGATGCGGAGCCGGGTACACGCCAGGGCAGCCGGGAGTCCACCGATGGGAGCACCAACAGCAACAGCTC
CGACGGCAGCTTCATCTTCCCCACTACCCGGCTAGGGGCTGAAAGCCAGTTTCAAGCATTCTGGATGGG
CTGGGACCAGCTCAGATTGTGGGGCAGACACTGGCAACACCACCCATGGGAGATGTGCACATTGCCA
TCATGGACCGGAGTGGCCAGCTGGAGGTGGAAGTGATTGAAGCTCGGGGCTGACCCCAAAACAGGCTC
CAAATCCCTCCCAGCCACCTATATCAAGGTTTACCTGCTGGAGAATGGGGCCTGCTTGGCCAAGAAGAAG
ACAAAGATGACCAAGAAGACCTGTGATCCCTGTACCAGCAGGCTCTGCTCTTTGACGAGGGACCCAGG
GCAAGGTGCTGCAGGTGATCGTCTGGGGAGACTATGGCCGCATGGACCACAAGTGCTTCATGGGCATGGC
CCAGATCATGCTGGACGAGCTGGACCTCAGCGCCGCGGTACCGGCTGGTACAACTCTTCCCCACCTCC
TCAGTGGCAGACTCCACACTCGGATCCCTACCAGGCGCTGTCCAGTCTTCCCTGGAGAGTGCCACCA
GCCCTCATGCTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG200540 representing NM_014747
Red=Cloning site Green=Tags(s)



MFNGEPGPASSGASRNVVRSSISGEICGSQQAGGAGTTTAKKRRSSLGAKMVAIVGLTQWSKSTLQLP
QPEGATKKLRNIRRTSTETGIAVEMRSRVTRQGSRESTDGTNSNSDGTFFIPTTRLGAESQFSDFLDG
LGPAQIVGRQTLATPPMGDVHIAIMDRSGQLEVEVIEARGLTPKPGSKSLPATYIKVYLLENGACLAKKK
TKMTKKTCDFLYQQALLFDEGPQGVKLVQIVWGDYGRMDHKCFMGMAQIMLDELDSLAAVTGWYKLFPTS
SVADSTLGSLTRRLSQSSLESATSPSCS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_014747

ORF Size: 924 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_014747.3](#)

RefSeq Size: 7233 bp

RefSeq ORF: 927 bp


Locus ID: 9783

UniProt ID: [Q9UJD0](#)

Cytogenetics: 1p34.2

Gene Summary: Regulates synaptic membrane exocytosis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG200540