

Product datasheet for **RG200561**

FBX09 (NM_033481) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | FBX09 (NM_033481) Human Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | FBX09 |
| Synonyms: | dj341E18.2; FBX9; NY-REN-57; VCIA1 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >RG200561 representing NM_033481 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTCCGAGCTCAGTGGATGTTTGAAGTCTGCTCCAGGTGTAAGCTCTAGCAATTTAGAAAATCGACCTT
GCAGAGCAGCAAGAGGCTCTCTCCAGAAAACATCGGCAGATACCAAAGGAAAACAAGAACAGGCAAAAAGA
AGAAAAGGCTCGAGAAGCTCTTCTAAAAGCAGTAGAAGAACAACAAAATGGAGCTCTCTATGAAGCCATC
AAGTTTTATCGTAGGGCTATGCAACTTGTACCTGATATAGAGTTCAAGTACTTATAACCGGTCTCCAG
ATGGTGTGGCGTTGAAACAGCTACATTGAAGATAATGATGATGACAGCAAAATGGCAGATCTCTTGTC
CTACTTCCAGCAGCAACTCACATTTTCAGGAGTCTGTGCTTAAACTGTGTCAGCCTGAGCTTGAGAGCAGT
CAGATTCACATATCAGTCTGCCAATGGAGTCTGATGTACATCTTCCGATGGGTGGTGTCTAGTGACT
TGGACCTCAGATCATTGGAGCAGTTGTCGCTGGTGTGTCAGAGGATTCTACATCTGTGCCAGAGACCTGA
AATATGGCGTCTGGCCTGCTTGAAGTTTGGGGCAGAAGCTGTATTAACCTGTTCCGTACACGTCCTGG
AGAGAGATGTTTTAGAACGGCTCGTGTTCGGTTTGTGGCGTGTATATCAGTAAAACCATATATTC
GTCAAGGGGAACAGTCTCTTGTGGTTTCTATAGAGCTGGCACCAGTGAATATTACAGGTACATAAG
ATCTTTTCTGATGGCCATGTGATGATGTTGACAACCCCTGAAGAGCCTCAGTCCATTGTTCCACGTTTA
AGAACTAGGAATACCAGGACTGATGCAATTTCTACTGGTCACTATCGCTTGTCAAGACACAGACAAT
AGACCAAAGTATTTGCTGTAATAACTAAGAAAAAGAAGAAAAACCACTTGACTATAAATACAGATATTT
TCGTGCTGTCCTGTACAAGAAGCAGATCAGAGTTTTTCATGTGGGGCTACAGCTATGTTCCAGTGGTCAC
CAGAGGTTCAACAACTCATCTGGATACATCATTCTTGTACATTACTTACAAATCAACTGGTGAGACTG
CAGTCAGTGCTTTTGGAGATTGACAAGATGTACACCCCTTGTCTTCGCCAGAGTAAGGAGCTACACAGC
TTTCTCAGAAAGGCTCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG200561 representing NM_033481
 Red=Cloning site Green=Tags(s)

MFRAQWMFELAPGVSSSNLENRPCRAARGSLQKTSADTKGKQEQAKEEKARELFLKAVEEEQNGALYEAI
 KFYRRAMQLVPDIEFKITYTRSPDGDGVGNSYIEDNDDSKMADLLSYFQQQLTFQESVLKLCQPELESS
 QIHISVLPMEVLMYIFRWVSSDLDLRSLEQLSLVCRGFYICARDPEIWRACLKLVWGRSCIKLVPYTSW
 REMFLERPRVRFDGVIYSKTTYIRQGEQSLDGFYRAWHQVEYYRYIRFFPDGHVMMLTTPPEEPQSI V PRL
 RTRNTRTDAILLGHYRLSQDTDNQTQVFVAVITKKKEEKPLDYKYRYFRRVVPVQEADQSFHVGLQLCSSGH
 QRFNKLIIWHHSCHITYKSTGETAVSAFEIDKMYTPLFFARVRSYTAFSERPL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_033481

ORF Size: 1209 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_033481.3](#), [NP_258442.2](#)

RefSeq Size: 4565 bp

RefSeq ORF: 1212 bp

Locus ID: 26268

UniProt ID: [Q9UK97](#)


Cytogenetics: 6p12.1

Domains: F-box

Protein Families: Druggable Genome

Gene Summary: This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class. Alternative splicing of this gene generates at least 3 transcript variants diverging at the 5' terminus. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG200561