

Product datasheet for RG200708

Myosin (MYL6) (NM_021019) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids Product Name: Myosin (MYL6) (NM_021019) Human Tagged ORF Clone Tag: TurboGFP Symbol: Myosin Synonyms: ESMLC; LC17; LC17-GI; LC17-NM; LC17A; LC17B; MLC-3; MLC1SM; MLC3NM; MLC3SM Mammalian Cell Neomycin Selection: Vector: pCMV6-AC-GFP (PS100010) E. coli Selection: Ampicillin (100 ug/mL) **ORF** Nucleotide >RG200708 representing NM_021019 Red=Cloning site Blue=ORF Green=Tags(s) Sequence: TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCGCGATCGCC ATGTGTGACTTCACCGAAGACCAGACCGCAGAGTTCAAGGAGGCCTTCCAGCTGTTTGACCGAACAGGTG ATGGCAAGATCCTGTACAGCCAGTGTGGGGGATGTGATGAGGGCCCTGGGCCAGAACCCTACCAACGCCGA GGTGCTCAAGGTCCTGGGGAACCCCCAAGAGTGATGAGATGAATGTGAAGGTGCTGGACTTTGAGCACTTT CTGCCCATGCTGCAGACAGTGGCCAAGAACAAGGACCAGGGCACCTATGAGGATTATGTCGAAGGACTTC GGGTGTTTGACAAGGAAGGAAATGGCACCGTCATGGGTGCTGAAATCCGGCATGTTCTTGTCACACTGGG TGAGAAGATGACAGAGGAAGAAGTAGAGATGCTGGTGGCAGGGCATGAGGACAGCAATGGTTGTATCAAC TATGAAGCGTTTGTGAGGCATATCCTGTCGGGG ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA >RG200708 representing NM_021019 **Protein Sequence:** Red=Cloning site Green=Tags(s) MCDFTEDQTAEFKEAFQLFDRTGDGKILYSQCGDVMRALGQNPTNAEVLKVLGNPKSDEMNVKVLDFEHF LPMLQTVAKNKDQGTYEDYVEGLRVFDKEGNGTVMGAEIRHVLVTLGEKMTEEEVEMLVAGHEDSNGCIN YEAFVRHILSG TRTRPLE - GFP Tag - V **Restriction Sites:** Sgfl-Mlul




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Cloning Scheme:




ACCN:	NM_021019
ORF Size:	453 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 021019.5</u>
RefSeq Size:	704 bp
RefSeq ORF:	456 bp
Locus ID:	4637
UniProt ID:	<u>P60660</u>

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

ORIGENE Myosin (MYL6) (NM_021019) Human Tagged ORF Clone – RG200708		
Cytogenetics:	12q13.2	
Domains:	EFh	
Protein Families:	Druggable Genome	
Protein Pathways:	Vascular smooth muscle contraction	
Gene Summary:	Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain that is expressed in smooth muscle and non- muscle tissues. Genomic sequences representing several pseudogenes have been described and two transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]	

Product images:



Circular map for RG200708

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US