

Product datasheet for **RG200754**

ACBD6 (NM_032360) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ACBD6
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG200754 representing NM_032360
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCTTCATCATTCTGCCCAGGGGGCCATACCGGCGACAGCGGTGGAGAGCTGAGCTCAGGGGACG
ACTCCGGGGAGGTGGAGTTCCCCATAGCCCTGAGATCGAGGAGACCAGTTGCCTGGCCGAGCTGTTTGA
GAAGGCTGCCGCGCACCTGCAAGGCCTGATTCAAGTGCCAGCAGGGAGCAGCTCTTGTACCTGTATGCC
AGGTACAAACAGGTCAAAGTTGAAATTGTAATACTCCTAAACCAAGCTTCTTTGATTTGAAGGAAAGC
AAAAATGGGAAGCTTGGAAAGCACTTGGTGATTCAAGCCCCAGCCAAGCAATGCAGGAATATATCGCAGT
AGTTAAAAAAGTAGATCCAGGTTGGAATCCTCAGATACCAGAGAAGAAAGGAAAGAAAGCAATACAGGT
TTTGGTGGGCCAGTTATTAGTTCTCTATATCATGAAGAAACCATCAGGGAAGAAGACAAAAATATATTTG
ATTACTGCAGGGAAAAACATTGACCATATAACCAAAGCCATCAAATCGAAAAATGTGGATGTGAATGT
GAAAGATGAAGAGGAGTAGGGCTCTACTTCACTGGGCTGTGATCGAGGACATAAGGAACTAGTCACAGTG
TTGCTGCAACATAGAGCTGACATTAAGTGCAGGACAATGAAGGCCAAACAGCTCTACATTATGCCTCTG
CCTGTGAGTTTCTGGATATTGTAGAGCTGCTGCTCCAGTCTGGTGTGACCCCACTCTCCGAGACCAGGA
TGGCTGCCTGCCAGAGGAGGTGACAGGCTGCAAAACAGTTTCTTTGGTGTGACGCGGCACACAAGTGGC
AAGGCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG200754 representing NM_032360
Red=Cloning site Green=Tags(s)


MASSFLPAGAITGDSGGELSSGDDSGEVEFPHSP EIEETSCLAELFEKAAAHLOGLIQVASREQLLYLYA
RYKQVKVGNCNTPKPSFFDFEGKQKWEAWKALGDSSPSQAMQEYIAVVKKLDPGWNPQIPEKKGKEANTG
FGGPVISSLYHEETI REEDKNIFDYCRENNIDHITKAISKNV DVNVKDEEGRALLHWACDRGHKELVTV
LLQHRADINCQDNEGQTALHYASACEFLDIVELL LQSGADPTLRDQDGCLPEEVTGCKTVSLVLQRHTTG
KA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



CTATAGGGCGGCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGGCGCGCAGATCT

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI

CAAGCTTAAGCTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC --- --- ---

HindIII NheI RsrII MluI NotI XhoI GFP Tag

T R T R P L E M E S D - - -

--- GAA GAA AGA GTT TAA ACGGCCGCCCGGAGCT

- - E E R V Stop

ACCN: NM_032360

ORF Size: 846 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_032360.4</u>
RefSeq Size:	1089 bp
RefSeq ORF:	849 bp
Locus ID:	84320
UniProt ID:	<u>Q9BR61</u>
Cytogenetics:	1q25.2-q25.3
Domains:	ACBP, ANK
Gene Summary:	Binds long-chain acyl-coenzyme A molecules with a strong preference for unsaturated C18:1-CoA, lower affinity for unsaturated C20:4-CoA, and very weak affinity for saturated C16:0-CoA. Does not bind fatty acids.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG200754