

Product datasheet for **RG200761**

HIC5 (TGFB111) (NM_015927) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HIC5 (TGFB111) (NM_015927) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	HIC5
Synonyms:	ARA55; HIC-5; HIC5; TSC-5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG200761 representing NM_015927
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCAAGGTCAGGGCTCCCAAAGAGCGCCCTGCGGAGCCTCTCACCCCTCCCCATCCTATGGCCACC
 AGCCACAGACAGGTCTGGGAGTCTTCAGGAGCCTCGGGGACAAGACCACCTGTACAGCACGGTATG
 CAAGCCTCGGTCCCAAAGCCTGCAGCCCGCGGCCCTCCATTCTCCTCTCCAGCGGTGCTTGGGT
 ACCGGGCTCTGTGAGCTAGATCGGTTGCTTCAGGAACCTAATGCCACTCAGTTCAACATCACAGATGAAA
 TCATGTCTCAGTCCCATCTAGCAAGGTGGCTTCAGGAGAGCAGAAGGAGGACCAGTCTGAAGATAAGAA
 AAGACCCAGCCTCCCTCCAGCCCGTCTCCTGGCCTCCCAAAGGCTTCTGCCACCTCAGCCACTCTGGAG
 CTGGATAGACTGATGGCCTCACTCTGACTTCCGCGTTCAAACCATCTCCAGCCTCTGGGCAACTC
 AGCCACCGGTGGTGAAGTCCACAATGAGGGCTCCCCATCCCACCAGAGCCGACTGGCAAGGGCAGCCT
 AGACACCATGCTGGGGCTGCTGCAGTCCGACCTCAGCCGCGGGGTGTTCCACCAGGCCAAAGGCTC
 TGTGGCTCCTGCAATAAACCTATTGCTGGGCAAGTGGTGACGGCTCTGGGCCGCGCTGGCACCCCGAGC
 ACTTCGTTTTCGGGAGGCTGTTCCACCGCCTGGGAGGCAGCAGCTTCTCGAGAAGGATGGAGCCCCCTT
 CTGCCCGAGTGTACTTTGAGCGCTTCTCGCAAAGATGTGGCTTCTGCAACCAGCCATCCGACACAAG
 ATGGTGACCGCTTGGGCACTACTGGCACCCAGAGCATTTCTGCTGCGTCAGTTGCGGGAGCCCTTCG
 GAGATGAGGGTTTCCACGAGCGCGAGGGCCGCCCTACTGCCCGCGGACTTCTGCAGCTGTTGCCCC
 GCGCTGCCAGGGTCCAGGGCCCCATCCTGGATAACTACATCTCGGCGCTCAGCGCGCTCTGGCACCCG
 GACTGTTTCGTCTGCAGGAATGCTTCGCGCCCTCTCGGGAGGCAGCTTTTTCGAGCACGAGGGCCGCC
 CGTTGTGCGAGAACCCTCCACGCACGACGCGGCTCGCTGTGCGCCACGTGTGGCCTCCCTGTGACCCG
 CCGCTGCGTGTCCGCCCTGGGTGCGCGCTTCCACCCGGACCCTTACATGCACCTTCTGCTGCGCCCG
 CTCACCAAGGGTCTTCCAGGAGCGCGCCGCAAGCCCTACTGCCAGCCCTGTTCTGAACTCTTCG
 GC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG200761 representing NM_015927
 Red=Cloning site Green=Tags(s)

MPRSGAPKERPAEPLTPPPSYGHQPQTGSGESSGASGDKDHLVSTVCKPRSPKPAAPAAPPFSSSSGVLG
 TGLCELDRLQLLNATQFNITDEIMSQFPSSKVASGEQKEDQSEDKRPSLPSSPSPGLPKASATSATLE
 LDRLMASLSDFRVQNHLPASGPTQPPVVSSTNEGSPSPPEPTGKGLDTMLGLLQSDLRRGVPTQAKGL
 CGSCNKP IAGQVVTALGRAWHPEHFVCGGCSTALGGSSFFEKDGAPFCPECYFERFSPRCGFCNQPIRHK
 MVTALGTHWHPEHFCCVSCGEPFGDEGFHEREGRPYCRRDFLQLFAPRCQGCQGPILDNYISALSALWHP
 DCFVCRECFAPFSGGSFFEHEGRPLCENHFHARRGSLCATCGLPVTGRCVSALGRRFHPDHFTCTFCLRP
 LTKGSFQERAGKPYCQPCFLKLF

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_015927

ORF Size: 1332 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_015927.4](#)

RefSeq Size: 1824 bp
 RefSeq ORF: 1335 bp
 Locus ID: 7041
 UniProt ID: [O43294](#)
 Cytogenetics: 16p11.2

Protein Families: Druggable Genome, Transcription Factors

Gene Summary: This gene encodes a coactivator of the androgen receptor, a transcription factor which is activated by androgen and has a key role in male sexual differentiation. The encoded protein is thought to regulate androgen receptor activity and may have a role to play in the treatment of prostate cancer. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]

Product images:



Circular map for RG200761