

Product datasheet for **RG200775**

ZSWIM1 (NM_080603) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZSWIM1
Synonyms:	C20orf162
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG200775 representing NM_080603
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCTTGAGAGACTCAAAGCCCCGTGGTCAGTGCCTGCAAAGAAAGTATTTGACCTTGGCATTGGA
 CAGCTCCCATCTCTCCCATGGCCCTGACAATGCTGAATGGGCTCCTGATTAAGGACTCAAGCCACCTAT
 GCTGCTGCACCAAGTTAAACAGACTGCCAGTTAGATACCTTCACTACCAAGAGCTGCTTTATGCAAAGT
 GTCTTTGACCATTTCCCTGAGATCTTATTTATCCACCGGACCTATAACCCAAGGGTAAGGTCTTATATA
 CCTTCTGGTGGATGGACCTCGGGTGCAGCTGGAGGGTCATCTTGCCCGAGCAGTCTACTTTGCCATCCC
 TGCCAAGGAGGACACTGAAGGCTGGCCCAAGATGTTCCAAGTATTCAAGAAGTTTAATCCAGCATGGGAG
 AGAGTCTGTACCATCCTGGTGGATCCTCATTTCTTCCACTGCCTATCCTAGCTATGGAGTTCCTCCACAG
 CTGAGGTCTTCTCTCAGCCTTCCACATTTGTAAGTTCCTCCAGGCCAAGTTCTATCAGCTGTCCCTTGA
 ACGGCCCGTGAAAGGCTGCTCCTGACCTCCCTGCAGAGCACAATGTGCTCAGCCACAGCAGGCAACCTG
 AGAAAGTTGTATACACTCCTGAGCAATGCATCCCTCCAGCCAAGCTGCCGAGCTTCACTCACACTGGC
 TGCTCAACGACCGCATCTGGCTGGCTCACCCTGGAGAAGCCGAGCTGAGAGCAGCCACTACTTCCAGAG
 CCTCGAGGTCAACCCACATCCTCAGCCAGTTCTTTGGTACCACCCCATCTGAGAAACAAGGTATGGCT
 TCTCTGTTCCGTTACATGCAGCAGAACTCTGCAGACAAGGCCAACTTCAACAGGGCCTGTGTGCCCAGA
 ACAATCATGCTCCCTCAGACACCATCCCCGAAAGCCCCAACTGGAGCAGCTGGTAGAATCCCACATCCA
 GCACTCCCTCAATGCCATCTGCACAGGGCCAGCAGCCCAACTGTGCCTGGGCGAGCTTGCTGTGGTCCAG
 AAATCCACACACCTCATTGGCTCTGGCTCAGAAAAGATGAACATACAGATCCTGGAAGATACCCATAAGG
 TGCAGCCCCAGCCCCCTGCCAGCTGCAGTGTCTACTTTAACCAGGCCTTCCACCTGCCCTGCCGCCACAT
 CCTAGCCATGCTCAGTGGCCGCGCCAGGTGCTCCAGCCGACATGCTGCCGGCTCAGTGGACGGCAGGC
 TGTGCTACCACTTAGACAGCATCCTGGGCAGCAAGTGGAGTGAGACCCTGGATAAGCACCTGGCAGTGA
 CTCACCTCACCGAGGAGGTGGGTGAGTGTGCAGCACTGCACCAAGGAGGAGTTTGAGCGGAGGTATAG
 CACCTGCGGGAAGTGGCCGACAGTGGATTGGCCTTATGAGCAGGTCCAATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG200775 representing NM_080603
 Red=Cloning site Green=Tags(s)

MLERLKAPWSAALQRKYFDLGIWTAPISPMALTMLNGLLIKDSSPPMLLHQVNKAQLDTFNYQSCFMQS
 VDFHFPEILFIHRTYNPRGKVLTYFLVDGPRVQLEGLARAVYFAIPAKEDTEGLAQMFQVFKFNPAWE
 RVCTILVDPHFLPLILAMEFPTAEVLLSAFHICKFLQAKFYQLSLERPVERLLLTSLQSTMCSATAGNL
 RKLYTLLSNCIPPAPKLPELHSHWLLNDRIWLAHRWRSRAESSHYFQSLEVTTHILSQFFGTTTPSEKQGM
 SLFRYMQQNSADKANFNQGLCAQNNHAPSDTIPESPKLEQLVESHIQHSLSNAICTGPAAQLCLGELAVVQ
 KSTHLIGSGSEKMNIQILEDTHKVQPPASCSCYFNQAFHLPCHILAMLARRQVLQPDMLPAQWTAG
 CATSLDSILGSKWSETLDKHLAVTHLTEEVGQLLQHCTKEEFERRYSTLRELADSWIGPYEQVQL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_080603

ORF Size: 1455 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_080603.5](#)

RefSeq Size: 2206 bp

RefSeq ORF: 1458 bp


Locus ID: 90204

UniProt ID: [Q9BR11](#)

Cytogenetics: 20q13.12

Domains: SWIM

Product images:



Circular map for RG200775