

Product datasheet for **RG200828**

KATNAL1 (NM_032116) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KATNAL1 (NM_032116) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	KATNAL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG200828 representing NM_032116
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAATTTGGCTGAGATTTGTGATAATGCAAAGAAAGGAAGAGAATATGCCCTTCTTGGAAATTACGACT
 CATCAATGGTATATTACCAGGGGTGATGCAGCAGATTACAGACATTGCCAGTCAGTCAGAGATCCAGC
 TATCAAAGGCAAATGGCAACAGGTTTCGGCAGGAATTATTGGAGGAATATGAACAAGTTAAAAAGTATTGTC
 AGCACTTTAGAAAAGTTTTAAAATTGACAAGCCTCCAGATTTCCCTGTGTCTGTCAAGATGAACCATTTA
 GAGATCCTGCTGTTTGGCCACCCCTGTTCTGCAGAACACAGAGCTCCACCTCAGATCAGGCGTCCCAA
 TCGAGAAGTAAGACCTCTGAGGAAAGAAATGGCAGGAGTAGGAGCCCGGGACCTGTAGGCCGAGCACAT
 CCTATATCAAAGAGTGAAAAGCCTTCTACAAGTAGGGACAAGGACTATAGAGCAAGAGGGAGAGATGACA
 AGGGAAGGAAGAATATGCAAGATGGTCAAGTGTGGTAAATGCCAAAATTTGATGGTGTGGTTATGA
 TAAGGATCTGGTGAAGCCCTTGAAGAGACATTGTATCCAGGAATCCTAGCATTATTGGGATGACATA
 GCAGATCTGGAAGAAGCTAAGAAGTTGCTAAGGGAAGCTGTTGTTCTTCCAATGTGGATGCCTGACTTTT
 TCAAAGGGATTAGAAGGCCATGGAAGGGTGTACTGATGGTTGGACCCCGAGCACTGGTAAAACATGCT
 AGCTAAAGCTGTTGCCACTGAATGTGGTACAACATTTCTCAACGTTTCGTCTTCTACACTGACATCTAAA
 TACAGAGGTGAATCTGAGAAGTTAGTTCGTCTGTTGTTTGGATGGCTAGATTTTATGCCCTACCACGA
 TCTTCATTGATGAGATAGATTCTATCTGCAGTCGAAGAGGAACCTCTGATGAACATGAGGCAAGTCGCAG
 GGTCAGTCTGAAGTCTCATTAGATGGATGGAGTTGGAGGAGCTTTAGAAAATGATGATCCTTCCAAA
 ATGGTTATGGTATTGGCTGCTACTAATTTCCCGTGGGACATTGATGAAGCTTTGCGAAGAAGGTTAGAAA
 AAAGGATATATACCTCTCCCAACAGCAAAGGAAGAGCTGAGCTTCTGAAGATCAACCTTCGTGAGGT
 CGAATTAGATCCTGATATTCAACTGGAAGATATAGCCGAGAAGATTGAGGGCTATTCTGGTGTGACATC
 ACTAATGTTTGCAGGGATGCCTCTTAAATGGCAATGAGACGGCGTATCAATGGCTTAAGTCCAGAAGAAA
 TCCGTGCACTTTCTAAAGAGGAACCTCAGATGCCTGTTACCAAAGGAGACTTTGAATTGGCCCTAAAGAA
 AATTGCTAAGTCTGTCTGCTGCAGACTTGGAGAAGTATGAAAATGGATGGTTGAATTTGGATCTGCT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG200828 representing NM_032116
 Red=Cloning site Green=Tags(s)

MNLAEICDNAKKGREYALLGNYDSSMVYQGVMMQIQRHQCQSVRDPAIKGGKQQVRQELLEEEYEQVKSIV
 STLESFKIDKPPDFPVSCQDEPFRDPAVWPPVPAEHRAPPQIRRPNREVRPLRKEMAGVGARGPVGRAH
 PISKSEKSTSRDKDYRARGRDDKGRKNMQDGASDGEMPKFDGAGYDKDLVEALERDIVSRNPSIHWDDI
 ADLEEAKLLREAVVLPMMWPDFFKGIRRPWKGVLVGGPPGTGKTM LAKAVATECGTTFFNVSSSTLTSK
 YRGESEKLVRLLFEMARFYAPTTIFIDEIDSICSRRTSDEHEASRRVKSELLIQMDGVGGALENDPDK
 MVMVLAATNFPWDIDEALRRRLEKRIYIPLPTAKGRAEL LKINLREVELDPDIQLEDIAEKIEGYSGADI
 TNVCRDASLMAMRRRINGLSPEIRALSKEELQMPVTKGFELALKKIAKSVSAADLEKYEKWMVEFGSA

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_032116

ORF Size: 1470 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032116.5](#)

RefSeq Size: 7536 bp

RefSeq ORF: 1473 bp

Locus ID: 84056


UniProt ID: [Q9BW62](#)

Cytogenetics: 13q12.3

Domains: AAA, AAA

Gene Summary: Regulates microtubule dynamics in Sertoli cells, a process that is essential for spermiogenesis and male fertility. Severs microtubules in an ATP-dependent manner, promoting rapid reorganization of cellular microtubule arrays (By similarity). Has microtubule-severing activity in vitro (PubMed:26929214).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG200828