

Product datasheet for **RG200875**

C1orf43 (NM_015449) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C1orf43 (NM_015449) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C1orf43
Synonyms:	HSPC012; NICE-3; NICE3; NS5ATP4; S863-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG200875 representing NM_015449 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGTCCGGCAGTAACTGGCTCTCCGGGTGAATGTCGTGCTGGTGATGGCCTACGGGAGCCTGGACT
TGAAAGAGGAGATTGATATTCGACTCTCCAGGGTTCAGGATATCAAGTATGAGCCCCAGCTCCTTGAGATG
TGATGATGCTAGACTACTACAACTGAAACCCAGGAAATCAAAGTTGCTACAATATCTGTATAGGATG
AAAGCTCTGGATGCCATTCGTACCTCTGAGATCCATTTTCATTCTGAAGCCGGCATCCCGTTCCTTAA
TGGGCAAGAATTTCCGCTCCTACCTGCTGGATCTGCGAAACACTAGTACGCCTTTCAAGGGTGTACGCAA
AGCACTCATTGATACCCTTTTGGATGGCTATGAAACAGCCCGCTATGGGACAGGGGTCTTTGGCCAGAAT
GAGTACCTACGCTATCAGGAGGCCCTGAGTGAGCTGGCCACTGCGGTTAAAGCACGAATTTGGGAGCTCTC
AGCGACATCACCAGTCAGCAGCCAAAGACCTAACTCAGTCCCCTGAGGTCTCCCCAACCAACCATCCAGGT
GACATACCTCCCCTCCAGTCAGAAGAGTAAACGTGCCAAGCACTTCTTGAATTGAAGAGCTTTAAGGAT
AACTATAACACATTGGAGAGTACTCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG200875 representing NM_015449
 Red=Cloning site Green=Tags(s)

MASGSNWLSGVNVVLVMAYGSLDLKEEIDIRLSRVQDIKYEPLLADDDARLLQLLETQGNQSCYNYLYRM
 KALDAIRTSEIPFHSEGRHPRSLMGKNFRSYLLDLRNTSTPFKGVKALIDTLLDGYETARYGTGVFGQN
 EYLRVQEALSELATAVKARIGSSQRHHQSAAKDLTQSPVSPPTTIQVTYLPSSQKSKRAKHFLKSFKD
 NYNTLESTL


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_015449

ORF Size: 657 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_015449.4
RefSeq Size:	1807 bp
RefSeq ORF:	660 bp
Locus ID:	25912
UniProt ID:	Q9BWL3
Cytogenetics:	1q21.3
Protein Families:	Transmembrane
Gene Summary:	General regulator of phagocytosis. Required to uptake Gram negative bacterium by macrophages.[UniProtKB/Swiss-Prot Function]