

Product datasheet for **RG200896**

TOR3A (NM_022371) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TOR3A (NM_022371) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TOR3A
Synonyms:	ADIR; ADIR2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG200896 representing NM_022371 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTTCGCGGTCCGTGGCGCCAGCTTTGGCTCTTTTTCCTGCTGCTGCTCCCGGGCGGCCTGAGCCCC
GCGGCGCTCCAGGCCGTGGGAGGGAACCGACGAGCCGGGCTCGGCCTGGGCTGGCCGGGCTTCCAGCG
CCTGCAGGAGCAGCTCAGGGCGGCGGTGCCCTCCAAGCGGTACTGGACGCTTTCAGCTGCCAGGTG
TGGCCCGACGACTGTGACGAGGACGAGGAGGCAGCCACGGGCCCCCTGGGCTGGCGCTTCTCTGTTGG
GCCAGCGTACCTGGACCTCCTGACCACGTGGTACTGCAGCTTCAAAGACTGCTGCCCTAGAGGGGATTG
CAGAATCTCCAACAACCTTACAGGCTTAGAGTGGGACCTGAATGTGCGGCTGCATGGCCAGCATTGGTC
CAGCAGCTGGTCCTAAGAACAGTGAGGGCTACTTAGAGACGCCCCAGCCAGAAAAGGCCCTTGCTCTGT
CGTTCCACGGCTGGTCTGGCACAGGCAAGAACCTCGTGGCACGGATGCTGGTGGAGAACCTGTATCGGGA
CGGGCTGATGAGTGACTGTGTCAGGATGTTTCATCGCCACGTTCCACTTTCCTCACCCAAATATGTGGAC
CTGTACAAGGAGCAGCTGATGAGCCAGATCCGGGAGACGACGAGCTCTGCCACCAGACCCTGTTTATCT
TCGATGAAGCGGAGAAGCTGCACCCAGGGCTGCTGGAGTCTTGGGCCACACTAGAACGCCGGGCCCC
TGAGGGCCACAGGGCTGAGTCTCCATGGACTATCTTTCTGTTTCTCAGTAATCTCAGGGGCGATATAATC
AATGAGGTGGTCTAAAGTTGCTCAAGGCTGGATGGTCCCGGGAAGAAATTACGATGGAACACCTGGGAC
CCACCTCCAGGCGGAGATTGTGGAGACCATAGACAATGGCTTTGGCCACAGCCGCTTGTGAAGGAAAA
CCTGATTGACTACTTCATCCCCTTCTGCTTTGGAGTACCGTACGTCAGTGAGGCTGTGTGCACGGGATGCC
TTCCTGAGCCAGGAGCTCCTGTATAAAGAAGAGACACTGGATGAAATAGCCAGATGATGGTGTATGTCC
CCAAGGAGGAACAACCTTTTCTCCAGGGCTGCAAGTCTATTTCCAGAGGATTAACACTTCTCTGTC
A

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >RG200896 representing NM_022371
 Red=Cloning site Green=Tags(s)


MLRGPWRQLWFLFLLLLPGAPEPRGASRPWEGTDEPGSAWAWPGFQRLQEQLRAAGALSKRYWTLFSCQV
 WPDDCDEDEEAATGPLGWRLPLLQRYLDLLTTWYCSFKDCCPRGDCRISNNFTGLEWDLNVRHLHGQHLV
 QQLVLRTRVRYLETPOPEKALALSFHGWSTGKNFVARMLENL YRDGLMSDCVRMFIATFHFPHPKYVD
 LYKEQLMSQIRETQQLCHQTLFIFDEAEKLPGLLEVLGPHLERRAPEGHRAESPWTIFLFLSNLRGDII
 NEVVLKLLKAGWSREEITMEHLEPHLQAEIVETIDNGFGHSRLVKENLIDYFIPFLPLEYRHRVRLCARDA
 FLSQELLYKEETLDEIAQMMVYVPKEEQLFSSQGCKSISQRINYFLS

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2569_a09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022371

ORF Size: 1191 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022371.4](#)

RefSeq Size: 2138 bp


RefSeq ORF: 1194 bp

Locus ID: 64222

UniProt ID: [Q9H497](#)

Cytogenetics: 1q25.2

Product images:



Circular map for RG200896