

Product datasheet for **RG200906**

SNX11 (NM_152244) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SNX11 (NM_152244) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SNX11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG200906 representing NM_152244 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGCTTTTGGTGTAGGATGTCGGAGAACCAAGAACAGGAGGAGGTGATTACAGTGCCTGTTTCAGGACC
 CCCGAGTGCAGAAATGAGGGCTCCTGGAACCTTATGTGGATTATAAGATATTCCTCCATACCAACAGCAA
 AGCCTTTACTGCCAAGACTTCCTGTGTGCGGCGCCGCTACCGTGAGTTCGTGTGGCTGAGAAAGCAGCTA
 CAGAGAAATGCTGGTTTGGTGCCTGTTCTGAACCTCCTGGGAAGTCAACCTTCTTCGGCACCTCAGATG
 AGTTCATTGAGAAGCGACGACAAGTCTGCAGCACTTCCTTGAAAAGTCTGCAGAGTGTGGTTCTCCT
 GTCAGACAGCCAGTTGCACCTATTCTGCAAAGCCAGCTCTCGGTGCCTGAGATAGAAGCCTGTGTCCAG
 GGCCGAAGTACCATGACTGTGTCTGATGCCATTCTTCGATATGCTATGTCAAAGTGTGGCTGGGCCAGG
 AAGAGAGGCAGAGCTCTTCTCACCTGGCTAAAGGAGACCAGCCTAAGAGTTGCTGCTTTCTTCCAAGATC
 GGGTAGGAGGAGCTCTCCCTCACCGCCTCCAGTGAAGAAAAGGACCATTTAGAAGTGTGGGCTCCAGTT
 GTTGACTCTGAGGTTCTTCTTGAAAGTCCCACTCTCCACCCCTCTCCTCACCATTATGCTGTGATT
 TTGGAAGACCCAAAGAGGAACCTCCACTTTCAGTCTGTGAGGAGGGCTGTGGGAGGAGATCATGCTGT
 GCCTTTGGACCCTGGTCAGTTAGAAACAGTTTGGAAAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >RG200906 representing NM_152244
 Red=Cloning site Green=Tags(s)

MGFWCRMSENQEVEITVRVQDPRVQNEGSWNSYVDYKIFLHTNSKAFTAKTSCVRRRYREFVWLRKQL
 QRNAGLVPVPELPGKSTFFGTSTDEFIEKRRQGLQHFLEKVLQSVLLSDSQLHLFLQSQLSVPEIEACVQ
 GRSTMTVSDAILRYAMSNCGWAQEERQSSSHLAKGDQPKSCCFLPRSGRRSSPSPPPSEEKDHLEVPAPV
 VDSEVPSLESPTLPPLSSPLCCDFGRPKEGTSTLQSVRRVAGGDHAVPLDPGQLETVLEK

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152244

ORF Size: 810 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_152244.1, NP_689450.1</u>
RefSeq Size:	2402 bp
RefSeq ORF:	813 bp
Locus ID:	29916
UniProt ID:	<u>Q9Y5W9</u>
Cytogenetics:	17q21.32
Domains:	PX
Gene Summary:	This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members. This gene encodes a protein of unknown function. This gene results in two transcript variants differing in the 5' UTR, but encoding the same protein. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG200906