

Product datasheet for **RG200920**

LRSAM1 (NM_001005373) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRSAM1 (NM_001005373) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LRSAM1
Synonyms:	CMT2P; RIFLE; TAL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG200920 representing NM_001005373
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCCTCTTCTCCGGAAGCGAAACCCAGTGAGGAGCTCGAAACGCCTGGAGTACCAGATGTGTT
 TGGCAAAAAGAAGCTGGGGCAGATGACATTCTCGACATCTCTAAATGTGAGCTCTCAGAGATTCCATTGG
 AGCTTTTGCAACATGCAAAGTTCTGCAGAAGAAGTGCTGATCGTCCACACGAATCACCTCACTTCCTG
 CTTCCAAATCTGCAGCCTCTGAGTCTGGCAACCATTAAGTTCTAGATCTCCACGATAATCAGCTGA
 CAGCCCTTCTGACGATCTGGGGCAGTGACTGCCCTCCAGGTCTTAAACGTGGAAGGAATCAACTGAT
 GCAGCTCCACGTTCCATTGGGAACCTGACCCAGCTCCAGACTCTCAATGTTAAAGACAACAAGCTGAAG
 GAGCTTCCAGACACCGTGGGGGAGCTTCGAAGCCTGCGTACCCTCAACATCAGTGGAAACGAGATCCAGA
 GATTGCCGCAGATGCTGGCTCACGTTCGAACCTGGAGATGCTGAGCCTTGACGCCTCGGCCATGGTCTA
 CCCGCCGGGAGGTGTGTGGTCCGGCCTGCGGCCATCTTGAGTTCTCTGCAAAGAGTCAGGGCTG
 GAATACTACCCCTTCTCAGTACTTGCTGCCAATTCGGAGCAAGATGGAATCGAGAACTCTCGGGACA
 GCCCTGATGGGCCACCGACAGATTCTCAAGGGAGGAGTTAGAGTGGCAGAACAGGTTCTCAGACTATGA
 GAAGAGGAAGGAACAGAAGATGCTGGAGAACTCGAGTTTGAACGGCGCTGGAAGTGGGGCAGCGGGAG
 CACACCCAGCTCCTTCAGCAGAGCAGCAGCCAGAAGGATGAGATCCTTCAGACGGTCAAGGAGGAGCAGT
 CCCGGCTGGAGCAGGGCTGAGTGAGCACCAGCGCCACCTCGACGCAGAGCGGCAGCGGCTGCAGGAGCA
 GCTGAAGCAGACGGAACAGAACATTTCCAGCCGGATCCAGAAGCTGCTGCAGGACAATCAGAGACAAAAG
 AAAAGTCCGAGATTTTGAATCGCTGGAAAATGAAAGAATAAGAATGGAACAGTTGATGTCCATAACCC
 AGGAGGAGACTGAGAGCCTGCGGGCAGTGACGTTGCCCTCCGCCATGCAGCAGATGCTGACTGAGAGCTG
 TAAGAACCGGCTCATCCAGATGGCCTACGAATCTCAGAGGCAGAACTTGGTCCAGCAGGCTGTTCCAGC
 ATGGCCGAAATGGATGAACGATTCCAGCAGATTCTGTCGTGGCAGCAAATGGATCAGAACAAGCCATCA
 GCCAGATCCTGCAGGAGAGCGCGATGCAGAAGGCTGCGTTCGAGGCACTCCAGGTGAAGAAAGACCTGAT
 GCATCGGCAGATCAGGAGCCAGATTAAGTTAATAGAACTGAGTTATTGCAGCTGACACAGCTGGAGTTA
 AAGAGGAAGTCCCTGGACACAGAGTCACTCCAGGAGATGATCTCGGAGCAGCGCTGGGCCCTCAGCTCCC
 TGCTCCAGCAGCTGCTCAAAGAGAAGCAGCAGCGAGAGGAAGAGCTCCGGGAAATCCTGACGGAGTTAGA
 AGCCAAAAGTGAACACAGGCAGGAAAATTAAGTGGCTGATTCAGTATCAACGGCTTTTGAACCAGAAGCCC
 TTGTCTTGAAGCTGCAAGAAGAGGGGATGGAGCGCCAGCTGGTGGCCCTCCTGGAGGAGCTGTCGGCTG
 AGCACTACCTGCCATCTTTGCGCACCACCGCCTCTCACTGGACCTGCTGAGCCAAATGAGCCCAGGGGA
 CCTGGCCAAAGTGGGCGTCTCAGAAGCTGGCCTGCAGCACGAGATCCTCCGGAGAGTCCAGGAAGTCTG
 GATGCAGCCAGGATCCAGCCAGAGCTGAAACCAATGGGTGAGGTGTCACCCCTACGGCCCCCAGG
 AGCCTCCTGAGTCTGTGAGGCCATCCGCTCCCCCTGCAGAGCTGGAGGTGCAGGCCTCAGAGTGTGTCG
 GTGCTGGAACGGGAGGCCAGATGATCTTCTCAACTGTGGCCACGTCTGCTGCTGCCAGCAGTGTCTGC
 CAGCCACTGCGCACCTGCCCGCTGTGCCCCAGGACATCGCCAGCGCCTCCGCATCTACCACAGCAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG200920 representing NM_001005373
 Red=Cloning site Green=Tags(s)

MPLFFRKRKPSEEARKRLEYQMCLAKEAGADDILDISKCELSSEIPFGAFATCKVLQKKVLIVHTNHLTSL
 LPKSCSLLSLATIKVLDLHDNQLTALPDDLQQLTALQVLNVERNQLMQLPRSIGNLTQLQTLNVKDNKLLK
 ELPDVTGELRSLRTLNI SGNEIQRLPQMLAHVRTLEMLSLDASAMVYPPREVCGAGTAAIQLFLCKESGL
 EYPPSQYLLPILEQDGIENSRDSPDGPTDRFSREELEWQNRFSDEYKKEQKMLEKLEFERRLELQRE
 HTQLLQQSSSQKDEILQTVKEEQSRLEQGLSEHQRHDAERQRLQEQLKQTEQNISSRIQKLLQDNQRQK
 KSSEILKSLENERIRMEQLMSITQEETESLRRRDVASAMQMLTESCKNRLIQMAYESQRQNLVQQACSS
 MAEMDERFQQILSWQMDQNKAISSQILQESAMQKAAFEALQVKKDLMHRQIRSQIKLIETELLQLTQLEL
 KRKSLDTESLQEMISEQRWALSSLLQQLLKEKQQRREEELREILTELEAKSETRQENYWLIIYQRLNQP
 LSLKLQEEGMRQLVALLEELSAEHYLPFAHRLSLDLLSQMSPGDLAKVGVSEAGLQHEILRRVQELL
 DAARIQPELPPMGEVVTPTAPQPPESVRPSAPPAELVQASECVVCLEREAQMIFLNCGHVCCCQCC
 QPLRTPCLCRQDIAQRLRIYHSS


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_001005373

ORF Size: 2169 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005373.1](#), [NP_001005373.1](#)

RefSeq Size: 3152 bp

RefSeq ORF: 2172 bp

Locus ID: 90678

UniProt ID: [Q6UWE0](#)

Cytogenetics: 9q33.3-q34.11

Protein Families: Druggable Genome

Gene Summary: This gene encodes a ring finger protein involved in a variety of functions, including regulation of signaling pathways and cell adhesion, mediation of self-ubiquitylation, and involvement in cargo sorting during receptor endocytosis. Mutations in this gene have been associated with Charcot-Marie-Tooth disease. Multiple transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jan 2012]