

## Product datasheet for **RG200927**

### ZDHHC16 (NM\_032327) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZDHHC16 (NM_032327) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZDHHC16
Synonyms:	APH2; DHHC-16
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG200927 representing NM_032327 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCGAGGCCAGCGGAGCCTGCTGCTGGGCCCGCCGCCTCTGCCTCCGCCTCCTTCTGCTGCTGGGT  
ACAGGGCCCGCTGTCCACCTCTACTCCGGGTCTAGTACAGCGCTGGCGCTACGGCAAGGTCTGCCTGCC  
CTCCCTGCTCTACAACCTCTTTGGGGCAGTGACACCGCTGTTGATGCTGCCTTTGAGCCTGCTACTGG  
CTGGTAGACAACGTGATCCGCTGGTTGGAGTGGTGTTCGTGGTCTGGTATCGTCTGACAGGCTCCA  
TTGTAGCTATCGCTACCTGTGTCTGCCTCTCATCCTCCGAACCTACTCAGTGCCACGACTCTGCTG  
GCATTTCTTCTATAGCCACTGGAATCTGATCCTGATTGTCTTCCACTACTACCAGGCCATCACCCTCCG  
CCTGGGTACCCACCCAGGGCAGGAATGATATCGCCACCGTCTCCATCTGTAAGAAGTGCAATTTACCCCA  
AGCCAGCCCAACACACCACTGCAGCATCTGCAACAGGTGTGTGCTGAAGATGGATCACCCTGCCCTG  
GCTAAACAATTGTGTGGGCCACTATAACCATCGGTACTTCTTCTTTCTGCTTTTTCATGACTCTGGGC  
TGTGTCTACTGCAGCTATGGAAGTTGGGACCTTTCCGGGAGGCTTATGCTGCCATTGAGAAAATGAAAC  
AGCTCGACAAGAACAACACTACAGGCGTTGCCAACCAAGGCTATCACCAGACCCACCACCTTCTC  
CTTTGAGAAAAGGATGACTCACAAAGTCTTGTCTACCTCTGGTTCCTGTGCAGTCTGTGGCACTTGCC  
CTGGGTGCCCTAACTGTATGGCATGCTGTTCTCATCAGTCGAGGTGAGACTAGCATCGAAAGGCACATCA  
ACAAGAAGGAGAGACGTCGGCTACAGGCCAAGGGCAGAGTATTTAGGAATCCTTACAACACTACGGCTGCTT  
GGACAACCTGGAAGGTATTCTGGGTGTGGATACAGGAAGGCACTGGCTTACTCGGGTGTCTTACCTTCT  
AGTCACTTGCCCATGGGAATGGAATGAGCTGGGAGCCCCCTCCCTGGGTGACTGCTCACTCAGCCTCTG  
TGATGGCAGTG

**ACGGT**ACGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG200927 representing NM\_032327  
 Red=Cloning site Green=Tags(s)

MRGQRSLLLGPARLCLRLLLLLGYRRRCPLLRGLVQRWRYGKVCRLSLLYNSFGGSDTAVDAAFEFVYW  
 LVDNVIRWFGVVFVVLVIVLTGSIVAIAYLCVLPILRTYVPRLCWHFFYSHWNLILIVFHYYQAITTP  
 PGYPPQGRNDIATVSICKKCIYPKPARTHHCSICNRCVLKMDHHCPLWLNCCVGHYNHRYFFSFCFFMTLG  
 CVYCSYGSWDLFREAYAAIEKMKQLDKNKLQAVANQTYHQTPPTFSFRERMTHKSLVYLWFLCSSVALA  
 LGALTVWHAVLISRGETSIERHINKERRRLQAKGRVFRNPYNYGCLDNWKVFLGVDTGRHWLTRVLLPS  
 SHLPHNGMSWEPPPPWTAHSASVMAV

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_032327

**ORF Size:** 1131 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_032327.3](#)

**RefSeq Size:** 1816 bp

**RefSeq ORF:** 1134 bp

**Locus ID:** 84287

**UniProt ID:** [Q969W1](#)


**Cytogenetics:** 10q24.1

**Domains:** zf-DHHC

**Protein Families:** Transmembrane

**Gene Summary:** Palmitoyl acyltransferase that mediates palmitoylation of proteins such as PLN and ZDHHC6 (PubMed:28826475). Required during embryonic heart development and cardiac function, possibly by mediating palmitoylation of PLN, thereby affecting PLN phosphorylation and homooligomerization (By similarity). Also required for eye development (By similarity). Palmitoylates ZDHHC6, affecting the quaternary assembly of ZDHHC6, its localization, stability and function (PubMed:28826475). May play a role in DNA damage response (By similarity). May be involved in apoptosis regulation (By similarity). Involved in the proliferation of neural stem cells by regulating the FGF/ERK pathway (By similarity).[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for RG200927