

Product datasheet for **RG200965**

TMEM51 (NM_018022) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM51 (NM_018022) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TMEM51
Synonyms:	C1orf72
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG200965 representing NM_018022 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGATGGCCAGTCCAAGGCCAATGGCTCGCACTATGCGCTGACCGCCATCGGCCTGGGGATGCTGGTCC
TTGGGGTGATCATGGCCATGTGGAACCTGGTACCCGGCTTCAGCGCGCCGAGAAGCCAACAGCTCAGGG
CAGCAACAAGACCGAGGTGGTGGCGGCATCCTCAAGAGCAAGACCTTCTGTGGCCTACGTGCTGGTC
GGGCGCGGGTGATGCTGCTGCTTTTCTATCTGCCTGAGTATCAGGGATAAGAGGAAGCAGCGGCAGG
GCGAGGACCTGGCCATGTCCAGCACCCGACAGGCGCTGGGCCTCACGCCAGGAGGAAGACAGCCAGGA
GGAAGAAGAGGAGGATGAGGAGGCTGCCTCAAGTACTATGTTCCAGCTACGAGGAAGTATGAACACA
AACTACTCAGAAGCAAGGGGAGAGGAGCAGAACCCGAGGTTGAGCATCTCTCTCCCGTCCATGAGTCAC
TGACGGGGCTCGACGAGACCACCCACATCCACCAGGGCTGACGTGGAGGCCAGCCCTGGGAACCCCC
TGACAGGCAGAACTCTAAGTTGGCCAAACGACTGAAACCACTGAAAGTTCGAAGGATTAATCTGAAAAG
CTTCACTCAAAGACTTTAGGATCAACCTCCAGACAAAACGTCCTCCTCCCTCGATAGAGCCTTTGA
CTCCTCCACCGCAGTATGATGAAGTCCAGGAGAAGGCCCGACCCCGCCGCCGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG200965 representing NM_018022
 Red=Cloning site Green=Tags(s)

MMAQSKANGSHYALTAIGLGMLVLGVIMAMWNLVPGFSAAEKPTAQGSNKTEVGGGILKSKTF SVAYVLV
 GAGVMLLLLSICLSIRDKRKQRQGEDLAHVQHPTGAGPHAQEEDSQEEEEDEEAASRYYYVPSYEEVMNT
 NYSEARGEEQNRLSISLPSYESLTGLDETTPTSTRADVEASPGNPPDRQNSKLAKRLKPLKVRRIKSEK
 LHLKDFRINLPDKNVPPPSIEPLTPPPQYDEVQEKAPDTRPPD

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_018022

ORF Size: 759 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018022.1](#), [NP_060492.1](#)

RefSeq Size: 1809 bp

RefSeq ORF: 762 bp


Locus ID: 55092

UniProt ID: [Q9NW97](#)

Cytogenetics: 1p36.21

Protein Families: Transmembrane

Product images:



Circular map for RG200965