

Product datasheet for **RG201029**

SLMO2 (PRELID3B) (NM_016045) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SLMO2 (PRELID3B) (NM_016045) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SLMO2
Synonyms: C20orf45; dj543j19.5; SLMO2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG201029 representing NM_016045
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGATCTGGACTTCGGAGCACGTCTTTGACCACCCGTGGGAACTGTTACAACAGCTGCAATGCAGA
AATACCCAAACCCTATGAACCAAGTGTGGTTGGAGTTGATGTGTTGGACAGACATATAGATCCCTCTGG
AAAGTTGCACAGCCACAGACTTCTCAGCACAGAGTGGGGACTGCCTTCCATTGTGAAGTCTTTATTGGT
GCAGCAAGAACGAAAACATATGTGCAAGAACATTCTGTAGTTGATCCTGTAGAGAAAACAATGGAACCTA
AATCTACTAATATTTTCAATTTACAACATGGTTTCAGTAGATGAGAGACTTATATACAACCACATCCTCA
GGATCCAGAAAAAAGTGTGACACAAGAAGCCATAATTACCGTGAAAGGAGTTAGCCTCAGCAGTTAC
CTTGAAGGACTGATGGCAAGTACGATATCCTCAAATGCTAGTAAAGGCCGAGAAGCAATGGAATGGGTAA
TACATAAATTAATGCTGAGATTGAAGAACTGACAGCCTCAGCAAGAGGAACCATAAGGACTCCAATGGC
AGCAGCAGCGTTTGCAGAGAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG201029 representing NM_016045
Red=Cloning site Green=Tags(s)

MKIWTSEHVFDPWETVTTAAMQKYPNPMNPSVVGVDVLDRIHIDPSGKLSHRLSTEWGLPSIVKSLIG
AARTKTYVQEHVVDPEKTMELKSTNISFTNMVSDERLIYKPHPQDPEKTVLTQEAIIIVKGVSLSSY
LEGLMASTISSNASKGREAMEWVIHKLNAEIEELTASARGTIRTPMAAAFAEK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_016045

ORF Size: 582 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016045.3](#)

RefSeq Size: 2581 bp

RefSeq ORF: 585 bp


Locus ID: 51012

UniProt ID: [Q9Y3B1](#)

Cytogenetics: 20q13.32

Domains: MSF1

Product images:



Circular map for RG201029