

Product datasheet for **RG201062**

MED18 (NM_017638) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MED18 (NM_017638) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: MED18
Synonyms: p28b; SRB5
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG201062 representing NM_017638
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGCACCTCCAGTCACCATGATGCCTGTCACTGGGGGCACCATTAACATGATGGAGTACCTGTTGC
 AGGGAAGTGTTTAGATCACAGTTTGAAAGCCTCATCCACCGCCTTCGTGGTTGTGTGACAACATGGA
 ACCTGAGACTTTCCTTGACCATGAGATGGTATTCTCCTTAAGGGCCAGCAAGCCAGCCCATTTGTTCTC
 AGGGCCCAGCCTCTATGGACAGGGCAGGGGCACCCTGGCATCTGCGCTACCTGGGACAGCCAGAAATGG
 GAGACAAGAACCAGCCATGCCCTGGTGCAGAACTGCGTGGACATTGCCACATCTGAGAACCTCACCAGCTT
 CTTGATGGAAATGGGCTTCCGCATGGACCATGAGTTTGTGCTAAGGGACATTTGTTCCGTAAGGGCATC
 ATGAAGATTATGGTGTACAAGATTTCCGCATCCTGGTGGCAGGGAACACAGACAGCACTGAGGCCTTGT
 CACTCTCTATCTCGTGGAAATTAAGTGTGGTAGCACCCGCTGGGCAGGACATGGTCTCTGATGACATGAA
 GAACCTCGCAGAACAGCTAAACCTCTGGTTCACCTAGAGAAAATAGACCCCAAGAGGCTCATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG201062 representing NM_017638
 Red=Cloning site Green=Tags(s)


MEAPPVTMPVTGGTINMMEYLLQGSVLDHSLES LIHRLRGLCDNMEPETFLDHEMVFLKGGQASPFVL
 RARRSMDRAGAPWHLRYLGQPEMGDKNRHALVRNCVDIATSENLTDFLMEMGFRMDHEFVAKGHLFRKGI
 MKIMVYKIFRILVPGNTDSTEALSLSYLVELSVVAPAGQDMVSDDMKNFAEQLKPLVHLEKIDPKRLM

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



View online »

Cloning Scheme:


ACCN: NM_017638

ORF Size: 624 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017638.3](#)

RefSeq Size: 1845 bp

RefSeq ORF: 627 bp


Locus ID: 54797

UniProt ID: [Q9BUE0](#)

Cytogenetics: 1p35.3

Gene Summary:

MED18 is a component of the Mediator complex, which is a coactivator for DNA-binding factors that activate transcription via RNA polymerase II (Sato et al., 2003 [PubMed 12584197]).[supplied by OMIM, Oct 2008]

Product images:


Circular map for RG201062