

## Product datasheet for **RG201230**

### PDHB (NM\_000925) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PDHB (NM_000925) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PDHB
Synonyms:	PDHBD; PDHE1-B; PDHE1B; PHE1B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG201230 representing NM_000925 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGCGGTGTCTGGCTTGGTGCAGAGACCCCTTCGGGAGGTCTCCGGCTGCTGAAGAGGCGCTTTC  
ACTGGACCGCGCCGGCTGCGCTGCAGGTGACAGTTCGTGATGCTATAAATCAGGGTATGGATGAGGAGCT  
GGAAAGAGATGAGAAGGTATTTCTGCTTGGAGAAGAAGTTGCCAGTATGATGGGCATACAAGGTTAGT  
CGAGGGCTGTGGAAGAAATATGGAGACAAGAGGATTATTGACACTCCCATATCAGAGATGGGCTTTGCTG  
GAATTGCTGTAGGTGCAGCTATGGCTGGTTGCGGCCATTTGTGAATTTATGACCTTCAATTTCTCCAT  
GCAAGCCATTGACCAGGTATAAACTCAGCTGCCAAGACCTACTACATGTCTGGTGGCCTTCAGCCTGTG  
CCTATAGTCTTCAGGGGGCCCAATGGTGCCTCAGCAGGTGTAGCTGCCAGCACTCACAGTGCTTTGCTG  
CCTGGTATGGGCACTGCCAGGCTTAAAGGTGGTCACTCCCTGGAATTCAGAGGATGCTAAAGGACTTAT  
TAAATCAGCCATTCGGGATAACAATCCAGTGGTGGTGTAGAGAAATGAATTGATGTATGGGGTCCTTTT  
GAATTTCTCCGGAAGCTCAGTCAAAGATTTCTGATTCCTATTGAAAAAGCCAAAATAGAAAAGCAAG  
GAACACATATAACTGTGGTTTCCATTCAAGACCTGTGGCCACTGCTTAGAAGCTGCAGCAGTGCATC  
TAAAGAAGGAGTTGAATGTGAGGTGATAAATATGCGTACCATTAGACCAATGGACATGGAACCATAGAA  
GCCAGTGCATGAAGACAAATCATCTTGAAGTGTGGAAGGAGGCTGGCCACAGTTTGGAGTAGGAGCTG  
AAATCTGTGCCAGGATCATGGAAGTCTGCGTTCAATTTCTGGATGCTCCTGCTGTTCTGCTACTGG  
TGCTGATGCCCTATGCCTTATGCAAAGATTCTAGAGGACAACCTTATACCTCAGGTCAAAGACATCATA  
TTTGCAATAAAGAAAACATTAATATT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG201230 representing NM\_000925  
 Red=Cloning site Green=Tags(s)

MAAVSGLVRRPLREVSGLLKRRFHWTAPAALQVTVRDAINQGMDEELERDEKVFLLGEEVAQYDGAYKVS  
 RGLWKKYGDKRIIDTPISEMGFAGI AVGAAMAGLRPICEFMTFNF SMQAIDQVINSAAKTYMSSGGLQPV  
 PIVFRGPNASAGVAAQHSQCFAAWYGHCPGLKVVSPWNSDAGLIKSAIRDNNPVVLENELMYGVPP  
 EFPPEAQSKDFLIPIGKAKIERQGTHTITVVSHSRPVGHCL EAAAVLSKEGVECEVINMRTIRPMDMETIE  
 ASVMKTNHLVTVEGGWPQFGVGAEICARIMEGPAFNFLDAPAVRVTGADVPMPYAKILEDNSIPQVKDII  
 FAIKKTLNI


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_000925

<b>ORF Size:</b>	1077 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_000925.3</a> , <a href="#">NP_000916.2</a>
<b>RefSeq Size:</b>	1544 bp
<b>RefSeq ORF:</b>	1080 bp
<b>Locus ID:</b>	5162
<b>UniProt ID:</b>	<a href="#">P11177</a>
<b>Cytogenetics:</b>	3p14.3
<b>Domains:</b>	transket_pyr, transketolase_C
<b>Protein Pathways:</b>	Butanoate metabolism, Citrate cycle (TCA cycle), Glycolysis / Gluconeogenesis, Metabolic pathways, Pyruvate metabolism, Valine, leucine and isoleucine biosynthesis
<b>Gene Summary:</b>	The pyruvate dehydrogenase (PDH) complex is a nuclear-encoded mitochondrial multienzyme complex that catalyzes the overall conversion of pyruvate to acetyl-CoA and carbon dioxide, and provides the primary link between glycolysis and the tricarboxylic acid (TCA) cycle. The PDH complex is composed of multiple copies of three enzymatic components: pyruvate dehydrogenase (E1), dihydrolipoamide acetyltransferase (E2) and lipoamide dehydrogenase (E3). The E1 enzyme is a heterotetramer of two alpha and two beta subunits. This gene encodes the E1 beta subunit. Mutations in this gene are associated with pyruvate dehydrogenase E1-beta deficiency. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2012]