

Product datasheet for **RG201239**

CKMT1A (NM_001015001) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CKMT1A (NM_001015001) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CKMT1A
Synonyms:	CKMT1; mia-CK; U-MtCK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG201239 representing NM_001015001 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGTCCCTTCTCCCGTCTGCTGTCCGCCGCCCGGGACTCAGGCTCCTGGCTTTGGCCGGAGCGG
GGTCTCTAGCCGCTGGGTTTCTGCTCCGACCGAACCTGTACGAGCTGCCAGTGAACGACGGAGGCTGTA
TCCCCCGAGCGCTGAGTACCCAGACCTCCGAAAGCACAACAACCTGCATGGCCAGTCACCTGACCCAGCA
GTCTATGCACGGCTCTGCGACAAGACCACCCACTGGTTGGACGCTAGATCAGTGTATCCAGACTGGCC
TGGACAACCCTGGCCACCCTTTCATCAAGACTGTGGGCATGGTGGCTGGAGATGAGGAGACCTATGAGGT
ATTTGCTGACCTGTTTGACCCTGTGATCCAAGAGCGACACAATGGATATGACCCCGGACAATGAAGCAC
ACCACGGATCTAGATGCCAGTAAATCCGTTCTGGCTACTTTGATGAGAGGTATGTATTGCTCCTCTAGAG
TCAGAACTGGCCGAAGCATCCGAGGACTCAGTCTGCCTCCAGCTTGCCTCGAGCAGAGCGACGAGAGGT
GGAACGTGTTGTGGTGGATGCACTGAGTGGCCTGAAGGGTGACCTGGCTGGACGTTACTATAGGCTCAGT
GAGATGACAGAGGCTGAACAGCAGCAGCTTATTGATGACCACTTTCTGTTTGATAAGCCTGTGCCCCGT
TGCTGACTGCAGCAGGAATGGCTCGAGACTGGCCAGATGCTCGTGAATTTGGCACAACAATGAGAAGAG
CTTCTGATCTGGGTGAATGAGGAGGATCATACCGGGTGATCTCCATGGAGAAGGGTGGTAACATGAAG
AGAGTGTGAAAGATTCTGCCGAGGCCCAAAGAGGTGGAGAGACTTATCCAAGAACGCTGGCTGGGAGT
TCATGTGGAATGAGCGTTTGGGATACATCTTGACCTGTCCATCTAACCTGGGCACCTGGACTTCGGGCAGG
AGTGACATCAAACCTGCCCTGCTAAGCAAAGATAGCCGCTTCCCAAAGATCCTGGAGAACCTAAGACTC
CAAAGCGTGGTACTGGAGGAGTGGACTGCTGCCACAGGCGGTGCTTTGATATTTCTAATTTGGACC
GACTAGGCAAAATCAGAGGTGGAGCTGGTGAACCTGGTCAATCGATGGAGTAAACTATTTGATTGATTGTGA
ACGGCGTCTGGAGAGAGGCCAGGATATCCGCATCCCCACACCTGTCATCCACACCAAGCAT

ACGGCTACGGCGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG201239 representing NM_001015001
 Red=Cloning site Green=Tags(s)

MAGPF S R L L S A R P G L R L L A L A G A G S L A A G F L L R P E P V R A A S E R R R L Y P P S A E Y P D L R K H N N C M A S H L T P A
 V Y A R L C D K T T P T G W T L D Q C I Q T G V D N P G H P F I K T V G M V A G D E E T Y E V F A D L F D P V I Q E R H N G Y D P R T M K H
 T T D L D A S K I R S G Y F D E R Y V L S S R V R T G R S I R G L S L P P A C T R A E R R E V E R V V D A L S G L K G D L A G R Y Y R L S
 E M T E A E Q Q L I D D H F L F D K P V S P L L T A A G M A R D W P D A R G I W H N N E K S F L I W V N E E D H T R V I S M E K G G N M K
 R V F E R F C R G L K E V E R L I Q E R G W E F M W N E R L G Y I L T C P S N L G T G L R A G V H I K L P L L S K D S R F P K I L E N L R L
 Q K R G T G G V D T A A T G G V F D I S N L D R L G K S E V E L V Q L V I D G V N Y L I D C E R R L E R G Q D I R I P T P V I H T K H


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001015001

ORF Size:	1251 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001015001.2 , NP_001015001.1
RefSeq Size:	1779 bp
RefSeq ORF:	1254 bp
Locus ID:	548596
UniProt ID:	P12532
Cytogenetics:	15q15.3
Protein Families:	Druggable Genome
Protein Pathways:	Arginine and proline metabolism, Metabolic pathways
Gene Summary:	Mitochondrial creatine (MtCK) kinase is responsible for the transfer of high energy phosphate from mitochondria to the cytosolic carrier, creatine. It belongs to the creatine kinase isoenzyme family. It exists as two isoenzymes, sarcomeric MtCK and ubiquitous MtCK, encoded by separate genes. Mitochondrial creatine kinase occurs in two different oligomeric forms: dimers and octamers, in contrast to the exclusively dimeric cytosolic creatine kinase isoenzymes. Many malignant cancers with poor prognosis have shown overexpression of ubiquitous mitochondrial creatine kinase; this may be related to high energy turnover and failure to eliminate cancer cells via apoptosis. Ubiquitous mitochondrial creatine kinase has 80% homology with the coding exons of sarcomeric mitochondrial creatine kinase. Two genes located near each other on chromosome 15 have been identified which encode identical mitochondrial creatine kinase proteins. [provided by RefSeq, Jul 2008]