

## Product datasheet for **RG201264**

### SSH3BP1 (ABI1) (NM\_001012752) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SSH3BP1 (ABI1) (NM_001012752) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SSH3BP1
Synonyms:	ABI-1; ABLBP4; E3B1; NAP1BP; SSH3BP; SSH3BP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG201264 representing NM\_001012752  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCAGAGCTGCAGATGTTACTAGAGGAGGAGATCCCGTCTGGCAAGAGGGCGCTGATCGAGAGTTACC  
 AGAACCTGACTCGGGTGGCAGACTACTGTGAAAACAACTACATACAGGCTACAGACAAGAGAAAAGCTTT  
 AGAGGAGACCAAAGCCTATACAACCCAATCTCTAGCTAGTGTGCTTATCAAATAAATGCATTGGCCAAC  
 AATGTAAGTCCAGTTGCTGGATATCCAAGCCTCTCAGCTTCGGAGAATGGAGTCTCCATCAATCATATCT  
 CACAGACTGTGGATATTCATAAGGAGAAAGTGGCAGCAAGAGAGATTGGTATTTTGACAACAAATAAGAA  
 TACATCAAGAAGTACAAAATAATAGCACCTGCGAATATGGAGCGCCCTGTAAGGTATATTCGGAACCT  
 ATCGATTACAGTTCGGATGATGTGGCCATGGTGTCAAGCATGGAATAACCCAGCCTGCAAGAAGT  
 GCACACTGTCGAGAACAAATCCTCCTACTCAGAAACCGCCAAGTCTCCATGTCAGGCCGGGGAACACT  
 GGGACGGAATACTCCTTATAAAACCTGGAACCTGTTAAACCCCAACAGTTCCTAATGACTATATGACC  
 AGTCTGTAGGCTTGAAGTCAGCATAGTCCAGGCAGGACAGCATCTTTAAATCAGAGACCAAGGACAC  
 ACAGTGAAGTAGTGGAGGAAGTGAAGTCGAGAAAAACAGTGGTAGCAGTAGTATTGGCATTCCCATTGC  
 TGTGCCTACACCTTCGCCACCCACTATTGGACCAGCAGCCCCGGGCTCAGCTCCTGGTCCCAGTATGGC  
 ACAATGACCAGGCAGATATCTCGACACAACCTCTACTCTTCTCGACATCTTCTGGTGGATACAGACGAA  
 CTCCCTCTGTGACTGCTCAATTTCTGCTCAGCCTCATGTTAATGGAGGTCCACTTTATTCTCAAAATTC  
 AATTTCTATTGCTCCACCCCTCCCTATGCCTCAGTTGACTCCACAGATACCTCTCACAGGCTTCGTG  
 GCCAGGGTGCAGGAAAACATTGCTGATAGTCCAACCTCCACCGCCACCCTCCACAGATGACATCCCA  
 TGTTGATGACTCTCACCTCCCCACCACCAGTGGATTATGAAGATGAGGAGGCTGCAGTACT  
 TCAGTATAATGATCCATATGCAGATGGGATCCTGCTTGGGCCCAAGAATTATATTGAGAAAGTTGTT  
 GCAATATATGATTATACAAAAGACAAGGATGATGAGCTGTCATTTATGGAGGGTGAATCATTTATGTTA  
 TAAAGAAGAATGATGATGGCTGGTATGAAGGAGTCTGCAATCGAGTACTGGTCTGTTCCCTGGGAATA  
 TGTTGAATCAATCATGCACTATACTGAT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:**

>RG201264 representing NM\_001012752  
 Red=Cloning site Green=Tags(s)

MAELQMLLEEEIPSGKRALIESYQNLTRVADYCENNYIQATDKRKALEETKAYTTQSLASVAYQINALAN  
 NVLQLLDIQASQLRRMESSINHISQTVDIHKEKVARREIGILTNNKNTSRTHKIIAPANMERPVRYIRKP  
 IDYTVLDDVGHGVKHGNNQPARTGTLSRNPPQKPPSPPMSSRGTGRNTPYKTLPEVKPPTVPNDYMT  
 SPARLGSQHSRGRTASLNQRPRTHSGSSGSGSRENSGSSSIGIPIAVPTSPPTIGPAAPGSAPGSQYG  
 TMTRQISRHNSTTSSTSSGGYRRTPSVTAQFSAQPHVNGGPLYSQNSISIAPPPPMPQLTPQIPLTGFV  
 ARVQENIADSPTPPPPPPDDIPMFDDSPPPPPPPVDYEDEEAADVQYNDPYADGDPWAWPKNYIEKVV  
 AIYDYTKDKDDELSFMEGAIYVIKKNDDGWYEGVGNRVTGLFPGNYVESIMHYTD

**TR**TRPLE - GFP Tag - V

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**


**ACCN:** NM\_001012752

**ORF Size:** 1428 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001012752.2](#)

**RefSeq Size:** 3572 bp

**RefSeq ORF:** 1431 bp


**Locus ID:** 10006

**UniProt ID:** [Q8IZP0](#)

**Cytogenetics:** 10p12.1

**Gene Summary:** This gene encodes a member of the Abelson-interactor family of adaptor proteins. These proteins facilitate signal transduction as components of several multiprotein complexes, and regulate actin polymerization and cytoskeletal remodeling through interactions with Abelson tyrosine kinases. The encoded protein plays a role in macropinocytosis as a component of the WAVE2 complex, and also forms a complex with EPS8 and SOS1 that mediates signal transduction from Ras to Rac. This gene may play a role in the progression of several malignancies including melanoma, colon cancer and breast cancer, and a t(10;11) chromosomal translocation involving this gene and the MLL gene has been associated with acute myeloid leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 14. [provided by RefSeq, Sep 2011]

**Product images:**



Circular map for RG201264