

Product datasheet for **RG201308**

MMTAG2 (C1orf35) (NM_024319) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MMTAG2 (C1orf35) (NM_024319) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MMTAG2
Synonyms:	MMTAG2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG201308 representing NM_024319 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTCCGGCTCCAGTCGTGGAGGCGTGC CGCGGGCAGGACCAGTTCAACTGGGAGGACGTGAAGACTG
ACAAGCAGCGGGAGAACCTGGGCAACTCGCTGATGGCGCCGGTAGGCCGCTGGCAGAAGGGCCGCGA
CCTCACCTGGTACGCCAAGGGCCGGCGCCATGCGCGGGCCGAGCCGCGAGGAGGAACCTGGCAGCCGTG
CGGGAGGCGGAGCGGAGGCGCTGCTGGCCGCCCTTGGCTACAAGAACGTGAAGAAGCAGCCACGGGCC
TGAGCAAGGAGGACTTCGCGGAGGTCTGCAAGCGGGAAGGAGGCGACCCCGAGGAGAAGGGCGTGGACCG
GCTGCTGGGGCTGGGAGCGCAAGTGGCTCCGTGGGCCGCGTGGCGATGTCCCGAGAGGACAAGGAGGCC
GCCAAACTGGGGCTGTCTGTGTTACGCATCACCGCTAGAGAGCGGCGGGCCCGGGACCTCGGCAGCCT
CGGCCAGGAGGAAGCCGCGGGCGGAGGATCAGACGGAAAGCAGTTGTGAGAGCCACAGGAAAAGCAAGAA
GGAGAAGAAGAAAAGAAAAGAGGAAACACAAGAAAGAGAAGAAGAAGAAAGACAAAGAGCACAGGCGG
CCAGCTGAGGCCACCTCCTCTCCACATCTCCTGAGAGGCCAGGCACCACCACCATGACTCCGACTCCA
ACTCCCCTGCTGTAAGAGGAGGAAGCGGGACACAGTGGGACAGGAGGAGCCCGTCTCGCAGTGGCA
TGACAGAGGCTCTGAGGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG201308 representing NM_024319
Red=Cloning site Green=Tags(s)

MFGSSRGGVRRGGQDFNWEDVKTDKQRENYLGNLMAPVGRWQKGRDLTWYAKGRAPCAGPSREEELAAV
 REAEREALLAALGYKNVKKOPTGLSKEDFAEVCKREGGDPEEKGVDRLLGLGSASGVSGRVAMSREDKEA
 AKLGLSVFTHHRVESGGPGTSAASARRKPRAEDQTESSCESHRKSKKEKKKKKKRKHKKKEKKKKDKEHRR
 PAEATSSPTSPERPRHHHSDSDSNPCKRRKRGHSGDRRSPRRWHDGSEA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024319

ORF Size: 789 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024319.4](#)

RefSeq Size: 1307 bp


RefSeq ORF: 792 bp

Locus ID: 79169

UniProt ID: [Q9BU76](#)

Cytogenetics: 1q42.13

Product images:



Circular map for RG201308