

Product datasheet for **RG201469**

CEP57 (NM_014679) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CEP57 (NM_014679) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CEP57
Synonyms:	MVA2; PIG8; TSP57
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG201469 representing NM_014679
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGCGGCGCGTCTGTCTCTGCGGCTTCTGGTTCTCACTTGTGCAACAGCTTTGCTGAGCCATCAAGGT
 CTAATGGAAGCATGGTTGCGCATTCTTCATCTCCATATGTAGTATATCCTTCGGATAAGCCTTTCCTTAA
 TAGTGATCTACGACGCTCCCAAGTAAGCTACACTTGCCTATCCAGAAAGCAACAGCAGAGCCATATTT
 TCTGCTCTTAAGAATCTTCAAGATAAGATTGACGCTTGAACCTGAGAGGATTGAGCAGAGAAAGT
 TGAAAACCTTGTCTAGAGAAACAATTGAATATAAGAAAGTACTGGATGAACAGATACAAGAAAGGAGAA
 TTCAAAGAATGAGGAATCAAAGCACAATCAAGAACTGACATCTCAGTTGTTAGCTGCAGAAAATAATGC
 AATCTATTAGAAAAACAATTGGAATACATGCGAAATATGATAAAGCATGCCGAAATGGAGAGGACATCTG
 TCTTAGAGAAACAAGTTTCCCTAGAAAGAGAACGACAACATGATCAAACACATGTTGAGGCCAACTTGA
 AAAATTGGATCTTCTTGAACAGGAGTATAACAACTTACCACAATGCAGGCCCTTGCAGAAAAAAATG
 CAAGAGTTGGAAGCAAACTCCATGAAGAAGAACAGGAAAGGAAACGCATGCAAGCTAAGGCAGCTGAGT
 TGCAGACTGGTCTAGAAACAAATAGACTTATCTTTGAAGATAAGGCAACTCCGTGTGTTCCCAATGCAAG
 AAGAATTAAGAAAGAAAGTCAAAACCACAGAAAAGAAAGTTCTAGGAACTATTTTGGTGCAACCA
 CATTATAGATTATGCTTGGGTGATATGCCATTTGTAGCTGGGAAGTCCACAAGCCCTAGCCATGCCGTGG
 TAGCCAATGTTGAGCTTGTCTTGCATCTAATGAAGCAACACAGTAAAGCTTTGTGCAATGATCGAGTCAT
 CAACAGTATTCCTTTGGCAAAGCAAGTATCTTACGAGGTGGTAAAGTAAGAAGTTGTCAGTAACACCT
 CCCTCCTCAACGGTATTAATGAGGAGTTGTCAGAAGTCTTACAGACTTTACAGGATGAATTTGGGCAAA
 TGAGCTTTGATCACCAGCAGCTTGCAAACTTATCCAGGAGTCGCAACCGTTGAACTGAAAGACAAGTT
 GGAGTGTGAATTGGAGGCATTAGTGGAAGGATGGAAGCAAAAGCCAACCAATAACTAAAGTTTCGAAAA
 TACCAAGCCCAGCTGGAGAAACAGAAGTTAGAGAAGCAGAAGAAGGAATTAAGCTACCAAAAAGACTC
 TTGATGAAGAAAGAAACAGCAGCAGCCGTTCTGGAATCACAGGGACCACAAATAAGAAAGATTTTATGAA
 ACTGAGACCTGGAGAAAAAGGAGAAAAATCTTCAGTTATTGAAGGACATGCAAGCATACAGAATTCA
 TTCAAAGCAGTAGTTTGTGTTGGGATTAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG201469 representing NM_014679
 Red=Cloning site Green=Tags(s)

MAAASVSAASGSHLSNSFAEPSRSNGSMVRHSSSPYVVYPSPDKPFLNSDLRRSPSKPTLAYPESNSRAIF
 SALKNLQDKIRRELERIQAEESVKTLSTETIEYKKVLDEQIQERENSKNEESKHNLQELTSQLLAENKC
 NLLEKQLEYMRNMIKHAEMERTSVLEKQVSLERERQHDQTHVQSLEKLDLLEQYNNKLTMTQALAEKMM
 QELEAKLHEEEQERKRMQAKAAELQTGLETNRLIFEDKATPCVPNARRIKKKSKPPEKKSSRNYFGAQP
 HYRLCLGDMPPFVAGKSTSPSHAVVANVQLVLHLMKQHSKALCNDRVINSIPLAKQVSSRGKSKKLSVTP
 PSSNGINEELSEVLQTLQDEFQMSFDHQLAKLIQESPTVELKDKLECELEALVGRMEAKANQITKVRK
 YQAQLEKQKLEKQKELKATKKTLDERNSSSRSGITGTTNKKDFMKLRPGEKRRKNLQLLKDMQSIQNS
 LQSSSLCWDY

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:




ACCN:	NM_014679
ORF Size:	1500 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_014679.5
RefSeq Size:	3158 bp
RefSeq ORF:	1503 bp
Locus ID:	9702

UniProt ID: [Q86XR8](#)

Cytogenetics: 11q21

Gene Summary: This gene encodes a cytoplasmic protein called Translokin. This protein localizes to the centrosome and has a function in microtubular stabilization. The N-terminal half of this protein is required for its centrosome localization and for its multimerization, and the C-terminal half is required for nucleating, bundling and anchoring microtubules to the centrosomes. This protein specifically interacts with fibroblast growth factor 2 (FGF2), sorting nexin 6, Ran-binding protein M and the kinesins KIF3A and KIF3B, and thus mediates the nuclear translocation and mitogenic activity of the FGF2. It also interacts with cyclin D1 and controls nucleocytoplasmic distribution of the cyclin D1 in quiescent cells. This protein is crucial for maintaining correct chromosomal number during cell division. Mutations in this gene cause mosaic variegated aneuploidy syndrome, a rare autosomal recessive disorder. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]

Product images:



Circular map for RG201469