

Product datasheet for **RG201581**

CFDP1 (NM_006324) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CFDP1 (NM_006324) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CFDP1
Synonyms:	BCNT; BUCENTAUR; CENP-29; CP27; p97; SWC5; Yeti
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG201581 representing NM_006324 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGAATTCGACTCCGAAGACTTCTCTACGTCGGAGGAGGACGAGGACTACGTGCCGTCGGGTGGAG
 AGTATAGTGAAGATGATGTAATGAATTAGTGAAGGAAGATGAAGTGGATGGTGAAGAGCAGACACAGAA
 AACCCAAAGGGAAGAAAGGAGCCAGAGCATTCCAGCCAGGAAGAGAAGACAAGGTGGCCTCTCATT
 GAAGAAGAGGAAGAGGAGGATGCCAATTCAGAATCTGAGGGAAGCAGTAGTGAGGAGGAAGATGACGCTG
 CAGAGCAGGAAAAAGGCATTGGATCAGAGGATGCCAGGAAAAAGAGGAGACGAACCTCTGGCCAGCTT
 CCTCAATGATGTGGGACCAAAATCAAAAGTGCCCAAGTACACAAGTTAAGAAAGGAGAGGAGACTGAA
 GAGACAAGTTCAAGTAAATTGTTGGTAAAAAGCAGAAGAGCTAGAGAAACCTAAAGAAACAGAAAAAGTTA
 AAATCACCAAGGTGTTTGATTTTGCTGGTGAAGAAGTAAGGGTAACTAAGGAAGTGGATGCTACATCTAA
 AGAGGCCAAATCCTTCTCAAGCAGAAATGAGAAAGAAAAACACAGGCTAATGTTCTTCAGCTCTGCCA
 TCACTCCCTGCCGGGTCAGGGTTAAAAAGATCAAGTGGCATGAGCAGCCTTTTGGGAAAAATTGGTGCCA
 AGAAGCAGAAAAATGAGCACCTTGAGAAGTCCAACTGGACTGGGAGAGCTTCAAGGAGGAAGAGGGGAT
 TGGTGAAGAACTGGCCATCCATAATCGAGGGAAGAGGGGTACATTGAACGGAAGCCTTCCTCGACCGA
 GTGGATCACAGGCAGTTTGAAATTGAGCGAGATCTCAGGCTGAGCAAAATGAAACCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >RG201581 representing NM_006324
 Red=Cloning site Green=Tags(s)

MEEFDSEDFSTSEEDDYVPSGGEYSEDDVNLVKEDEVGEEQTQKTQGKKRKAQSIPARKRRQGGLSL
 EEEEEEDANSESEGSSEEDDAAEQEKIGSEDARKKKEDELWASFLNDVGPKSKVPPSTQVKKGEETE
 ETSSSKLLVKAEELEKPKETEKVKITKVFDFAGEEVRVTKEVDATSKAKSFFKQNEKEKPQANVPSALP
 SLPAGSGLKRSSGMSSLLGKIGAKKQKMSTLEKSKLDWESFKEEEGIGEELAIHNRGKEGYIERKAFLDR
 VDRHQFEIERDLRLSKMKP

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_006324

ORF Size: 897 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_006324.1](#)

RefSeq Size: 1281 bp

RefSeq ORF: 900 bp


Locus ID: 10428

UniProt ID: [Q9UEE9](#)

Cytogenetics: 16q23.1

Gene Summary: May play a role during embryogenesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG201581