

## Product datasheet for **RG201589**

### HMGN4 (NM\_006353) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** HMGN4 (NM\_006353) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** HMGN4  
**Synonyms:** HMG17L3; NHC  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG201589 representing NM\_006353  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCCCAAGAGAAAGGCAAAAGGAGATGCTAAAGGTGATAAGCAAAGGTGAAGGATGAGCCACAGAGGA  
 GATCAGCTCGTTGTCTGCTAAACCAGCTCCTCCAAACCAGAGCCAGGCCTAAAAAGGCTCTGCAAA  
 GAAGGGAGAGAAGCTTCCCAAAGGGAGAAAGGGAAAGCAGATGCTGAAAGGATGGAAACAACCTGCA  
 AAAAACCGAGATGCCTCTACACTCCAGTCCAGAAAGCGGAAGGCACTGGGGATGCCAAG

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >RG201589 representing NM\_006353  
 Red=Cloning site Green=Tags(s)  
 MPKRKAKGDAKGDKAKVKDEPQRRSARLSAKPAPPKPEPRPKKASAKKGEKLPKGRKGKADAGKDGNP  
 KNRDASTLQSQAEGTGD

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_006353

**ORF Size:** 270 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM\\_006353.1](#), [NP\\_006344.1](#)

RefSeq Size: 1980 bp

RefSeq ORF: 273 bp

Locus ID: 10473

UniProt ID: [O00479](#)

Cytogenetics: 6p22.2

**Gene Summary:** The protein encoded by this gene, a member of the HMGN protein family, is thought to reduce the compactness of the chromatin fiber in nucleosomes, thereby enhancing transcription from chromatin templates. [provided by RefSeq, Mar 2013]

## Product images:



Circular map for RG201589