

## Product datasheet for **RG201601**

### **SULT1A1 (NM\_001055) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SULT1A1 (NM_001055) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SULT1A1
Synonyms:	HAST1/HAST2; P-PST; P-PST 1; PST; ST1A1; ST1A3; STP; STP1; ts-PST; TSPST1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG201601 representing NM_001055 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGAGCTGATCCAGGACACCTCCCGCCGCCACTGGAGTACGTGAAGGGGTCCCGCTCATCAAGTACT  
TTGCAGAGGCACTGGGGCCCTGCAGAGCTTCCAGGCCCGCCTGATGACCTGCTCATCAGCACCTACCC  
CAAGTCCGGCACTACCTGGTAAGCCAGATTCTGGACATGATCTACCAGGGTGGTACCTGGAGAAGTGT  
CACCGAGCTCCCATCTTCATGCGGGTGCCTTCCCTTGAGTTCAAAGCCCCAGGGATCCCTCAGGGATGG  
AGACTCTGAAAGACACACCGGCCCCACGACTCCTGAAGACACACCTGCCCTGGCTCTGCTCCCCAGAC  
TCTGTTGGATCAGAAGGTCAAGGTGGTCTATGTTGCCCGCAACGCAAAGGATGTGGCAGTTTCTACTAC  
CACTTCTACCACATGGCCAAGGTGCACCCTGAGCCTGGGACCTGGGACAGCTTCCCTGGAGAAGTTCATGG  
TCGGAGAAGTGTCTACGGATCCTGGTACCAGCACGTGCAGGAGTGGTGGGAGCTGAGCCGACCCACCC  
TGTTCTCTACCTCTTCTATGAAGACATGAAGGAGAACCCGAAAAGGGAGATTCAAAGATCTCTGGAGTTT  
GTGGGGCACTCCCTGCCAGAGGAGACCGTGGACTTCGTGGTTCAGCACACGTCGTTCAAGGAGATGAAGA  
AGAACCCTATGACCAACTACACCACCGTCCCCAGGAGTTCATGGACCACAGCATCTCCCCCTTCATGAG  
GAAAGGCATGGCTGGGGACTGGAAGACCACCTTACCCTGGCGCAGAATGAGCGCTTCGATGCGGACTAT  
GCGGAGAAGATGGCAGGCTGCAGCCTCAGCTTCCGCTCTGAGCTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG201601 representing NM\_001055  
Red=Cloning site Green=Tags(s)

MELIQDTSRPPLEYVKGVPLIKYFAEALGPLQSFQARPDDL ISTYPKSGTTWVSQILDMIYQGGDLEKC  
 HRAPIFMRVPFLEFKAPGIPSGMETLKDTPAPRLLKTHLPLALLPQTLLDQKVKVYVARNADVAVSYY  
 HFYHMAKVHPEPGTWDSFLEKFMVGEVSYGWSYQHVQEWELSRTHPVL YLFYEDMKENPKREIQILEF  
 VGHSLPEETVDFVVQHTSFKEMKKNPMTNYTTVPQEFMDHSISPFMRKGMAGDWKTTFTVAQNERFDADY  
 AEKMAGCSLSFRSEL

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001055

**ORF Size:** 885 bp


**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001055.2</a> , <a href="#">NP_001046.2</a>
<b>RefSeq Size:</b>	1222 bp
<b>RefSeq ORF:</b>	888 bp
<b>Locus ID:</b>	6817
<b>UniProt ID:</b>	<a href="#">P50225</a>
<b>Cytogenetics:</b>	16p11.2
<b>Domains:</b>	Sulfotransfer
<b>Protein Pathways:</b>	Sulfur metabolism
<b>Gene Summary:</b>	Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. This gene encodes one of two phenol sulfotransferases with thermostable enzyme activity. Multiple alternatively spliced variants that encode two isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

### Product images:



Circular map for RG201601