

## Product datasheet for **RG201718**

### **TMEM106C (NM\_024056) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TMEM106C (NM\_024056) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** TMEM106C  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG201718 representing NM\_024056  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGGTCTCAGCATTCCGCTGCTGCTCGCCCCTCCTCCTGCAGGCGAAAGCAAGAAGATGACAGGGACG  
 GTTTGCTGGTGAACGAGAGCAGGAAGAAGCCATTGCTCAGTCCCATATGTGGAATTCACCGGGAGAGA  
 TAGCATCACCTGTCTCACGTGCCAGGGGACAGGCTACATCCAACAGAGCAAGTAAATGAGTTGGTGGCT  
 TTGATCCCACACAGTGATCAGAGATTGCGCCCTCAGCGAACTAAGCAATATGTCCTCCTGCCATCCTGC  
 TTTGTCTCCTGGCATCTGGTTTGGTGGTTTTCTTCCTGTTCCGCATTAGTCCTTGGGATGATGACCG  
 CATCAAAGTGGTAAAGTCACATTTAATAAGCAAGACTCCCTTGTAAATTCACCATCATGGCCACCCTG  
 AAAATCAGGAACTCCAATTCTACACGGTGGCAGTGACCAGCCTGTCCAGCCAGATTCAGTACATGAACA  
 CAGTGGTCAGTACATATGTGACTACTAACGTCTCCCTATTCCACCTCGGAGTGAGCAACTGGTGAATTT  
 TACCGGGAAGGCCGAGATGGGAGGACCCTTTTCCATGTGTACTTCTTCTGCACGGTACCTGAGATCCTG  
 GTGCACAACATAGTGATCTTCATGCGAACTTCAGTGAAGATTTACATATTGGCCTCATGACCCAGAGCT  
 CCTTGGAGACACATCACTATGTGGATTGTGGAGAAATTCACAGCTATT

**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA**

**Protein Sequence:** >RG201718 representing NM\_024056  
 Red=Cloning site Green=Tags(s)

MGSQHSAAARPSSCRKQEDDRDGLLAEREQEEIAQFPYVEFTGRDSITCLTCQGTGYIPTEQVNELVA  
 LIPHSQRLRPQRTKQYVLLSILLCLLASGLVFFLPHSVLVDDGKIKVVKVTFNKQDSLVIITIMATL  
 KIRNSNFYTVAVTSLSSQIQYMNTVVSTYVTTNVSLIPRSEQLVNFTGKAEMGGPFYVYFFCTVPEIL  
 VHNIVIFMRTSVKISYIGLMTQSSLETHHYVDCGGNSTAI


**TRTRPLE - GFP Tag - V**



[View online »](#)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_024056

**ORF Size:** 750 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_024056.3](#)

**RefSeq Size:** 1449 bp

**RefSeq ORF:** 753 bp


**Locus ID:** 79022

**UniProt ID:** [Q9BVX2](#)

Cytogenetics: 12q13.11

Protein Families: Transmembrane

**Product images:**



Circular map for RG201718