

Product datasheet for **RG201743**

IKK gamma (IKBKG) (NM_003639) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IKK gamma (IKBKG) (NM_003639) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IKK gamma
Synonyms:	AMCBX1; EDAID1; FIP-3; FIP3; Fip3p; IKK-gamma; IKKAP1; IKKG; IMD33; IP; IP1; IP2; IPD2; NEMO; ZC2HC9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG201743 representing NM_003639
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAATAGGCACCTCTGGAAGAGCCAAGTGTGTGAGATGGTGCAGCCAGTGGTGGCCCGCAGCAGATC
 AGGACGTACTGGGCGAAGAGTCTCCTCTGGGGAAGCCAGCCATGCTGCACCTGCCTCAGAACAGGGCGC
 TCCTGAGACCCTCCAGCGCTGCCTGGAGGAGAATCAAGAGCTCCGAGATGCCATCCGGCAGAGCAACAG
 ATTCTGCGGGAGCGCTGCGAGGAGCTTCTGCATTTCCAAGCCAGCCAGAGGGAGGAGAAGGAGTTCTCA
 TGTGCAAGTTCCAGGAGGCCAGGAACTGGTGGAGAGACTCGGCCTGGAGAAGCTCGATCTGAAGAGGCA
 GAAGGAGCAGGCTCTGCGGGAGGTGGAGCACCTGAAGAGATGCCAGCAGCAGATGGCTGAGGACAAGGCC
 TCTGTGAAAGCCCAGGTGACGTCTTGTCTCGGGAGCTGCAGGAGAGCCAGAGTCGCTTGGAGGCTGCCA
 CTAAAGGAATGCCAGGCTCTGGAGGTCGGGCCCGGGCCAGCAGCAGCGCGCCAGCTGGAGAGTGA
 GCGCGAGGCGCTGCAGCAGCAGCACAGCTGCAGGTGGACCAGCTGCGCATGCAGGGCCAGAGCGTGGAG
 GCCGCGCTCCGCATGGAGCGCCAGGCCGCTCGGAGGAGAAGAGGAAGCTGGCCAGTTGCAGGTGGCCT
 ATCACCAGCTCTTCCAAGAATACGACAACCACATCAAGAGCAGCGTGGTGGGCGAGTGAGCGGAAGCGAGG
 AATGCAGCTGGAAGATCTCAAACAGCAGCTCCAGCAGGCCAGGAGGCCCTGGTGGCCAAACAGGAGGTG
 ATCGATAAGCTGAAGGAGGAGCCGAGCAGCACAAGATTGTGATGGAGACCGTTCCGGTGTGAAGGCC
 AGGCGGATATCTACAAGCGGACTTCCAGGCTGAGAGGCGAGCCCGGGAGAAGCTGGCCGAGAAGAAGGA
 GCTCCTGCAGGAGCAGCTGGAGCAGCTGCAGAGGAGTACAGCAAACCTGAAGGCCAGCTGTCAGGAGTCG
 GCCAGGATCGAGGACATGAGGAAGCGGCATGTCGAGGTCTCCAGGCCCTTGCCTCCCGCCCTGCCT
 ACCTCTCTCTCCCTGGCCCTGCCAGCCAGAGGAGGAGCCCCCGGAGGAGCCACCTGACTTCTGTG
 TCCAAGTGCCAGTATCAGGCCCTGATATGGACACCCTGCAGATACATGTCATGGAGTGCATTGAG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG201743 representing NM_003639
 Red=Cloning site Green=Tags(s)

MNRHLWKSQ LCEMVQPSGGPAADQDVLGEESPLGKPAMLHLPSEQGAPETLQRCLLENQELRDAIRQSNQ
 ILRERCEELLHFQASQREEKEFLMCKFQEARKLVERLGLEKLDLKRQKEQALREVEHLKRCQQQMAEDKA
 SVKAQVTSLLGELQESQSRLEAATKECQALEGRARAASEQARQLESEREALQQQHSVQVDQLRMQGSVE
 AALRMRQAASEEKRLAQLQVAYHQLFQEYDNHIKSSVVGSRKRMQLEDLKQQLQQAEEALVAKQEV
 IDKLKEEAQHKIVMETVPVLKAQADIYKADFQAERQAREKLAEKKELLQEQLQREYSKLGKASCQES
 ARIEDMRKRHVEVSQAPLPPAPAYLSSPLALPSQRRSPPEEPPDFCCPKCQYQAPMDTLQIHVMECIE

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_003639

ORF Size: 1257 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_003639.4](#)


RefSeq Size: 2121 bp

RefSeq ORF: 1260 bp

Locus ID: 8517

UniProt ID: [Q9Y6K9](#)
Cytogenetics: Xq28
Domains: zf-C2H2
Protein Families: Druggable Genome, Transcription Factors
Protein Pathways: Acute myeloid leukemia, Adipocytokine signaling pathway, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Cytosolic DNA-sensing pathway, Epithelial cell signaling in Helicobacter pylori infection, MAPK signaling pathway, NOD-like receptor signaling pathway, Pancreatic cancer, Pathways in cancer, Primary immunodeficiency, Prostate cancer, RIG-I-like receptor signaling pathway, Small cell lung cancer, T cell receptor signaling pathway, Toll-like receptor signaling pathway
Gene Summary: This gene encodes the regulatory subunit of the inhibitor of kappaB kinase (IKK) complex, which activates NF-kappaB resulting in activation of genes involved in inflammation, immunity, cell survival, and other pathways. Mutations in this gene result in incontinentia pigmenti, hypohidrotic ectodermal dysplasia, and several other types of immunodeficiencies. A pseudogene highly similar to this locus is located in an adjacent region of the X chromosome. [provided by RefSeq, Mar 2016]

Product images:



Circular map for RG201743